



Feb 02, 2025

SeV fluorescent reporter virus titration by TCID50

 Journal of Virology

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1d73yvr2/v1>

Carolina Lopez¹

¹Washington University



Carolina Lopez

Washington University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1d73yvr2/v1>

External link: <https://doi.org/10.1128/jvi.00894-25>

Protocol Citation: Carolina Lopez 2025. SeV fluorescent reporter virus titration by TCID50. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1d73yvr2/v1>

Manuscript citation:

Pye SE, Achouri E, Yang Y, Musa A, López CB Validation of diverse and previously untraceable Sendai virus copyback viral genomes by direct RNA sequencing. Journal of Virology 99(8). doi: [10.1128/jvi.00894-25](https://doi.org/10.1128/jvi.00894-25)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 14, 2025

Last Modified: February 02, 2025

Protocol Integer ID: 118345

Keywords: sev fluorescent reporter virus titration by tcid50 protocol, sev fluorescent reporter virus titration, tcid50 protocol, virus, rsev

Abstract

Protocol to titer rSeV carrying a fluorescent reporter

Materials

1. Cultured LLCMK2 cells (ATCC, #CCL-7)
2. Tissue culture media and Infection media (See media prepare protocol)
3. 0.05% Trypsin (Thermo, #25300054)
4. PBS (Phosphate-buffered saline)
5. Cell counter and Trypan blue
6. 96-well plates
7. Pipettes, 1.5 ml tubes
8. TPCK-Trypsin (for Sendai virus infection, see SeV infection protocol)

Seed cells to 96-well plate

- 1
 1. Check LLCMK2 cells to ensure they are in good condition, approximately 90% confluent.
 2. Wash cells twice with PBS.
 3. Add  3 mL 0.05% trypsin-EDTA to the T75 flask, incubate at room temperature for 3-4 minutes.
 4. Add  7 mL TCM to the flask to stop trypsinization and collect cells.
 5. Mix  10 μ L cells and  10 μ L Trypan blue, then add  10 μ L mixture to a cell counter slide and count the cells.
 6. Plate 20,000 cells in  100 μ L of TCM media to each well. 3 rows of a 96-well plate are needed for one test.
 7. Incubate plates at  37 °C overnight.

Notes:

- Usually, LLCMK2 cells are used for SeV titration since most subsequent experiments are done with LLCMK2. However, if you plan to infect A549 cells or other cell lines with the virus stocks, it is recommended to do the titration on the same cell line you intend to use.
- The seeding number for LLCMK2 and A549 cells is similar, usually around 20,000 cells per well. This results in 95–100% confluence overnight, though this can vary depending on the condition of the cells. In my experiments, seeding numbers between 12,000 and 20,000 cells per well have been successfully used.

Titration

- 2
 1. Prepare infection media: Add TPCK-trypsin to the infection media at a final concentration of 2 μ g/mL. Prepare 30 mL of the media for each plate, ensuring sufficient volume for dilutions and replacing the TCM from the plate. (This step is for LLCMK2 cells, if you are using A549 cells, TPCK-trypsin is not needed).
 2. Use the prepared media to do a 10-fold serial dilution from 1:10¹ to 1:10⁸: Prepare 8 tubes. Add 450 μ L of prepared media to each tube. Add 50 μ L of virus to the first tube, vortex to mix thoroughly, then transfer 50 μ L to the second tube. Repeat this process sequentially until the 8th tube.
 3. Check the cells under a microscope to confirm they are healthy and evenly distributed.
 4. Remove the TCM from the plate using a multichannel pipette.
 5. Wash the cells by adding 200 μ L of PBS to each well, then discard the PBS.
 6. Add 100 μ L of the prepared infection media (containing TPCK-trypsin) to each well.



7. Add 100 μ L of the diluted virus to the cells, with each dilution to three wells (e.g., columns A-C as below).

8. Incubate the plate at 37 °C for 3–5 days.

	A	B	C	D	E	F	G	H	I	J	K	L
	Stoc k1 1:10 ¹	Stoc k1 1:10 ¹	Stoc k1 1:10 ¹	Stoc k2 1:10 ¹	Stoc k2 1:10 ¹	Stoc k2 1:10 ¹	...					
	Stoc k1 1:10 ²	Stoc k1 1:10 ²	Stoc k1 1:10 ²	Stoc k2 1:10 ²	Stoc k2 1:10 ²	Stoc k2 1:10 ²						
	Stoc k1 1:10 ³	Stoc k1 1:10 ³	Stoc k1 1:10 ³	Stoc k2 1:10 ³	Stoc k2 1:10 ³	Stoc k2 1:10 ³						
	Stoc k1 1:10 ⁴	Stoc k1 1:10 ⁴	Stoc k1 1:10 ⁴	Stoc k2 1:10 ⁴	Stoc k2 1:10 ⁴	Stoc k2 1:10 ⁴						
	Stoc k1 1:10 ⁵	Stoc k1 1:10 ⁵	Stoc k1 1:10 ⁵	Stoc k2 1:10 ⁵	Stoc k2 1:10 ⁵	Stoc k2 1:10 ⁵						
	Stoc k1 1:10 ⁶	Stoc k1 1:10 ⁶	Stoc k1 1:10 ⁶	Stoc k2 1:10 ⁶	Stoc k2 1:10 ⁶	Stoc k2 1:10 ⁶						
	Stoc k1 1:10 ⁷	Stoc k1 1:10 ⁷	Stoc k1 1:10 ⁷	Stoc k2 1:10 ⁷	Stoc k2 1:10 ⁷	Stoc k2 1:10 ⁷						
	Stoc k1 1:10 ⁸	Stoc k1 1:10 ⁸	Stoc k1 1:10 ⁸	Stoc k2 1:10 ⁸	Stoc k2 1:10 ⁸	Stoc k2 1:10 ⁸						

Read titration results

3 Read the results at 4dpi as below:



	A	B	C
	Positive wells of last row with signal	Dilution of last row with signal	TCID ₅₀ / 100uL
	+ + +	10-X	10 ^X + 0.7
	+ + -	10-X	10 ^X + 0.4
	+ - -	10-X	10 ^X - 0.1

X = row in which last signal occurs

+ = Fluorescent signal positive

Notes:

- The TCID₅₀ result is reported per 100 µL. To convert it to per mL, plus 1 to the X. For example, a TCID₅₀ of 10^{7.7}/100 µL is equivalent to 10^{8.7} TCID₅₀/mL.
- For rSeV-C-eGFP, the TCID₅₀ results remain consistent from 3 dpi to 5 dpi. So, the results can be read on any of these days.