

Nov 18, 2025

SeV titration using FFU

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzjxj5lx1/v1>

Carolina Lopez¹

¹Washington University



Carolina Lopez

Washington University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzjxj5lx1/v1>

Protocol Citation: Carolina Lopez 2025. SeV titration using FFU. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzjxj5lx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: November 18, 2025

Protocol Integer ID: 100402

Keywords: ffu protocol for titration, sev titration, protocol for titration, using ffu protocol, titration, sev

Abstract

Protocol for titration of SeV

Materials

Infection medium

- DMEM 500 mL
- Pen/Strep 5 mL
- 35% BSA 5 mL
- 5% NaHCO₃ 12 mL

2X Infection medium

- DMEM F12 (Powder) 1 sachet
- Pen/Strep 10 mL
- L/Glut10 mL
- 35% BSA 6 mL
- HEPES 1M 10 mL
- NaHCO₃ 18 mL
- Sterile dH₂O to 500 mL

Overlay medium

- 2X Infection media
- 2.5% Sterile Avicel
- TPCK-Trypsin

Reagents for Staining

- Anti-SeV primary antibody (example: ms SeV NP antibody, Clone M73/2)
- Goat anti-mouse H+L Biotinylated Antibody (Southern Biotech #1036-08)
- Streptavidin AP (ThermoFisher #21324)
- 1-StepTM NBT/BCIP Substrate Solution (ThermoFisher #34042)



Preparing Overlay Medium

- 1 To prepare the overlay medium, first, add the correct amount of TPCK to the 2X infection medium. Then mix the 2X infection media/TPCK with the same volume of 2.5% Avicel to prepare a 1:1 solution. Alternately, 0.6% Agarose can be used in a similar 1:1 manner for a final concentration of 0.3% Agarose. All media must be prepared in sterile conditions.

Example for  50 mL of Overlay medium using Avicel:

-  25 mL of 2x infection medium
-  50 μ L of TPCK (2 mg/mL)
-  25 mL of 2.5% Avicel

Day 0 - SEEDING CELLS

- 2
 - Seed 2×10^5 LLCMK-2 cells/well in a 24-well plate (or in a way they are 100% confluent next day).
 - *Suggestion:*
 - Trypsinize a whole, 80% confluent T75 flask;
 - Resuspend cells in  10 mL (regular 10% DMEM) and seed  450 μ L each well - Homogenize and check if cells are evenly distributed in the wells.
 - If scaling up to 48-well plates, seed 10^5 cells/well.
 - Incubate plates at  37 °C .

Day 1 - DILUTION AND INFECTION

1h

- 3
 - Check cells and ensure they are 100% confluent.
 - Prepare 10fold dilutions of your samples in Infection Media;
 - Change tips for every dilution;
 - Wash cells twice with PBS to remove any FBS trace;
 - Inoculate  200 μ L of your serially diluted samples on each well (if using 48-well plates inoculate  180 μ L).
 - This step can be done in duplicates or triplicates.
 - Inoculate in the opposite direction from the dilution so there is no need for changing tips;
 - Save wells for the negative control (Infection media only) and the 2ry antibody negative control (Also infection media only).
 - Incubate at  37 °C for  01:00:00 , gently rock the plate every 10';

1h



- Remove the inoculum and add  500 μL of Overlay medium + TPCK 2 $\mu\text{g}/\text{mL}$ (250 μL if using 48-well plates).
- Incubate the plates for 2 days at  37 $^{\circ}\text{C}$.

Day 3 (48hpi) - FIXATION

2h

- 4
 - Aspirate the overlay medium from the wells and fix cells with 10% Formaldehyde (in PBS) for  02:00:00 . If using overlay medium containing Agarose, directly add 500 μL 4% PFA over the wells and incubate as follows.
 - 200 $\mu\text{L}/\text{well}$
 - Use the plate rocker to avoid drying
 - Remove the fixative agent and wash twice with PBS.
 - Store at  4 $^{\circ}\text{C}$ or move directly to the staining steps.

2h

Day 3 - STAINING

1h

- 5
 - Remove the PBS and replace with  150 μL of PBS/1% BSA 1%/0.1% Saponin.
 - Incubate for 30' to 1hr to block/permeabilize cells

1h

OBS.: All the following steps are done with 150 to 180 $\mu\text{L}/\text{well}$ of reagent. Scale up according to the number of wells. To achieve better results also keep the plates on the rocker for all the following incubation steps.

- Aspirate the blocking/permeabilization solution and incubate with primary antibody;
 - Unlabeled **Ms anti-SeV NP (1.3 mg/mL) M73/2**
 - 1:2000 diluted in PBS/1% BSA 1%/0.1% Saponin;

OBS.: Don't add the primary antibody to the 2ry-only negative control well.

- Incubate for  01:00:00 - rocking

- Wash 3x with PBS ( 180 μL of PBS/well, aspirate carefully for not scraping the monolayers).
- Incubate with secondary antibody for 30';
 - Biotinylated Gt anti-Mouse (H+L) diluted 1:1000 in PBS/1% BSA 1%/0.1% Saponin.
- Wash 3x with PBS
- Incubate with Streptavidin for 30'
 - AP-Conjugated Streptavidin, diluted 1:1000 in PBS/1% BSA 1%/0.1% Saponin.
- Wash 3x with PBS
- Incubate with AP substrate for up to 20' (or until the FFUs are pretty stained).



- 1-Step NBT/BCIP Substrate solution for AP

- Wash 1X with dH₂O.

Day 3 – FOCUS COUNTING

- 6
 - In the last dilution with visible focus, count the stained focus in all the wells from that particular sample and divide by the number of replicates.
 - The titer will be expressed by the average number of focus counted times the dilution factor, per volume of initial inoculum.




Figure 1. Representative titration plate. 10-fold dilutions of the initial virus inoculum are indicated in red. Infection media was used as negative control (NC), and is indicated in the bottom-right. Fixed NC monolayer was either incubated with the secondary (2ry) antibody only, or both primary and secondary antibodies as indicated.

- Example: In the plate above, there are 5 FFUs in the -7 dilution. Considering the dilution factor, the titer would be 5×10^7 FFUs. Finally considering the initial volume inoculated (200 uL), the final titer of that particular sample is 5×10^7 FFU/200uL.

Protocol references

Annotated by IAC