

Aug 26, 2023

shRNA knockdown

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwdjj7lmk/v1>

wusj¹, Nancy C. Hernandez Villegas², schekman¹

¹Department of Molecular and Cell Biology, Howard Hughes Medical Institute, University of California, Berkeley, Berkeley, United States;

²Helen Wills Neuroscience Institute, University of California, Berkeley, Berkeley, United States



Nancy C. Hernandez Villegas

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwdjj7lmk/v1>

Protocol Citation: wusj , Nancy C. Hernandez Villegas, schekman 2023. shRNA knockdown. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwdjj7lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 04, 2023

Last Modified: August 26, 2023

Protocol Integer ID: 85962

Keywords: dnajc5 knockdown, shrna knockdown this protocol, shrna knockdown, using shrna, stable cell line

Funders Acknowledgements:

Nancy C. Hernandez Villegas

Grant ID: NIH training program T32GM139780

Abstract

This protocol describes a standard procedure to generate a stable cell line for a DNAJC5 knockdown using shRNA.

Materials

Reagents

	Item	Catalog number	Manufacturer
	HEK 293T cells		Cell Culture Facility, UC Berkeley
	pIKO.1-Hygroplasmids		
	DNAJC5-shRNA		
	pMD2.G plasmid		Addgene
	PsPAX2 plasmid		Addgene
	Lenti-X Concentrator	631231	Takara Bio
	SH-SY5Y cells		Cell Culture Facility, UC Berkeley
	Hygromycin B (50 mg/mL)	10687010	Thermo Scientific



shRNA knockdown

- 1 HEK293T cells were plated at a density such that on the day of transfection they are no more than 50% confluent.
- 2 pIKO.1-Hygroplasmids-containing shRNA targeting DNAJC5 (ccggGCAACCTCAGATGACATTAACTCGAGTTTAATGTCATCTGAGGTTGCTTTTGG) together with pMD2.G and PsPAX2 were transfected into HEK293T cells to produce lentiviral particles for 72 hr.
- 3 Lentivirus particles were concentrated with Lenti-X Concentrator follow manufacturer protocol (Takara Bio).
- 4 SH-SY5Y cells were transduced by lentivirus before differentiation.
- 5 Three days post transduction, cells were selected with 250 µg/ml hygromycin for 10 days.
- 6 The selected cells were differentiated, and the knockdown was verified with immunoblot.