

Nov 02, 2018 Version 1

Single-cell mapping of lineage and identity via CellTagging V.1

DOI

dx.doi.org/10.17504/protocols.io.u9yez7w

Samantha A Morris¹, Brent A. Biddy¹, Wenjun Kong¹, Kenji Kamimoto¹, Chuner Guo¹, Sarah Waye¹, Tao Sun¹

¹Washington University, Saint Louis

Morris Lab



Samantha A Morris

Brigham and Women's Hospital; Harvard Medical School

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.u9yez7w

Protocol Citation: Samantha A Morris, Brent A. Biddy, Wenjun Kong, Kenji Kamimoto, Chuner Guo, Sarah Waye, Tao Sun 2018. Single-cell mapping of lineage and identity via CellTagging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.u9yez7w>

Manuscript citation:

<https://www.biorxiv.org/content/early/2017/04/28/127860>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: November 02, 2018

Last Modified: November 02, 2018

Protocol Integer ID: 17432

Keywords: cell barcoding, lineage tracing, single-cell analysis, direct lineage reprogramming



Abstract

Single-cell technologies are offering unprecedented insight into complex biology, revealing the behavior of rare cell populations that are typically masked in bulk population analyses. One current limitation of single-cell approaches is that lineage relationships are lost as a result of cell processing, restricting interpretations of the data collected. Elegant computational approaches have been developed in an effort to infer these missing observations, but it remains a challenge to reconstruct true reprogramming trajectories using these tools. Although sophisticated lineage tracing solutions to connect cell history with fate are emerging, these protocols are either not compatible with high-throughput scRNA-seq, or require genome editing strategies that are not readily deployed in some systems. Here, our protocol describes a single-cell resolution clonal tracking approach, 'CellTagging', based on combinatorial cell indexing, permitting the parallel capture of lineage information and cell identity. CellTagging integrates with high-throughput single-cell RNA-sequencing, where iterative rounds of cell labeling enable the construction of multi-level lineage trees. This straightforward lentiviral-labeling approach can be applied to an array of cell biological applications to simultaneously profile lineage and identity, at single-cell resolution.

Materials

MATERIALS

⊗ Stellar Competent Cells **Takara Bio Inc. Catalog #636763**

⊗ CellTag pooled library V1 **addgene Catalog #115643**





⊗ CellTag pooled library V2 **addgene Catalog #115644**

⊗ CellTag pooled library V3 **addgene Catalog #115645**

Safety warnings

! For generation of lentivirus, the appropriate institutional regulations must be satisfied.

Amplification of CellTag library

- 1
Thaw Stellar Competent Cells in an ice bath just before use.
- 2
After thawing, mix 100 μ L of cells with 10-50 ng of CellTag V1 DNA in a 1.5-mL microcentrifuge tube.
- 3
Place transformation mixture on ice for 30 minutes.
- 4
Heat shock the cells for exactly 60 seconds at  42 °C .
- 5
Place tube on ice for 1-2 minutes.
- 6
Add SOC medium to bring the final volume to 1000 μ L.
- 7
Incubate by shaking (~250 rpm) for 1 hour at  37 °C .
- 8
Take 6 μ L of the recovery. Use 1 μ L of this to prepare a 1:10 and then a 1:100 dilution. Plate both dilutions and the rest of the 5 μ L of recovery onto three separate LB/Amp plates. Grow overnight at  37 °C .
- 9
Add the rest of the recovery to 500 mL of LB+Amp. Grow overnight while shaking (~250 rpm) at  37 °C .
- 10
After counting the CFUs from the LB/Amp plates the next day, calculate the transformation efficiency.
- 11
Harvest the cells from the liquid culture and use multiple Maxiprep columns (Sigma or Zymogen) to purify the library.

12

