



Aug 01, 2024

Version 1

Spore based infection assay on *Pinus sylvestris* seedlings with *Diplodia sapinea* V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyw9revx9/v1>

Anne Oostlander¹, Laura Brodde², Miriam von Bargaen¹, Rasmus Enderle³, Marco Leiterholt¹, Dagmar Trautmann³, Malin Elfstrand², Jan Stenlid², André leißner¹

¹Institute of Genetics, Technical University Braunschweig, Braunschweig, Germany;

²Department of Forest Mycology and Plant Pathology, Swedish University of Agricultural Sciences, Uppsala, Sweden;

³Institute of Forest Protection, Julius Kuehn Institute (JKI), Braunschweig



Anne Geertje Oostlander

Technische Universität Braunschweig

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyw9revx9/v1>



Protocol Citation: Anne Oostlander, Laura Brodde, Miriam von Bargen, Rasmus Enderle, Marco Leiterholt, Dagmar Trautmann, Malin Elfstrand, Jan Stenlid, André leißner 2024. Spore based infection assay on *Pinus sylvestris* seedlings with *Diplodia sapinea*. protocols.io <https://dx.doi.org/10.17504/protocols.io.eq2lyw9revx9/v1>

Manuscript citation:

Oostlander AG, Brodde L, Bargen M von, Leiterholt M, Trautmann D, Enderle R et al. A reliable and simple method for the production of viable pycnidiospores of the Pine pathogen *Diplodia sapinea* and a spore-based infection assay on Scots Pine. Plant Dis 2023; 107(11):3370–7.
doi.org/10.1094/PDIS-01-23-0107-RE

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 26, 2024

Last Modified: August 01, 2024

Protocol Integer ID: 104147

Keywords: *Diplodia sapinea*, *Diplodia* tip blight, infection assay, *Pinus sylvestris* infection, pycnidiospores, *Sphaeropsis sapinea*, *diplodia sapinea* pathogenicity in scots pine seedling, infection assay on *pinus sylvestris* seedling, assessing *diplodia sapinea* pathogenicity, *pinus sylvestris* seedling, scots pine seedling, *diplodia sapinea*, *diplodia sapinea* this protocol, *pinus sylvestri*, spore, infection assay

Abstract

This protocol describes a spore based method for assessing *Diplodia sapinea* pathogenicity in Scots pine seedlings (*Pinus sylvestris*), including inoculation, symptom evaluation, and reisolation.

Materials

Material:

Petri dishes

D. sapinea (e.g. ex-type strain CBS 138184)

Light shelf with cold white daylight (e.g. Osram Lumilux 18W/865), intensity from 5000 to 6500 lux

Cleanbench

Centrifuge

Sterile scalpel

Ethanol (70%)

Sodium hypochlorite (3% NaOCl)

Sterile water

Plant Material:

Scots pine (*P. sylvestris*) container seedlings

Potting soil

Plastic pots

Minimal medium (VMM):

as described in Vogel, 1964

20 ml Vogel's solution

20 g sucrose

15 g agar

add 1 l water autoclave

Trace element solution:

50 g citric acid

50 g zinc sulfate (J.T. Baker)

10 g ammonium iron(II) sulfate

2,5 g copper sulfate

0,5 g manganese sulfate

0,5 g boric acid

0,5 g sodium molybdate

add 1 l water

Chloroform (1 ml) is added as a preservative, store at room temperature.

Vogel's solution:

125 g sodium citrate

250 g potassium dihydrogen phosphate

100 g ammonium nitrate

10 g magnesium sulfate

5 g calcium chloride



5 ml trace element solution

2,5 ml biotin solution

add 1 l water

Chloroform (2 ml) is added as a preservative, store at room temperature.

Biotin solution:

0,1 g biotin

add 1 l water

Store at -20°C

0,01 % (v/v) Tween:

0,1 ml Tween

add 1 l water

aliquod and autoclave

Plant Material and Greenhouse Conditions

- 1 Use 2 year old container seedlings of Scots pine (*Pinus sylvestris*).
- 2 Ensure seedlings have no visible symptoms and have undergone two annual cycles.
- 3 Apply the last fungicide treatment 8 weeks before the first inoculation.
- 4 Replant seedlings in plastic pots filled with potting soil
- 5 Maintain plants at 20–25°C, 16 h of light per day, and high humidity (>90% RH) for 4 d after inoculation, then switch to moderate humidity (60% RH).

Fungal Inoculum Preparation

- 6 Grow *Diplodia sapinea* on VMM (Vogels Minimal Medium) for 21 d at approximately 27°C under constant light (5000–6500 lx) to induce sporulation.
- 7 Harvest the spores from the plate by adding approximately 2 ml of 0.01% (v/v) Tween to the plate and rinse the surface of the colony several times by pipetting. Repeat the procedure for a higher yield.
- 8 Dilute the spore suspension to 2×10^6 spores/ml for inoculation.
- 9 Assess spore viability by spreading 600 µl of the spore suspension onto three VMM plates. Incubate for 7 h at 27°C and record the germination status of 200 spores per plate. Aim for an average germination rate of 88%.
- 10 Confirm the absence of hyphal fragments microscopically.

Inoculation

- 11 Wounded or not wounded plants can be inoculated. Wounding is likely to influence symptom development.



12 Wounding: Use a sterile scalpel to make a 5 mm long cut down to the cambium on dormant seedlings and on the last year's growth segments of actively growing seedlings. Be careful that the cut does not become too deep.

13 Inoculation with Spores:

13.1 Pipetting: Apply 2 μ l of spore suspension (approx. 4000 spores) directly onto the wound or unwounded stem.

13.2 Spraying: Spray approximately 360 μ l (approx. 720,000 spores) from a distance of 10 cm onto the wound or unwounded stem.

13.3 Control Treatments: Use sterile 0.01% Tween 20 for mock inoculation.

Symptom Assessment

14 Evaluate the symptoms 4- and 6-weeks post-inoculation based on the classification into previously determined symptom classes, for example:

1. No symptoms
2. Necrotic needles, no symptoms on stem
3. Upper third of shoot necrotic
4. Upper two thirds necrotic
5. Seedling dead

Sampling for reisolation

15 After 6 weeks, cut off and discard the shoot tip.

16 Remove an approx. 1.5 cm long piece of the stem above the wound (if applicable), remove needles, and use for reisolation.

Surface Sterilization

17 Sterilize stem pieces:



- 17.1 Immerse in 70% ethanol for 30 s.
- 17.2 Immerse in 3% NaOCl for 60 s.
- 17.3 Rinse twice in sterile water for 30 s each.
- 17.4 Dry on sterile filter paper.

Isolation and Identification

- 18 Cut four approx. 2 mm long pieces from the center of the stem samples and place on MYP agar. Incubate for 27 d at room temperature and ambient daylight.
- 19 Sort isolates into morphological groups from day 4 of cultivation.
- 20 Extract fungal DNA of a representational number of isolates for molecular identification.
- 21 Use species-specific primers to identify *D. sapinea* (for example as described in Adamson et al. 2021).

Protocol references

Adamson, Kalev/Laas, Marili/Blumenstein, Kathrin/Busskamp, Johanna/Langer, Gitta J./Klavina, Darta/Kaur, Anu/Maaten, Tiit/Mullett, Martin S./Müller, Michael M./Ondrušková, Emília/Padari, Allar/Pilt, Enn/Riit, Taavi/Solheim, Halvor/Soonvald, Liina/Tedersoo, Leho/Terhonen, Eeva/Drenkhan, Rein (2021). Highly clonal structure and abundance of one haplotype characterise the *Diplodia sapinea* populations in Europe and Western Asia. Journal of fungi (Basel, Switzerland)
<https://doi.org/10.3390/jof7080634>.