



Feb 15, 2024

# Stable Transfection of Plasmid DNA into Adherent Rodent Cell Lines Using Calcium Phosphate

DOI

[dx.doi.org/10.17504/protocols.io.rm7vzxz5rgx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzxz5rgx1/v1)

Vipasha Dwivedi<sup>1</sup>, Jen Davis<sup>1</sup>, Ella Wilson<sup>1</sup>, Markelle Scott<sup>1</sup>, Madison Zelan<sup>1</sup>, Nicholas DeFeo<sup>1</sup>, Victoria Huffman<sup>1</sup>, Rees Cooke<sup>1</sup>, Kaitlyn O'Daniel<sup>1</sup>, David J Riese II<sup>1,2</sup>

<sup>1</sup>Auburn University; <sup>2</sup>University of Alabama-Birmingham



David J Riese II

Auburn University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzxz5rgx1/v1>

**Protocol Citation:** Vipasha Dwivedi, Jen Davis, Ella Wilson, Markelle Scott, Madison Zelan, Nicholas DeFeo, Victoria Huffman, Rees Cooke, Kaitlyn O'Daniel, David J Riese II 2024. Stable Transfection of Plasmid DNA into Adherent Rodent Cell Lines Using Calcium Phosphate. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzxz5rgx1/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 15, 2024

**Last Modified:** February 15, 2024

**Protocol Integer ID:** 95298

**Keywords:** stable transfection of plasmid dna, plasmid dna, dna into rodent fibroblast cell line, puromycin resistance gene, transfecting dna, transfected dna, such as the neomycin resistance gene, hygromycin resistance gene, neomycin resistance gene, antibiotic resistance gene, expression vector for an antibiotic resistance gene, rodent fibroblast cell line, plasmid, stable transfection, adherent rodent cell line, using calcium phosphate, calcium phosphate precipitate, dna, gene, calcium

**Funders Acknowledgements:**

NIH

Grant ID: R15CA280767

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

A traditional strategy for stably transfecting DNA into rodent fibroblast cell lines features calcium phosphate precipitates. Here we describe our laboratory protocol for this strategy. We have assumed that the transfected DNA contains an expression vector for an antibiotic resistance gene, such as the neomycin resistance gene, the hygromycin resistance gene, the puromycin resistance gene, and the like.

## Troubleshooting

## Introduction

- 1 A traditional strategy for stably transfecting DNA into rodent fibroblast cell lines features calcium phosphate precipitates. Here we describe our laboratory protocol for this strategy. We have assumed that the transfected DNA contains an expression vector for an antibiotic resistance gene, such as the neomycin resistance gene, the hygromycin resistance gene, the puromycin resistance gene, and the like.

## Methods

### 2 ***Cell Lines and Cell Culture***

- 2.1 NIH/3T3 [1] and C127 [2] adherent rodent cell lines are a generous gift from Daniel DiMaio. PA317 [3], Psi-2 [4], and Psi-Cre [5] recombinant retrovirus packaging cell lines are derived from NIH/3T3 cells and are a generous gift from Daniel DiMaio [6].
- 2.2 The cell line to be transfected is typically maintained in T-75 (nominally 75 cm<sup>2</sup>) cell culture flasks according to published or vendor recommendations [1-5, 7]. The cell line is transfected no less than two passages and no greater than five passages after being established from a frozen vial of archived cells.
- 2.3 Cells to be transfected should be seeded at three densities to account for experiment-to-experiment variations in cell growth rate and plating efficiency. Thus, for each transfection,  $1 \times 10^5$ ,  $2 \times 10^5$ , and  $3 \times 10^5$  cells are seeded into a 60 mm (diameter) cell culture dish. In other words, three dishes of cells, each at a different density, are prepared for each transfection. Prepare enough plates of cells to allow for mock and control transfections. Cells are incubated under standard conditions for 24-48 hours to allow cell recovery. Cells should be at 50% confluence at the time of transfection.

### 3 ***Preparation for DNA Transfection***

- 3.1 Two hours prior to transfection, the set of plates that are closest to 50% confluence are selected, while the other sets of plates are discarded. The medium should be changed on the selected set of plates using 5 mL of standard (pre-warmed to 37°C) culture medium and the cells should be returned to the 37°C incubator.
- 3.2 During this two-hour incubation, a sterile 5 mL snap-cap tube (Falcon 352063; VWR 60819-728 [8]) should be labeled for each transfection, including mock and control transfections.

- 3.3 In a biosafety cabinet, 200  $\mu$ L of 2x HEPES-buffered saline (2x HEBS – see recipe below) should be added to each 5 mL snap-cap tube.
- 3.4 In a biosafety cabinet, add 10  $\mu$ g of uncut plasmid DNA to the appropriate 5 mL snap-cap tube(s).
- 3.5 Using sterile technique, make a fresh solution of 250 mM  $\text{CaCl}_2$  by diluting sterile 2 M  $\text{CaCl}_2$  with deionized water. Each transfection requires 200  $\mu$ L of 250 mM  $\text{CaCl}_2$ . The 250 mM  $\text{CaCl}_2$  solution should be formulated in a sterile 1.5 mL microcentrifuge tube or a sterile 50 mL centrifuge tube to facilitate sterile pipetting.
- 3.6 Prepare the calcium phosphate precipitates of DNA in a biosafety cabinet. Approximately 90 minutes after changing the medium on the cells to be transfected, use a micropipetter to add 200  $\mu$ L of 250 mM  $\text{CaCl}_2$  to each HEBS/DNA solution while simultaneously blowing air bubbles into the HEBS/DNA solution using a pipetter and a 1 mL sterile serological pipette. This facilitates formation of calcium phosphate precipitates of DNA.
- 3.7 Vortex to mix the precipitates and incubate for 20 minutes at room temperature in the biosafety cabinet.

#### 4 ***DNA Transfection and Shock***

- 4.1 Drip each precipitate onto a 60 mm dish of cells with swirling. Confirm the presence of precipitates by microscopic inspection. Incubate 4-6 hours at 37°C.
- 4.2 Shock cells as follows:
  - Prepare the glycerol/HEBS shock solution using a 1:1 mixture of 30% glycerol and 2x HEBS
  - Prepare complete (pre-warmed to 37°C) culture medium containing 5 mM NaButyrate.
  - Aspirate medium from each plate.
  - Add 0.7 mL glycerol/HEBS solution and incubate 45 seconds.
  - Aspirate glycerol/HEBS solution and wash twice with PBS.
  - Add 4 mL complete culture medium supplemented with 5 mM NaButyrate.
  - Incubate 12-16 hours at 37°C.

#### 5 ***Passage Cells and Select for Transfected Cells***



- 5.1 Use standard cell culture techniques to subculture each 60 mm dish of transfected cells into three 100 mm (diameter) cell culture dishes.
- 5.2 Incubate 16-24 hours at 37°C.
- 5.3 Select for transfected cells using the appropriate antibiotic.
- 5.4 Depending on the antibiotic, colonies of stably transfected (antibiotic-resistant) cells typically become apparent following 7-12 days of selection.

## Buffer Recipes

### 6 **2x HEBS (250 mL recipe)**

- 6.1 Add the following to a 500 mL Erlenmeyer flask on the benchtop
  - 118 mg  $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$  (62.5 mg anhydrous)
  - 0.5 g Dextrose (D-glucose)
  - 4 g NaCl
  - 185 mg KCl
  - 2.5 g HEPESAdd 230 mL deionized water; stir until the solids are dissolved  
Adjust pH to 7.05 w/10 N NaOH  
Repeat pH next day (use concentrated HCl if pH is too high)  
Adjust volume to 250 mL using deionized water  
Sterile filter using a 0.2  $\mu\text{m}$  filter and aliquot 10 mL into twenty-five 15 mL screw top tubes  
Store at -20°C until use

### 7 **1M NaButyrate ( $\text{NaC}_4\text{H}_7\text{O}_2$ – 250 mL recipe)**

- 7.1 Add the following to a 500 mL Erlenmeyer flask in a fume hood
  - 200 mL deionized water
  - 27.5 g NaButyrateStir until the solid is dissolved  
Adjust volume to 250 mL using deionized water  
Sterile filter using a 0.2  $\mu\text{m}$  filter into a sterile 250 mL media bottle  
Store at 4°C until use



## 8 **30% Glycerol (250 mL recipe)**

- 8.1 Add the following to a 500 mL Erlenmeyer flask in a fume hood
- 150 mL deionized water
  - 75 mL Glycerol
- Stir until the two liquids are thoroughly mixed
- Adjust the volume to 250 mL using deionized water
- Sterile filter using a 0.2 um filter into a sterile 250 mL media bottle
- Store at 4° C until use

## References

- 9
1. Jainchill, J.L., S.A. Aaronson, and G.J. Todaro, *Murine sarcoma and leukemia viruses: assay using clonal lines of contact-inhibited mouse cells*. J Virol, 1969. **4**(5): p. 549-53.
  2. Lowy, D.R., E. Rands, and E.M. Scolnick, *Helper-independent transformation by unintegrated Harvey sarcoma virus DNA*. J Virol, 1978. **26**(2): p. 291-8.
  3. Miller, A.D. and C. Buttimore, *Redesign of retrovirus packaging cell lines to avoid recombination leading to helper virus production*. Mol Cell Biol, 1986. **6**(8): p. 2895-902.
  4. Mann, R., R.C. Mulligan, and D. Baltimore, *Construction of a retrovirus packaging mutant and its use to produce helper-free defective retrovirus*. Cell, 1983. **33**(1): p. 153-9.
  5. Danos, O. and R.C. Mulligan, *Safe and efficient generation of recombinant retroviruses with amphotropic and ecotropic host ranges*. Proc Natl Acad Sci U S A, 1988. **85**(17): p. 6460-4.
  6. Daniel DiMaio - Yale University. [Accessed February 3, 2024]; Available from: <https://medicine.yale.edu/profile/daniel-dimaio/>.
  7. Dwivedi, V. and D.J. Riese 2nd. *Cell Culture Reagents and Sources - 2023-08-30*. protocols.io 2023 [Accessed on February 3, 2024]; Available from: <https://dx.doi.org/10.17504/protocols.io.261ged55wv47/v1>.
  8. *Falcon 352063 5 mL snap-cap tubes - VWR*. [Accessed February 3, 2024]; Available from: <https://pr.vwr.com/store/product/4948460/falcon-round-bottom-high-clarity-polypropylene-tube-disposable-corning>.