

Mar 28, 2025

Standard Uniplex PCR SOP for *Escherichia coli* and *Salmonella*

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp9km1vzp/v1>



Justine Condon¹, Lisa Durso²

¹USDA-ARS; ²USDA - ARS



Justine C. Condon

USDA-ARS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp9km1vzp/v1>

Protocol Citation: Justine Condon, Lisa Durso 2025. Standard Uniplex PCR SOP for *Escherichia coli* and *Salmonella*. protocols.io <https://dx.doi.org/10.17504/protocols.io.dm6gp9km1vzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 18, 2025

Last Modified: March 28, 2025

Protocol Integer ID: 121937

Keywords: Gel Electrophoresis , PCR, Primers, invA, uidA, E. coli, Salmonella, Escherichia coli, Polymerase Chain Reaction, Pathogen Indicators, standard uniplex pcr sop for escherichia coli, standard uniplex pcr sop, gel electrophoresis for the detection, escherichia coli, pcr, salmonella, gel electrophoresis, polymerase chain reaction

Disclaimer

The use of trade, firm, or corporation names in this publication (or page) is for the information and convenience of the reader. Such use does not constitute an official endorsement or approval by the United States Department of Agriculture or the Agricultural Research Service of any product or service to the exclusion of others that may be suitable.

Abstract

To describe the procedures on how to set-up polymerase chain reaction (PCR) and process gel electrophoresis for the detection of *Escherichia coli* (*E. coli*) and *Salmonella*.



Guidelines

Scope

E. coli is a gram-negative rods genus belonging to the Enterobacteriaceae family. They are commonly found in the lower intestinal regions of humans and animals. Most are part of the normal microbiota and can be harmless or beneficial. Some serotypes are pathogenic and can cause severe gastrointestinal illness as a food-borne disease. (WHO, 2018). The *uidA* gene, which encodes the beta-glucuronidase enzyme, was detected in 97.7% of 435 *Escherichia coli* isolates (Martins, 1993) and is used as an indicator of *E. coli*.

Salmonella is a gram-negative rods genus belonging to the Enterobacteriaceae family. *Salmonella* is a ubiquitous and hardy bacterium that can survive several weeks in a dry environment and several months in water. (WHO, 2018). Many *Salmonella* serovars have a broad host range and can infect a wide variety of animals, including mammals, birds, reptiles, amphibians, and insects. In addition, *Salmonella* can grow in plants and can survive in protozoa, soil, and water. (Silva, 2014). They are a group of bacteria that can cause gastrointestinal illness and fever called salmonellosis. It's a common food-borne illness often spread through contaminated food items, or direct contact with infected animals and feces. *invA* is a common molecular target for *Salmonella*-specific detection methods and is recommended by the U.S. Food and Drug Administration Bacteriological Analytical Manual as a target for PCR confirmation of putative *Salmonella* isolates (Buehler, 2019).

Pathogens and antimicrobial resistance are a global public health concern, and tracking indicators assist in surveillance and mitigation efforts. This protocol describes the procedures used to detect the presence of the *uidA* gene from *E. coli* and *invA* gene from *Salmonella* using standard PCR, with detection via gel electrophoresis.

Responsibility

This SOP applies to all staff members and students. These individuals must be knowledgeable about the requirements set forth within this document. The lab manager or designee shall ensure that all staff and students know the proper techniques.

Materials

PCR

- Forward and Reverse Primers (IDT, Coralville, IA) for:
 1. *invA* - **Salmonella** (Sequences in **Table 2**)
 2. *uidA* – **E. coli** (Sequences in **Table 4**)
-  JumpStart™ REDTaq® ReadyMix™ Reaction Mix **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P0982-800RXN**, St. Louis, MO)
-  Molecular Biology Grade Water **HyClone Catalog #SH30538.01**, St. Louis, MO)
- Positive Control – *Salmonella typhimurium* (for *invA*)
- Positive Control – STEC O157 (for *uidA*)
- PCR Prep Cabinet (3620804 Labconco, Kansas City, MO)
- 0.2mL x 8 PCR reaction strips (T320-2N, Simport, Quebec, Canada)
- 2mL sterile microcentrifuge tube (02.681.299 Fisher Scientific, Pittsburg, PA)
- Aerosol-barrier, no-retention pipette tips for various volumes from 1-1000µL
- Pipettors for various volumes from 1-1000µL
- Thermal cycler (T100, BioRad, Hercules, CA)
- Mini Centrifuge (C1301P, Labnet International, Big Flats, NY)

Gel Electrophoresis

1.  Agarose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9539** 500g, St. Louis, MO)
2. 1X Sodium Borate (SB) Buffer (SB20-4 20X SB Buffer FasterBetter Media, Hunt Valley, MD) (Recipe **Appendix A - Step 34**)
3. Microwave
4. Silicone Hot-Pad gloves
5.  SYBR™ Safe DNA Gel Stain **Thermo Fisher Scientific Catalog #S33102**, Carlsbad, CA)
6. PCR amplicons
7.  DirectLoad™ PCR 100 bp Low Ladder **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D3687-1VL**, St. Louis, MO)
8. 500mL or 1-Liter sterile Pyrex glass bottle
9. Gel casting tray, gel combs, and gel boat
10. Electrophoresis machine (PowerEase 300W, Life Technologies, Carlsbad, CA)
11. Micro-Pipettors and appropriate tips (10-20µL)
12. Gel imaging apparatus (Gel Doc XR+, BioRad, Hercules, CA)
13. Lighter

Procedure - *invA* for *Salmonella* PCR

15m 20s

- 1 If frozen, thaw all reagents and samples and keep on an ice block or ice for stability. If kept in the refrigerator, keep there until ready for use, and keep on an ice block or ice for stability when in use.
- 2 Label tubes for the samples and controls.
- 3 Prepare the PCR Prep area for setting up a PCR reaction to eliminate contamination – UV PCR tubes, PCR water, and centrifuge tube for at least  00:15:00 .
- 4 Combine the Master Mix components for the number of reactions (allow for 10% extra Master Mix) using the recipe from **Table 1** with the correct primer sets in **Table 2** for  25 μL reactions.
 - Before adding to the Master Mix - mix the primers by aspirating and expelling with the pipette tip.
 - Vortex the completed Master Mix for at least  00:00:20 . Centrifuge briefly.

15m

20s

A	B	C	D
Component	Description	Volume/RXN (μL)	Final Concentration
TAQ	JumpStart RedTaq ReadyMix Reaction Mix (P0982)	12.5	1X
PCR Water	PCR Water (HyClone SH30538.01)	7.3	N/A
Primer 1	Forward Primer (100 μM)	0.1	2 μM
Primer 2	Reverse Primer (100 μM)	0.1	2 μM
Sample	Template	5	1 μL

Table 1. PCR Master Mix recipe. 25 μL total reaction volume.



	A	B	C	D	E	F	G
	Gene	Type	Sequence 5' to 3'	TM (°C)	TC Conditions	Amplicon Size (bp)	Sequence Reference
	invA	Forward (invA139) Reverse (invA141)	GTGAAATTA TCGCCACGT TCGGGCAA TCATCGCAC CGTCAAAGG AACC	64	1 cycle at 94°C for 10 min; 30 cycles at 94°C for 30 sec, 64°C for 30 sec, 72°C for 30 sec; 1 cycle at 72°C for 5 min.	284	Rahn et al., 1992

Table2. Primer sequences, melting temperature, and thermal cycling conditions.

5

Note

Check the pipette's volume as you go along, making sure it is exactly at the correct volume.

6 Carefully pipette  5 µL of PCR Water to the non-template control (NTC) tube. 

7 Add  20 µL of Master Mix to each labeled PCR tube. 

8 Next, gently vortex to mix then pipette  5 µL of sample ( 5 µL of working stock; or  1 µL straight DNA +  4 µL PCR water) into its same labeled tube. 

Note

- **Use a new pipette tip for each sample**
- To make a working stock dilution, see **SOP: IDT Primer and gBlock Hydration and Aliquots**

9 Add  1 µL of positive control (mix first by aspirating and expelling with the pipette tip) to its same labeled tube. 



10 Centrifuge the strip tubes so that the entire volume is at the bottom.

Note

Make sure the centrifuge is balanced with both sides containing tubes.

11 Load the tubes into the Thermocycler, keeping them in order.

12 Find the programmed gene conditions and double check that they match the conditions listed in **Table 2**.

	A	B	C	D	E	F	G
	Gene	Type	Sequence 5' to 3'	TM (°C)	TC Conditions	Amplicon Size (bp)	Sequence Reference
	invA	Forward (invA139) Reverse (invA141)	GTGAAATTA TCGCCACG TTCGGGCA A TCATCGCAC CGTCAAAG GAACC	64	1 cycle at 94°C for 10 min; 30 cycles at 94°C for 30 sec, 64°C for 30 sec, 72°C for 30 sec; 1 cycle at 72°C for 5 min.	284	Rahn et al., 1992

Table2. Primer sequences, melting temperature, and thermal cycling conditions.

13 Run PCR

13.1 When cycle is complete:

- Gel electrophoresis (See Section **Procedure Gel Electrophoresis** [⇒ go to step #17](#))
OR
- Store the tubes in a labeled rack at -20 °C until ready to run a gel.
- **Turn off the Thermocycler.**

Procedure - *uidA* for *E. coli* PCR

20s

14 Follow **Steps 1–3** [go to step #1](#) to prep for PCR.

15 Combine the Master Mix components for the number of reactions (allow for 10% extra Master Mix) using the recipe from **Table 3** with the correct primer sets in **Table 4** for  25 µL reactions.

20s

- Before adding to the Master Mix - mix the primers by aspirating and expelling with the pipette tip.
- Vortex the completed Master Mix for at least  00:00:20 .

A	B	C	D
Component	Description	Volume/RXN (µL)	Final Concentration
TAQ	JumpStart RedTaq ReadyMix Reaction Mix (P0982)	12.5	1X
PCR Water	PCR Water (HyClone SH30538.01)	7.3	N/A
Primer 1	Forward Primer (100µM)	0.1	2µM
Primer 2	Reverse Primer (100µM)	0.1	2µM
Sample	Template	5	1µL

Table 3. PCR Master Mix recipe. 25µl total reaction volume

A	B	C	D	E	F	G
Gene	Type	Sequence 5' to 3'	TM (°C)	TC Conditions	Amplicon Size (bp)	Sequence Reference
uidA	Forward (uidA241) Reverse (uidA383)	CAGTCTGG ATCGCGAA AACTG ACCAGACG TTGCCAC ATAATT	63	1 cycle at 95°C for 1 min; 40 cycles at 94°C for 10 sec, 63°C for 40 sec, 72°C for 30 sec; 1 cycle at 72°C for 5 min.	142	USDA SOP: MDP-MTH-11

Table 4. Primer sequences, melting temperature, and thermal cycling conditions.



16 Follow **Steps 5 – 13**  [go to step #5](#) using **Table 4** to run the PCR. 

- Since the total reaction volume is  25 μL , pipette  20 μL of Master mix into each tube then  5 μL of sample or  1 μL control.

Procedure - Gel Electrophoresis

2h 11m 30s

17 Using a spoon, measure out  2.6 g of powdered agarose into a weigh boat.

18 Measure  260 mL of **1X SB** buffer into a graduated cylinder.

19 Add powdered agarose and **1X SB** into 1L bottle and swirl to mix. 

20 Microwave **SB** and agarose mixture for  00:01:00 with bottle lid LOOSE. 1m

21 Remove bottle from microwave (**using hot pad gloves**), tighten the lid and gently swirl.

22 LOOSEN the lid and put in microwave for another minute.

23 Check mixture visually, if flecks are present repeat  [go to step #21](#) , loosen lid and microwave again for  00:00:30 . 30s

24 Add  26 μL of SYBR safe to the molten gel. 

25 Swirl gently with a loose lid to distribute the SYBR.

26 Place bottle in 55°C water bath for  00:10:00 to temper. 10m



- 27 While the mixture is tempering, set up gel cast by putting two plastic stoppers into the ends of the cast, with the orange rubber side of the stopper facing out.
- 28 Put two (or three) 36-well combs into the first and fourth notches of the gel cast (gel cast is right side up if first notch is close to the top of the cast).
- 29 After the gel has tempered, slowly pour it (**using hot pad gloves**) into the gel cast.
- 30 If there are any noticeable bubbles, take a lighter and hold it over the bubble, which should make it pop.
- 31 Cover gel loosely with foil and allow it to set (up to an hour) before use.
- 32 Make sure **SB gel** goes into a gel boat with **SB buffer**.

33 Electrophoresis Conditions

33.1  10 μL of PCR amplicon will go into a well, flanked by  5 μL of ladder.

33.2 Run at 120 V for  02:00:00 in 1X SB buffer.

2h

33.3 Remove processed gel and analyze using a gel doc (BioRad Gel Doc).

33.4 Samples with a positive band will be confirmed on qPCR (See **SOP: 5-Gene Uniplex Probe-Based qPCR SOP**).

Appendix A – Recipes

34 **1x SB Buffer Stock – 10 liters**



34.1 **Directions to make 10L:**

Prepare a clean,  10 L carboy.

34.2 In a graduated cylinder, measure out  500 mL of 20x SB, and pour into the carboy.

34.3 In a graduated cylinder, measure out  9500 mL of dH₂O, and pour into the carboy.

34.4 Put the cap on tight, and swirl gently to mix.

34.5 Autoclave the carboy.

- Place autoclave tape on the cap, and write SB on the tape
- Put carboy in a large autoclave bin
- Loosen the cap
- Autoclave at  121 °C
- Label and store at  Room temperature .

Protocol references

1. Buehler , Ariel J., Martin Wiedmann, Zeina Kassaify, Rachel A. Cheng. Evaluation of *invA* Diversity among *Salmonella* Species Suggests Why Some Commercially Available Rapid Detection Kits May Fail To Detect Multiple *Salmonella* Subspecies and Species. *Journal of Food Protection*, Volume 82, Issue 4, 2019, Pages 710-717. ISSN 0362-028X. <https://doi.org/10.4315/0362-028X.JFP-18-525>.
2. Martins MT, Rivera IG, Clark DL, Stewart MH, Wolfe RL, Olson BH. Distribution of *uidA* gene sequences in *Escherichia coli* isolates in water sources and comparison with the expression of beta-glucuronidase activity in 4-methylumbelliferyl-beta-D-glucuronide media. *Appl Environ Microbiol*. 1993 Jul;59(7):2271-6. doi: 10.1128/aem.59.7.2271-2276.1993. PMID: 8357258; PMCID: PMC182268.
3. Rahn, K., S. A. De Grandis, R. C. Clarke, S. A. McEwen, J. E. Gala'n, C. Ginocchio, R. Curtiss III, and C. L. Gyles. 1992. Amplification of an *invA* gene sequence of *Salmonella typhimurium* by polymerase chain reaction as a specific method of detection of *Salmonella*. *Mol. Cell. Probes* 6:271–279.
4. Schwesig, D., Borchers, U., Chancerelle, L., Duffek, A., Eriksson, U., Goksoyr, A., Lamoree, M., Lepom, P., Leonards, P., Leverett, D., McLachlan, M., Poulsen, V., Robinson, R., Silharova, K., Tolgyessy, P., Wegener, JW., Westwood, D. 2009. NORMAN Network of reference laboratories and related organisations for monitoring and bio-monitoring of emerging environmental pollutants. European Commission, Contract no018486, V4.1.
5. Silva, C., Calva, E., & Maloy, S. (2014). One Health and Food-Borne Disease: *Salmonella* Transmission between Humans, Animals, and Plants. *Microbiology Spectrum*, 2(1). <https://doi.org/10.1128/microbiolspec.oh-0020-2013>
6. USDA, MDP-MTH-11. 2011. Real-time PCR Detection of Shiga-toxin Producing *Escherichia coli* (STEC), Serotype O157 and non-O157 in Fresh Produce and Food with non-O157 Isolation and Identification. Real-time PCR Detection of STEC (usda.gov)