

Feb 02, 2024

STC-1 Passaging and PD Associated Microbe Stimulation Protocol

DOI

dx.doi.org/10.17504/protocols.io.261gedmn7v47/v1

Malu G Tansey¹

¹College of Medicine | University of Florida



Senthilkumar Karuppagounder

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gedmn7v47/v1>

Protocol Citation: Malu G Tansey 2024. STC-1 Passaging and PD Associated Microbe Stimulation Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261gedmn7v47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 02, 2024



Last Modified: May 31, 2024

Protocol Integer ID: 94578

Keywords: ASAPCRN, pd associated microbe stimulation protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) Collaborative Research Network

Abstract

STC-1 Passaging and PD Associated Microbe Stimulation Protocol

Materials

Complete Medium

DMEM, high glucose, no glutamine

4mM L-glutamine

1mM sodium pyruvate

5% FBS

Other Reagents

0.05% Trypsin-EDTA

1X D-PBS without calcium and magnesium

75cm² Flask

12 well plates

Troubleshooting



Passaging & Plating Protocol

- 1 Cells are passaged when they reach approximately 70% confluence

Warm trypsin, 1XPBS, and complete medium in water bath for 15 min.
- 2 Aspirate off culture medium
- 3 Rinse cells with 10mL D-PBS
- 4 Add 3mL Trypsin to flask. Place flask back into 37°C/5% CO₂ incubator for 7 minutes.
- 5 Observe cells under microscope to ensure they have detached from the bottom of the flask.

a. These cells are particularly sticky, so giving them a nice tap may be necessary to detach them.
- 6 Add 9mL of complete medium to flask and transfer cell suspension to a 15mL tube.
- 7 Spin down at 130g for 5 min
- 8 Aspirate off supernatant and resuspend cell pellet in 1mL culture medium.
- 9 Count cells
- 10 Prepare cell suspension at a concentration of 3×10^5 cells/mL
- 11 Add 1mL cell suspension per well in 12 well plate(s)
- 12 For PD associated microbe stimulation, allow cells to rest for approximately 24 hours,



- 13 Treat STC-1s with 10 colony forming units (CFU) of PD associated microbe per STC-1 cell
- 14 Incubate for desired time and then collect conditioned media and/or cells for downstream analyses

LPS Pre-Treatment

- 15 Following plating, allow cells to rest for 24 hours
- 16 Treat cells with 10ng/mL LPS for 24 hours
- 17 24 hr after LPS treatment, stimulate cells with 10 CFU of PD associated microbes per STC-1 cell
- 18 Incubate for desired time and then collect conditioned media and/or cells for downstream analyses