

Dec 08, 2023

TFome CRISPRko screens

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbxk8olpk/v1>

Sean R. McCutcheon¹, Adam M. Swartz¹, Michael C. Brown¹, Alejandro Barrera¹, Christian McRoberts Amador¹, Keith Siklenka¹, Lucas Humayun¹, Maria A. ter Weele¹, James M. Isaacs¹, Andrea Daniel¹, Timothy E. Reddy¹, Andrew Allen¹, Smita K. Nair¹, Scott J. Antonia¹, Charles A. Gersbach¹

¹Duke University

Andrea Daniel: This protocol was adapted from work of Sean McCutcheon and colleagues in the Gersbach lab at Duke University.

Gersbach Lab



Andrea R Daniel

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbxk8olpk/v1>



Protocol Citation: Sean R. McCutcheon, Adam M. Swartz, Michael C. Brown, Alejandro Barrera, Christian McRoberts Amador, Keith Siklenka, Lucas Humayun, Maria A. ter Weele, James M. Isaacs, Andrea Daniel, Timothy E. Reddy, Andrew Allen, Smita K. Nair, Scott J. Antonia, Charles A. Gersbach 2023. TFome CRISPRko screens. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbxk8olpk/v1>

Manuscript citation:

McCutcheon, S.R., Swartz, A.M., Brown, M.C. *et al.* Transcriptional and epigenetic regulators of human CD8⁺ T cell function identified through orthogonal CRISPR screens. *Nat Genet* (2023). <https://doi.org/10.1038/s41588-023-01554-0>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 08, 2023

Last Modified: December 08, 2023

Protocol Integer ID: 92033

Keywords: CRISPR, screen, knockout, BATF3, T cells, methods for transcription factor crispr knockout screen, transcription factor crispr knockout screen, cofactors of batf3, tfome crisprko screen, tfome crisprko screens this protocol, batf3, transcription

Funders Acknowledgements:

NIH

Grant ID: HG012053

Abstract

This protocols describes methods for transcription factor CRISPR knockout screens to reveal cofactors of BATF3 in T cells.

gRNA library construction

- 1 The Brunello genome-wide KO⁸³ library (four gRNAs per gene) was subset for 1,612 TFs⁴⁵ and *IL7R*.
- 2 A total of 550 NT gRNAs were included in the library for a total of 7,000 gRNAs (available in Supplementary Table 6 of McCutcheon et al. Nature Genetics, 2023. <https://doi.org/10.1038/s41588-023-01554-0>).
- 3 This gRNA library was cloned into SpCas9 gRNA lentiviral plasmids with either mCherry or BATF3.

gRNA library cloning

- 4 Oligonucleotide pools containing variable gRNA sequences and constant regions for polymerase chain reaction (PCR) amplification were synthesized by Twist Bioscience.
- 5 gRNA amplicons were gel extracted, PCR purified and input into 20 µl Gibson reactions (5:1 molar ratio of insert to backbone) with 200 ng of Esp3I digested and 1 × solid-phase reversible immobilization (SPRI)-selected (Beckman Coulter) plasmid backbone.
- 6 Gibson reactions were purified using ethanol precipitation and transformed into Lucigen's Endura ElectroCompetent Cells.
- 7 Transformed cells were cultured overnight and plasmids were isolated using Qiagen Midi Kits.

TFome CRISPRko screens

- 8 A total of 20×10^6 CD8⁺ T cells from two donors were activated with CD3/CD28 dynabeads at a 1:1 ratio.
- 9 At 24 h post-activation, CD8⁺ T cells were split evenly and transduced in parallel with TFome CRISPRko gRNA libraries with mCherry or BATF3.
- 10 At 48 h post-activation, cells were electroporated with Cas9 protein. Briefly, the cells were collected, spun down at 90g for 10 min, resuspended in 100 µl of Lonza P3 Primary Cell buffer with 3.2 µg Cas9 (Thermo) per 10^6 cells, and electroporated with the pulse code EH115.

- 11 After electroporation, warm medium was immediately added to each cuvette and cells were recovered at 37 °C for 20 min before being transferred into a six-well plate.
- 12 On day 3 post transduction, cells were selected with 2 µg ml⁻¹ of puromycin for 3 days. On day 9 post transduction, cells were stained for CD8, IL7R and a viability dye.

12.1

	Antibody Target	Fluorophore/Sequence	Clone	Isotype	Dilution	Application	Manufacturer	Catalog #
	IL7RA	PE-Cy5	eBioR DR5	Mouse / IgG1, kappa	1:100	Flow cytometry	Thermo	15-1278-42
	CD8	bv-421	HIT8a	Mouse IgG1, κ	1:50	Flow cytometry	BD Biosciences	740078

- 13 An SH800 FACS Cell Sorter (Sony Biotechnology) was used for cell sorting and analysis.
- 14 Cells were then washed with flow buffer, spun down at 300g for 5 min and resuspended in flow buffer for cell sorting.
- 15 Viable CD8⁺ T cells in the lower and upper 10% tails of IL7R expression were sorted for subsequent gRNA library construction and sequencing.
- 16 All replicates were maintained and sorted at a minimum of 75× coverage.

Protocol references

45. Lambert, S. A. et al. The human transcription factors. *Cell* **172**, 650–665 (2018).
83. Doench, J. G. et al. Optimized sgRNA design to maximize activity and minimize off-target effects of CRISPR-Cas9. *Nat. Biotechnol.* **34**, 184–191 (2016).