


Jun 01, 2023

🌐 Transient CRISPR-Cas9 Coupled with Electroporation Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l64x2kvqe/v1>



Amy Gladfelter¹

¹Duke University



Amy Gladfelter

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l64x2kvqe/v1>

Protocol Citation: Amy Gladfelter 2023. Transient CRISPR-Cas9 Coupled with Electroporation Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l64x2kvqe/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 02, 2018

Last Modified: June 01, 2023

Protocol Integer ID: 17448

Keywords: electroporation protocol transient crispr, cas9 transformation of cryptococcus neoforman, transient crispr, cryptococcus neoforman, cas9 transformation, cas9

Abstract

Transient CRISPR-Cas9 transformation of *Cryptococcus neoformans*.

Attachments



CRISPR Electroporati...

70KB

Guidelines

References

1. Lin, X., et al., Generation of stable mutants and targeted gene deletion strains in *Cryptococcus neoformans* through electroporation. *Med Mycol*, 2015. **53**(3): p. 225-234.
2. Fan, Y. and X. Lin, Multiple applications of a transient CRISPR-Cas9 coupled with electroporation (TRACE) system in the *Cryptococcus neoformans* species complex. *Genetics*, 2018.



1 PCR amplification of CAS9, sgRNA, your construct.

- CAS9: Use plasmid pXL1-CAS9-HYG as template with primers CAS9-F and CAS9-R. (6985 bp)
- sgRNA: U6P and sgRNA scaffold
 1. U6pPromoter is PCR amplified using serotype D genomic DNA as template with primers U6P-F and GOI-sgRNA-R. ~295 bp
 2. SgRNA scaffold is PCR amplified using pYF515 as template with primers GOI- sgRNA-F and sgRNA-R. ~108 bp
 3. sgRNA construct is PCR amplified using above two PCR products as template with primers U6P-F and sgRNA-R. ~383 bp

2 Mix 2 µg your construct DNA, 100 ng sgRNA, and 170 ng CAS9 in an Eppendorf tube.

3 Vacuum dry the DNA and elute in 5 µL DNase/RNase free water.

Note

Notes: use the combination of 2 µg construct DNA, 1 µg CAS9 DNA, and 700 ng sgRNA can increase the transformation efficiency, but low dose is sufficient to obtain transformants.

4 Inoculate recipient strain in 5 mL YPD liquid medium, culture overnight at 30 °C with shaking at 250 rpm.



- 5 Use the overnight culture to inoculate  100 mL fresh YPD medium at an initial inoculum of OD600=0.2. Grow the cells for additional  04:00:00 to  05:00:00 until the cell density reached OD600 between 0.6-1.0.
From this step on, everything on ice and centrifugation at  4 °C
- 6 Collect cells by centrifugation at 3200g for  00:05:00 at  4 °C .
- 7 Wash cells with ice-cold water (EB Buffer instead of water in 2015 paper). (wash 1/2)
- 8 Wash cells with ice-cold water (EB Buffer instead of water in 2015 paper). (wash 2/2)
- 9 Suspend cells in  10 mL ice-cold EB buffer (10 mM Tris-HCl, pH 7.5, 1mM MgCl₂, 270 mM Sucrose) with 1mM DTT.
- 10 Incubate the cells on ice for an hour ( 01:00:00).

**Note**

(30 to 60 mins in 2015 paper)

- 11 (Optional) Wash cells with  10 mL ice-cold EB buffer once (2015 paper).
- 12 Collect the cells by centrifugation and resuspend in  250 μ L EB buffer.
- 13 Mix  45 μ L cells with  5 μ L DNA mix from step 2 in a pre-cooled 2 mm gap electroporation cuvette.
- 14 Transform the DNA by electroporation using the BioRad gene pulser with settings of 0.45 kV, 125 μ F, 600 Ω . (If using an Eppendorf multiporator, use the bacterial mode with V=2 kV with τ optimized for 5 ms)
- 15 Suspend electroporated cells in  1 mL of YPD medium and culture at  30 $^{\circ}$ C for  02:00:00 ( 01:30:00 in 2015 paper) before plating onto the appropriate selective agar medium.

