

May 13, 2023

Version 2

TS Procure 812 - primary fixation Karnovsky's - transwells (TM - 013) V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwjwpwlmk/v2>

sandra.crameri¹

¹CSIRO ACDP



sandra.crameri

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwjwpwlmk/v2>

Protocol Citation: sandra.crameri 2023. TS Procure 812 - primary fixation Karnovsky's - transwells (TM - 013). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwjwpwlmk/v2> Version created by [sandra.crameri](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 13, 2023



Last Modified: May 13, 2023

Protocol Integer ID: 81816

Keywords: Procure 812, formalin, epon, procure 812 resin block, transwells to procure, fixed transwell, resin block, ts procure, transwell, resin

Abstract

The method is used to Karnovsky's EM fixed transwells to Procure 812 resin blocks. The presence of 4% paraformaldehyde is only used to satisfy microsecurity requirements.

Guidelines

Sample is not optimal for EM, best to dissect and use only the outer 2mm of tissue mass. Sample is usually derived from histology.

Protocol materials

⊗ 25% Glutaraldehyde

⊗ 16% Paraformaldehyde ProSciTech Catalog #C004

⊗ Osmium Tetroxide ProSciTech

⊗ Ethanol

⊗ Propylene oxide ProSciTech Catalog #C156

⊗ Propylene oxide ProSciTech Catalog #C156

⊗ Propylene oxide ProSciTech Catalog #C156



HEADER

1 SAN:

SPEC No:

OPERATOR & STEPS:

OPERATOR & STEPS:

Conventional

5h 15m

2 All solutions are exchanged above and below the membrane, up to step 10.

3 Karnovsky's fixative: [M] 2.5 % volume 25% Glutaraldehyde and
[M] 4 % volume 16% Paraformaldehyde **ProSciTech Catalog #C004** in
[M] 0.1 Molarity (M) Sorenson's Phosphate buffer (pH 7.2, 300mosmol/kg) for at least
 00:40:00 or overnight

40m

4 [M] 0.1 Molarity (M) buffer 00:15:00

15m

5 [M] 1 % volume Osmium Tetroxide **ProSciTech** in buffer 01:00:00

1h

6 [M] 70 % volume Ethanol for at least 00:20:00 or overnight

20m



- 7 [M] 95 % volume Ethanol 00:20:00 20m
- 8 [M] 100 % volume Ethanol 00:20:00 20m
- 9 [M] 100 % volume Ethanol dried with molecular sieve 00:20:00 20m

Remove membrane

- 10 Remove membrane from transwell scaffold using a 45o scalpel and transfer to a glass vial with [M] 100 % volume Ethanol dried with molecular sieve.

Conventional

5h 15m

- 11 [M] 50 % volume Ethanol [M] 50 % volume Propylene oxide ProSciTech Catalog #C156 00:20:00 20m
- 12 [M] 100 % volume Propylene oxide ProSciTech Catalog #C156 00:20:00 20m
- 13 [M] 50 % volume Propylene oxide ProSciTech Catalog #C156 [M] 50 % volume Procure 812
- 14 [M] 100 % volume Procure 812 00:20:00 20m
- 15 [M] 100 % volume Procure 812 01:00:00 1h
- 16 Embed & polymerise overnight in oven at 60 °C , section with diamond knife as usual.

