

Oct 08, 2025

Use of cholera toxin subunit B to label neural projections to major pelvic ganglia

 Forked from [Use of cholera toxin subunit B to label neural projections to lower urinary tract organs](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgq71jovk5/v1>

Janet R Keast¹, Peregrine B Osborne¹, Nicole Wiedmann¹, John-Paul Fuller-Jackson¹

¹University of Melbourne


SPARC

Tech. support email: info@neuinfo.org



Janet R Keast

University of Melbourne



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgq71jovk5/v1>

Protocol Citation: Janet R Keast, Peregrine B Osborne, Nicole Wiedmann, John-Paul Fuller-Jackson 2025. Use of cholera toxin subunit B to label neural projections to major pelvic ganglia. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.36wgg71jovk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 24, 2022

Last Modified: October 08, 2025

Protocol Integer ID: 65218

Keywords: lower urinary tract, retrograde tracing, neural tracing, major pelvic ganglia, spinal cord, neural projections to major pelvic ganglia, major pelvic ganglia, innervating pelvic visceral organ, preganglionic neurons in the spinal cord, pelvic visceral organ, preganglionic neuron, lower urinary tract, spinal cord, neuron, neural projection, methods of anesthesia, anesthesia

Funders Acknowledgements:

NIH SPARC

Grant ID: OT2OD023872

Abstract

This protocol is used to visualise preganglionic neurons in the spinal cord innervating pelvic visceral organs (e.g., the lower urinary tract) in an experimental adult male or female rat. The protocol is performed under anesthesia and should incorporate all local requirements for standards of animal experimentation, including methods of anesthesia, surgical environment, and post-operative monitoring and care.

Materials

MATERIALS

-  Parafilm
-  Isoflurane **Zoetis Catalog #10015516**
-  Lacrilube **Ellar Laboratories**
-  Cholera toxin subunit B **List Labs Catalog #104**
-  Glass capillaries **Warner instruments Catalog ##GC100F-15**
-  Sterile Saline (0.9% NaCl)
-  Evans Blue Dye **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E2129**

Equipment

Picospritzer III Intracellular Microinjection Dispense System	NAME
Injection system	TYPE
Picospritzer	BRAND
052-0500-900	SKU
https://ph.parker.com/us/12051/en/picospritzer-iii-intracellular-microinjection-dispense-systems-picospritzer-micro-dispense-system/052-0500-900	LIN K
100 psi, 2 channel	SPECIFICATIONS

Preparation for surgery

- 1 Prepare cholera toxin subunit B solutions: low salt formulation with 0.05% Evans Blue.
- 2 Prepare glass pipettes for surgery by prefilling each pulled glass pipette with cholera toxin subunit B
- 3 Anesthetise animal (2.5% isoflurane in oxygen, or as required for maintenance)
- 4 Apply eye lubricant and place animal on heated pad.
- 5 Shave and clean the ventral abdomen.

Surgery

- 6 Perform a midline incision in the skin and then the muscle, then gently move organs to visualise the required injection site.

The major pelvic ganglia is located on the dorsal lateral lobe of prostate (male) or cervix (female)
- 7 Using fine angled forceps, gently blunt dissect underneath the pelvic ganglia. Once the ganglia is separated from the underlying tissue, slide a sterile 2mm x 2mm piece of parafilm between it and the tissue.
- 8 Microinject sterile tracer solution at the selected injection site using a glass pipette attached to a picospritzer. At each injection site, hold the glass pipette in place for ~5 seconds after ejection of the dye, to enable the dye to spread to the underlying tissue. This also minimises leakage. Continue injections until the desired volume is reached.
- 9 Wash all injection sites with sterile saline.
- 10 For a bilateral injection, perform steps 6-8 on the alternate side.



- 11 Close the muscle and skin using approved procedures. Administer analgesics and monitor animal during postoperative period as per local approved procedure.

Tissue harvesting

- 12 To analyse tracer distribution in spinal cord or the injection site, 3 days after surgery, anaesthetise animals as per local ethical requirements, and perform intra-cardiac perfusion with fixative, then dissect tissues of interest for further study.