

Feb 11, 2026

VarChamp Cell Painting

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqnwxgk5/v1>

Tanisha Teelucksingh¹, Florent Laval²

¹University of Toronto; ²Center for Cancer Systems Biology (CCSB), Dana-Farber Cancer Institute, Boston, MA



Tanisha Teelucksingh

University of Toronto

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqnwxgk5/v1>

Protocol Citation: Tanisha Teelucksingh, Florent Laval 2026. VarChamp Cell Painting . protocols.io

<https://dx.doi.org/10.17504/protocols.io.36wgqnwxgk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: October 07, 2024

Last Modified: February 11, 2026

Protocol Integer ID: 109318

Keywords: varchamp cell painting, missense variants on protein, subcellular localization, missense variant

Funders Acknowledgements:

National Human Genome Research Institute (NHGRI)

Grant ID: UM1HG011989

Abstract

Assessment of the impact of missense variants on protein subcellular localization



Virus Production - Seeding transfection plates (96-well format)

- 1 Seed HEK293T Cells - 96-well format,  100 μL per well. Media used: DMEM

4h 30m

Note

All media contains 10% FBS and 5% penicillin streptomycin

1. if transfecting packaging plasmids on the same day, seed 2.4×10^4 cells/well. Allow at least  04:00:00 for cells to adhere, incubating at  37 °C , 5% CO₂
2. if transfecting packaging plasmids the next day, seed 1.2×10^4 cells/well. Incubate  Overnight at  37 °C , 5% CO₂

- 2 Remember to harvest mycoplasma sample for testing downstream



Virus Production - Preparing packaging/transfection mix (96-well format)

- 3 Prepare packaging mix -  10 μL /well
1. 50ng/well psPAX2
 2. 5ng/well pVSV-G
 3. Top up to  10 μL /well with OptiMEM

- 4 Prepare transfection mix -  10 μL /well
1.  9.4 μL OptiMEM/well
 2.  0.6 μL TransIT/well

Note

Make sure to add the TransIT drop wise and mix well

- 5 Prepare 96-well plate - 33.33 μL total volume per well
1. 10 μL packaging mix
 2. 10 μL transfection mix
 3. 10 μL sterile water
 4. 3.33 μL lentiviral expression plasmid (normalized to 15ng/ μL)

30m





5. Incubate at room temperature for 00:30:00

Note

Mix well before incubation

6 Transfer entire mix (33.33uL) to HEK293T cells

1. Incubate for 48:00:00 at 37°C, 5% CO₂

2d

Transduction of U2OS cells

7 Seed U2OS Cells - 384-well format, 40uL per well. Media used: McCoy's 5A with 10% FBS and 1% penicillin streptomycin

1. Seed 1000 cells/well → include 10ug/mL polybrene

2. Allow cells to settle for ~ 1 hour at room temperature before incubating at 37°C, 5% CO₂ for an additional 3 hours

8 Remember to harvest mycoplasma sample for testing downstream



9 Spinfection - 384-well format

1. Add 6uL of virus per well

2. Spin plates at 1200 x g, 37°C for 30 minutes

3. Incubate for 48 hours at 37°C, 5% CO₂



Puromycin Selection of Transduced U2OS cells

10 Dilute puromycin in media (McCoy's 5A with 10% FBS and 1% penicillin streptomycin) to a final concentration of 1.5ug/mL

1. Flick media from 384-well plates

2. Add 30uL of diluted puromycin to each well

3. Incubate for 48:00:00 at 37°C, 5% CO₂

4. Make sure to monitor untransduced control for complete selection

2d



Mycoplasma Testing

11 Before proceeding with staining and fixation of cells, samples taken on steps 2 and 8 should be tested for mycoplasma as this will impact downstream applications





Mitotracker Staining

- 12 Dilute MitoTracker Deep Red in McCoy's 5A (with 10% FBS and 1% penicillin streptomycin) to a final concentration of 500nM
1. Flick media from 384-well plates
 2. Add 20uL of diluted MitoTracker Deep Red to each well
 3. Incubate for 1 hour at 37°C, 5% CO₂

Cell Fixation

- 13 Dilute 32% paraformaldehyde to a final concentration of 4% in McCoy's 5A (with 10% FBS and 1% penicillin streptomycin) 
1. Flick media from 384-well plates
 2. Add 20uL of diluted PFA to each well
 3. Cover with tin foil - PFA is light sensitive
 4. Incubate for 15-30 minutes at room temperature
 5. Wash each plate 4 times with 1x HBSS
 6. If not staining right away, store plates with 50uL 1x HBSS per well, covered in tin foil at 4°C

Cell Staining

- 14 Prepare Staining and Permeabilization Buffer in 1x HBSS:
1. 0.1% Triton (vol/vol)
 2. 1% BSA (wt/vol)
- 15 Dilute dyes to the following concentrations:
1. Phalloidin 568 - 0.125uL/mL (stock 6.6uM)
 2. Hoechst 33342 - 1ug/mL (stock 10mg/mL)
 3. WGA 555 - 1.5ug/mL (stock 1mg/mL)

Note

All dyes are spun down before use according to Cell Painting Protocol V3

- 16 Flick out 1x HBSS from cells after fixation
1. Add 20uL of staining solution per well
 2. Centrifuge at 500 x g for 1 min to remove air bubbles
 3. Incubate, covered in tin foil for 30 minutes
 4. Wash each plate 4 times with 1x HBSS
 5. Leave 50uL 1x HBSS per well

Imaging

17 All confocal images are captured on a Perkin Elmer Opera Phenix Microscope with the following conditions:

1. 20X water objective
2. Single plane (-7um)
3. 384 wells, 9 fields
4. Separate all channels:

Channel	Excitation	Emission	Exposure	Power
DAPI	405	435-480	100 ms	60%
Alexa488	488	500-550	60 ms	30%
Alexa568	561	571-596	60 ms	30%
Mitotracker Deep Red	640	650-760	20 ms	20%
Brightfield High	Transmissio n		100 ms	50%
Brightfield	Transmissio n		100 ms	50%
Brightfield Low	Transmissio n		100 ms	50%

Note

Current run time is 2hr 12min for 1 plate