

Mar 25, 2025

Variant painting via fluorescence

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8erydl2w/v1>

Chloe Reno¹

¹University of Toronto

IGVF



Chloe Reno

University of Toronto

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8erydl2w/v1>

Protocol Citation: Chloe Reno 2025. Variant painting via fluorescence. protocols.io
<https://dx.doi.org/10.17504/protocols.io.5jyl8erydl2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 25, 2025

Last Modified: March 25, 2025

Protocol Integer ID: 124928

Keywords: variant painting via immunofluorescence, cell painting protocol v3, variant painting, variation on protein localization, immunofluorescence, protein localization, protein, imaging, fluorescence, throughput imaging platform, cell, coding variation, additional cellular compartment, mitochondria

Funders Acknowledgements:

NIH IGVF

Grant ID: HG011989

Abstract

High-throughput imaging platform to assay the impact of coding variation on protein localization. Proteins are fused to mNeonGreen for visualization, alongside additional cellular compartments including the nucleus, mitochondria, actin, golgi, and plasma membrane. This protocol is adapted from Variant Painting via immunofluorescence and Cell Painting Protocol v3.



Day 1: Cell Seeding

- 1 Seed 1.2×10^4 HEK293T cells/well in 96-well plate, 100 μL per well. Media is high glucose DMEM supplemented with 10% FBS, 1% pepenicillin-streptomycin. Incubate Overnight at 37 $^{\circ}\text{C}$, 5% CO_2 .

Note

Harvest a sample of HEK293T cells for mycoplasma testing.

Note

If performing same-day transfection, double the seeding density and allow cells to adhere for ~4 hrs.

Day 2: Viral Packaging

2d

- 2 Prepare packaging master mix:
 1. 50 ng /well psPAX2
 2. 5 ng /well pVSV-G
 3. Top up to 10 μL /well with OptiMEM
- 3 Prepare transfection reagent master mix:
 1. 9.4 μL /well OptiMEM
 2. 0.6 μL /well TransIT-LT1 Transfection Reagent
- 4 Prepare 96-well transfection complexes:
 1. Add 10 μL /well molecular grade water
 2. Add 3.33 μL /well lentiviral expression plasmid (normalized to 15 ng/ μL)
 3. Add 10 μL /well packaging master mix
 4. Add 10 μL /well reagent master mixIncubate at room temperature 00:30:00 .



- 5 Transfer transfection complexes to seeded HEK293T cells. Incubate  48:00:00 at  37 °C , 5% CO₂.

2d

Note

Make sure to monitor EGFP transfection control for transfection efficiency.

Day 4: Transduction

2d 4h 30m

- 6 Seed 1000 U2OS cells/well,  40 µL /well in 384-well black optically clear flat-bottom PhenoPlates. McCoy's 5A media is supplemented with 10 µg/mL polybrene + 10% FBS + 1% penicillin-streptomycin. Allow cells to adhere for  01:00:00 at room temperature before incubating at  37 °C , 5% CO₂ for  03:00:00 .

4h

Note

Harvest a sample of U2OS cells for mycoplasma testing. Test before Day 8.

- 7 Perform Spinfection:
1. Add  6 µL virus /well
 2. Spin plates at  1200 x g, 37°C, 00:30:00
 3. Incubate  48:00:00 at  37 °C , 5% CO₂.

2d 0h 30m

Day 6: Puromycin Selection

2d

- 8 Prepare 1.5 µg/mL puromycin in McCoy's 5A media supplemented with 10% FBS + 1% penicillin-streptomycin.
1. Flick media out of 384-well imaging plates into a bleach bucket.
 2. Add  30 µL puromycin media /well
 3. Incubate  48:00:00 at  37 °C , 5% CO₂.

2d

**Note**

Make sure to monitor untransduced (EGFP transfection) control for complete selection

Day 8: Fixation & Staining**2h 16m**

- 9 Stain live cells with MitoTracker DeepRed:
1. Dilute MitoTracker DeepRed to **1mM 500 nanomolar (nM)** in McCoy's 5A + 10% FBS + 1% penicillin-streptomycin.
 2. Flick off media from imaging plates.
 3. Add **20 µL** /well of prepared MitoTracker media.
 4. Incubate **01:00:00** at **37 °C**, 5% CO₂.

1h

- 10 Fix cells:
1. Dilute 32% Paraformaldehyde to 4% in McCoy's 5A + 10% FBS + 1% penicillin-streptomycin.
 2. Flick off media from imaging plates.
 3. Add **20 µL** /well of diluted PFA.
 4. Wrap plates in foil to protect from light.
 5. Incubate **00:15:00** - **00:30:00** at room temperature.
 6. Wash with plate washer 4x with 1xHBSS, including final aspiration.

45m**Note**

If you are not staining & imaging the same day, do not include the final aspiration. Instead leave ~ **50 µL** of HBSS behind and store plates at **4 °C** wrapped in foil.

- 11 Prepare your blocking & permeabilization buffer:
0.1% Triton-X(vol/vol) + 1% BSA (wt/vol) in 1xHBSS.
- 12 Prepare your staining solution. Dilute the following dyes in your blocking & permeabilization buffer:
1. Spin down all dyes
 2. Dilute Phalloidin AlexaFluor 568 to 8.25 nM

- 13 Stain cells:
1. Flick off media from imaging plates.
 2. Add  20 μ L /well of staining solution.
 3. Centrifuge  500 x g, 00:01:00 to remove bubbles.
 5. Wrap plates in foil to protect from light. Incubate  00:30:00 at room temperature.
 6. Wash with plate washer 4x with 1xHBSS, excluding final aspiration t leave ~  50 μ L /well.

- 14 Image:
Confocal images are captured on a PerkinElmer Opera Phenix microscope with the following conditions:

1. 20x water objective
2. Single plane
3. 384-wells, 9 fields
4. Separate Channels:

	Channel	Excitation	Emission	Exposure	Power
	DNA	405	435-480	100 ms	60%
	Protein	488	500-550	60 ms	30%
	AGP	561	571-596	60 ms	30%
	Mitochondria	640	760	20 ms	20%