

Mar 24, 2020 Version 1

Viral Sequencing, from Gunk to Graph V.1

DOI

dx.doi.org/10.17504/protocols.io.bd3yi8pw

David A Eccles¹


¹Malaghan Institute of Medical Research (NZ)

Coronavirus Method De...



David A Eccles

Malaghan Institute of Medical Research (NZ)



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bd3yi8pw

Protocol Citation: David A Eccles 2020. Viral Sequencing, from Gunk to Graph. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bd3yi8pw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 23, 2020

Last Modified: March 24, 2020

Protocol Integer ID: 34648

Keywords: SARS-CoV-2, COVID-19, nanopore, sequencing,



Abstract

This is a fast "gunk to graph" protocol for analysing viral RNA from nasopharyngeal swabs. The approach involves swab lysis and inactivation at the point of sampling, uses a cellulose binding / wash protocol to reduce extraction cost, incorporates sample-specific barcodes during first-strand synthesis, nanopore rapid-attachment primers during PCR amplification, and nanopore sequencing with parallel RAMPART analysis for fast assembly and phylogenetics.

Materials

MATERIALS

⊗ Q5 Hot Start High-Fidelity 2X Master Mix - 500 rxns **New England Biolabs Catalog #M0494L**

⊗ MinION sequencer **Oxford Nanopore Technologies**

⊗ ONT MinION Flow Cell R9.4.1 **Oxford Nanopore Technologies Catalog #FLO-MIN106D**

Additional materials TBA.







Safety warnings


⚠ This protocol is UNTESTED, and is in the early stages of development. Do not trust the protocol; question everything.








Assume samples are potentially infectious during extraction, and make sure to use proper sterile technique to avoid cross-contamination.









Swab Lysis

- 1 Prepare a  1.5 mL centrifuge tube with heated lysis buffer and a cellulose disc
- 1.1 Add  500 μ L lysis / RNase inactivation buffer ([Twitter reference](#)) to 1.5ml centrifuge tube:
 -  10 millimolar (mM) Tris
 -  10 millimolar (mM) EDTA
 -  0.5 % volume SDS
 -  150 millimolar (mM) NaCl



OR  500 μ L extraction buffer #2 (see [paper](#)):






 -  800 millimolar (mM) guanidine hydrochloride
 -  50 millimolar (mM) Tris [pH 8]
 -  0.5 % volume Triton X100
 -  1 % volume Tween-20
- 1.2 Add a  3 mm diameter punched disc from [Whatman #1 filter paper](#) (see [paper](#))
- 1.3 Preheat  1.5 mL tube to  60 °C
- 2 Collect sample using a sterile polystyrene swab with a 30mm breakpoint (e.g. [Puritan 25-3606-U; PurFlock Ultra 6" Sterile Elongated Flock Swab w/Polystyrene Handle, 30mm Breakpoint](#)).

RNA Wash

- 3 Transfer disc to a new  1.5 mL tube containing  200 μ L wash buffer using a pipette tip to remove contaminants:
 -  10 millimolar (mM) Tris [pH 8.0]
 -  0.1 % volume Tween-20
- 4 Incubate tube at  Room temperature for  00:01:00




cDNA Synthesis


- 5 Add the following additional components into the  200 μL PCR tube (see the [Nanopore protocol for Sequence-specific cDNA-PCR Sequencing \(SQK-PCS109\)](#)) in a  11 μL reaction:






-  1 μL x  2 micromolar (μM) reverse primers
-  1 μL x  10 millimolar (mM) dNTPs
-  9 μL RNase-free water

Reverse primers should be prefixed with sample-specific barcode sequences (if used) and the ONT reverse anchor sequence, i.e. [5' - ACTTGCCTGTCGCTCTATCTTC - [barcode] - [sequence-specific] - 3']

- 6 Mix gently **by flicking the tube** and spin down  00:00:05

- 7 Denature RNA and anneal primers at  65 $^{\circ}\text{C}$ for  00:05:00 and then snap cool on a pre-chilled freezer block for  00:01:00







- 8 In a separate tube, mix together the following in an  8 μL reaction:

-  4 μL 5X RT Buffer
-  1 μL RNaseOUT
-  1 μL Nuclease-free water
-  2 μL x  10 micromolar (μM) Strand-switching primer (SSP)

Note: It might be possible to instead carry out only the first-strand synthesis (i.e. excluding SSP), then use a forward primer, tailed with sample-specific barcode sequences (if used) and the ONT forward anchor sequence, i.e. [5' - TTTCTGTTGGTGCTGATATTGC - [barcode] - [sequence-specific] - 3']. One-Step RT-PCR sequencing kits may help with this (e.g. [OneTaq One-Step RT-PCR Kit](#)). For more



details about the reverse anchor sequence and four-primer amplicon sequencing, see the Nanopore protocol for [Four-primer PCR \(SQK-PSK004 or SQK-PBK004\)](#).








- 9 Mix gently **by flicking the tube** and spin down  00:00:05
- 10 Add the strand-switching buffer to the snap-cooled, annealed RNA, mix by **flicking the tube** and spin down
- 11 Incubate at  42 °C for  00:02:00
- 12 Add  1 µL of Maxima H Minus Reverse Transcriptase, to a total volume of  20 µL
- 13 Mix gently by **flicking the tube** and spin down  00:00:05
- 14 Incubate using the following protocol:

Cycle step	Temperature	Time	No. of cycles
Reverse transcription and strand-switching	42° C	90 mins	1
Heat inactivation	85° C	5 mins	1
Hold	4° C	∞	

Thermal cycler settings for reverse transcription and strand switching

PCR amplification








- 15 In four new  200 μL PCR tubes, prepare the following reaction at  Room temperature in a  50 μL reaction:
-  25 μL 2X Q5 Hot Start High-Fidelity Master Mix
 -  1.5 μL cDNA primer (cPRM)
 -  18.5 μL Nuclease-free water
 -  5 μL Reverse-transcribed cDNA from the previous step (pool, or single sample)

- 16 Amplify using the following cycling conditions:

















Cycle step	Temperature	Time	No. of cycles
Initial denaturation	95 °C	30 secs	1
Denaturation	95 °C	15 secs	10-18*
Annealing	62 °C	15 secs	10-18*
Extension	65 °C	50 secs per kb	10-18*
Final extension	65 °C	6 mins	1
Hold	4 °C	∞	

Thermal cycler settings for PCR amplification

* The recommended starting point is 14 cycles - adjust this depending on experimental needs.

- 17 Add  1 μL of NEB Exonuclease 1 (20 units) directly to each PCR tube to remove unextended primers. Mix by ***pipetting***.
- 18 Incubate the reaction at  37 °C for  00:15:00 , followed by  80 °C for  00:15:00






Bead Cleanup

- 19 Add 160 μ L of resuspended AMPure XP beads to the  1.5 mL tube and mix by ***pipetting***
- 20 Incubate on a gentle agitator (e.g. hula mixer or rotator mixer) for  00:05:00 at  Room temperature
- 21 Spin down  00:00:05 the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.
- 22 Keep the tube on the magnet and wash the beads with  200 μ L of freshly-prepared  70 % volume ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.
- 23 Repeat the previous step: wash with  200 μ L  70 % volume ethanol , and discard the ethanol / wash liquid.
- 24 Spin down  00:00:05 and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for  00:00:30 [at most] but do not dry the pellet to the point of cracking (the magnetic beads should just start to lose their shiny sheen).
- 25 Remove the tube from the magnetic rack and resuspend pellet in  12 μ L of Elution Buffer (EB).
- 26 Incubate at  Room temperature for  00:10:00
- 27 Pellet beads on magnet  00:05:00 until the eluate is clear and colourless
- 28 While still on the magnet, carefully remove and retain  12 μ L of eluate into a clean  1.5 mL Eppendorf DNA LoBind tube




- 29 Quantify 1 µl of the amplified cDNA library using the Quantus Fluorometer using the ONE dsDNA assay (see [ncov 2019 sequencing protocol, step 16](#))

Adapter Addition

- 30 Add  1 µL of Rapid Adapter (RAP) to the amplified cDNA library
- 31 Mix by ***pipetting*** and spin down  00:00:05
- 32 Incubate the reaction for  00:05:00 at  Room temperature
- 33 Store the prepared library  On ice until ready to load onto a flow cell.

Nanopore Sequencing

- 34 Load  20 ng sequencing library onto a MinION flow cell (see [ncov 2019 sequencing protocol, step 21](#))
- 35 Start the sequencing run using MinKNOW, using SQK-PCS109 as the sample preparation protocol (see [ncov 2019 sequencing protocol, step 22](#))

RAMPART Analysis

- 36 Analyse the run results using RAMPART (see <https://artic.network/ncov-2019/ncov2019-using-rampart.html>)