

Mar 23, 2020

Viral Sequencing, from Gunk to Graph (One-Step four-primer PCR alternative)

 Forked from [Viral Sequencing, from Gunk to Graph](#)

DOI

dx.doi.org/10.17504/protocols.io.bd5ni85e

David A Eccles¹


¹Malaghan Institute of Medical Research (NZ)

Coronavirus Method De...



David A Eccles

Malaghan Institute of Medical Research (NZ)



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bd5ni85e

Protocol Citation: David A Eccles 2020. Viral Sequencing, from Gunk to Graph (One-Step four-primer PCR alternative). [protocols.io https://dx.doi.org/10.17504/protocols.io.bd5ni85e](https://dx.doi.org/10.17504/protocols.io.bd5ni85e)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 23, 2020

Last Modified: March 23, 2020

Protocol Integer ID: 34702

Keywords: SARS-CoV-2, COVID-19, nanopore, sequencing,



Abstract

This is a fast "gunk to graph" protocol for analysing viral RNA from nasopharyngeal swabs. The approach involves swab lysis and inactivation at the point of sampling, uses a cellulose binding / wash protocol to reduce extraction cost, incorporates sample-specific barcodes during first-strand synthesis, nanopore rapid-attachment primers during PCR amplification, and nanopore sequencing with parallel RAMPART analysis for fast assembly and phylogenetics.

Note: this is a one-step RT-PCR reaction protocol. For an alternative method that allows more multiplexing flexibility, see the protocol that this has been forked from:

<https://dx.doi.org/10.17504/protocols.io.bd3yi8pw>

Materials

MATERIALS

⊗ Q5 Hot Start High-Fidelity 2X Master Mix - 500 rxns **New England Biolabs Catalog #M0494L**

⊗ MinION sequencer **Oxford Nanopore Technologies**

⊗ ONT MinION Flow Cell R9.4.1 **Oxford Nanopore Technologies Catalog #FLO-MIN106D**







Additional materials TBA.


Safety warnings








⚠ This protocol is UNTESTED, and is in the early stages of development. Do not trust the protocol; question everything.

Assume samples are potentially infectious during extraction, and make sure to use proper sterile technique to avoid cross-contamination.







Swab Lysis

- 1 Prepare a  1.5 mL centrifuge tube with heated lysis buffer and a cellulose disc
- 1.1 Add  500 μL lysis / RNase inactivation buffer ([Twitter reference](#)) to 1.5ml centrifuge tube:
 -  10 millimolar (mM) Tris
 -  10 millimolar (mM) EDTA
 -  0.5 % volume SDS
 -  150 millimolar (mM) NaCl

OR  500 μL extraction buffer #2 (see [paper](#)):




 -  800 millimolar (mM) guanidine hydrochloride
 -  50 millimolar (mM) Tris [pH 8]
 -  0.5 % volume Triton X100
 -  1 % volume Tween-20
- 1.2 Add a  3 mm diameter punched disc from [Whatman #1 filter paper](#) (see [paper](#))
- 1.3 Preheat  1.5 mL tube to  60 °C
- 2 Collect sample using a sterile polystyrene swab with a 30mm breakpoint (e.g. [Puritan 25-3606-U; PurFlock Ultra 6" Sterile Elongated Flock Swab w/Polystyrene Handle, 30mm Breakpoint](#)).








RNA Wash

- 3 Transfer disc to a new  1.5 mL tube containing  200 μL wash buffer using a pipette tip to remove contaminants:
 -  10 millimolar (mM) Tris [pH 8.0]
 -  0.1 % volume Tween-20
- 4 Incubate tube at  Room temperature for  00:01:00







cDNA Synthesis setup

- 5 Transfer disc to a new  200 μL PCR tube using a pipette tip
- 6 Add the following additional components into the  200 μL PCR tube (see the [TaqMan Fast Virus 1-Step Master Mix Product Sheet](#) and [Nanopore protocol for Sequence-specific cDNA-PCR Sequencing \(SQK-PCS109\)](#)) in a  50 μL reaction:

-  1 μL x  2 micromolar (μM) reverse primers
-  1 μL x  2 micromolar (μM) forward primers
-  1.5 μL cDNA primer (cPRM)
-  12.5 μL TaqMan Fast Virus 1-StepMaster Mix
-  34 μL RNase-free water

Reverse primers should be prefixed with sample-specific barcode sequences (if used) and the ONT reverse anchor sequence, i.e. [5' - ACTTGCCTGTCGCTCTATCTTC - [barcode] - [sequence-specific] - 3']

Forward primers should be prefixed with sample-specific barcode sequences (if used) and the ONT forward anchor sequence, i.e. [5' - TTTCTGTTGGTGCTGATATTGC - [barcode] - [sequence-specific] - 3']

- 7 Mix gently **by flicking the tube** and spin down  00:00:05
- 8 Denature RNA and anneal reverse primers at  65 $^{\circ}\text{C}$ for  00:05:00 and then snap cool on a pre-chilled freezer block for  00:01:00

[Note: this step may not be needed]

- 9 Mix gently **by flicking the tube** and spin down  00:00:05

RT and PCR amplification

- 10 In four new 200 μL PCR tubes, prepare the following reaction at Room temperature in a 50 μL reaction:
- 25 μL 2X Q5 Hot Start High-Fidelity Master Mix
 - 1.5 μL cDNA primer (cPRM)
 - 18.5 μL Nuclease-free water
 - 5 μL Reverse-transcribed cDNA from the previous step (pool, or single sample)

- 11 Amplify using the following cycling conditions:

Cycle step	Temperature	Time	No. of cycles
Reverse transcription	50° C	5 mins	1
RT inactivation /initial denaturation	95° C	30 secs	1
Denaturation	95 °C	15 secs	10-40*
Anneal / Extended	62 °C	60 secs	10-40*
Final extension	65 °C	6 mins	1
Hold	4 °C	∞	

Thermal cycler settings for PCR amplification

* The recommended starting point is 20 cycles (i.e. 14 cycles + 6 for direct RNA amplification penalty, see [here](#)) - adjust this depending on experimental needs.

- 12 Add 1 μL of NEB Exonuclease 1 (20 units) directly to each PCR tube to remove unextended primers. Mix by ***pipetting***.








- 13 Incubate the reaction at 37 °C for 00:15:00 , followed by 80 °C for 00:15:00 to








Bead Cleanup

- 14 **Pool liquid from amplified samples** into 1.5 µL Eppendorf DNA LoBind tubes, with no more than 500 µL in each tube
- 15 Add 2X resuspended AMPure XP beads (e.g. 1000 µL XP beads to 500 µL pooled sample) to the 1.5 mL tube and mix by ***pipetting***
- 16 Incubate on a gentle agitator (e.g. hula mixer or rotator mixer) for 00:05:00 at Room temperature
- 17 Spin down 00:00:05 the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.
- 18 Keep the tube on the magnet and wash the beads with 500 µL of freshly-prepared 70 % volume ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.
- 19 Repeat the previous step: wash with 500 µL 70 % volume ethanol , and discard the ethanol / wash liquid.
- 20 Spin down 00:00:05 and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for 00:00:30 [at most] but do not dry the pellet to the point of cracking (the magnetic beads should just start to lose their shiny sheen).
- 21 Remove the tube from the magnetic rack and resuspend pellet in 12 µL of Elution Buffer (EB).




- 22 Incubate at  Room temperature for  00:10:00
- 23 Pellet beads on magnet  00:05:00 until the eluate is clear and colourless
- 24 While still on the magnet, quantify 1 µL of the amplified cDNA library using the Quantus Fluorometer using the ONE dsDNA assay (see [ncov 2019 sequencing protocol, step 16](#))
- 25 While still on the magnet, carefully remove and retain  11 µL of eluate **from each pooled sample** into a clean  1.5 mL Eppendorf DNA LoBind tube

Adapter Addition

- 26 Transfer  11 µL from the total pool into a clean  1.5 mL Eppendorf DNA LoBind tube
- 27 Add  1 µL of Rapid Adapter (RAP) to the amplified cDNA library
- 28 Mix by ***pipetting*** and spin down  00:00:05
- 29 Incubate the reaction for  00:05:00 at  Room temperature
- 30 Store the prepared library  On ice until ready to load onto a flow cell.

Nanopore Sequencing

- 31 Load  20 ng sequencing library onto a MinION flow cell (see [ncov 2019 sequencing protocol, step 21](#))
- 32 Start the sequencing run using MinKNOW, using SQK-PCS109 as the sample preparation protocol (see [ncov 2019 sequencing protocol, step 22](#))



RAMPART Analysis

- 33 Analyse the run results using RAMPART (see <https://artic.network/ncov-2019/ncov2019-using-rampart.html>)