


Apr 17, 2020

Viral Titration of SARS-COV-2 by Plaque Assay (Semi-Solid Agarose)

DOI

<https://dx.doi.org/10.17504/protocols.io.be4zjgx6>



Björn eyer¹

¹Institut Pasteur

Coronavirus Method De...



Bjoern Meyer

Institut Pasteur

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.be4zjgx6>

Protocol Citation: Björn eyer 2020. Viral Titration of SARS-COV-2 by Plaque Assay (Semi-Solid Agarose). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.be4zjgx6>

Manuscript citation:

Gordon, D.E., Jang, G.M., Bouhaddou, M.*et al.* A SARS-CoV-2 protein interaction map reveals targets for drug repurposing. *Nature* (2020). <https://doi.org/10.1038/s41586-020-2286-9>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 15, 2020

Last Modified: April 17, 2020

Protocol Integer ID: 35705

Keywords: plaque assay, viral titration, Covid-19, Covid 19, Coronavirus, viral titration of sar, viral titration, process of plaque assay, plaque assay, sar, titration, assay,

Abstract

This protocol outlines the process of plaque assay for the viral titration of SARS-CoV-2.

Attachments



SARSCoV2_SemiSolid_P

...

277KB

Materials

MATERIALS

⊗ Ethyl Alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**

⊗ MEM **Thermo Fisher Catalog #11095080**

⊗ UltraPure® Agarose **Thermo Fisher Catalog #16500500**

⊗ DMEM, high glucose, GlutaMAX® Supplement, pyruvate **Thermo Fisher Catalog #31966047**

⊗ Gibco™ Trypsin-EDTA (0.05%) phenol red **Fisher Scientific Catalog #11580626**

⊗ Gibco™ Fetal Bovine Serum qualified One Shot™ format **Fisher Scientific Catalog #A3160802**

⊗ Corning™ Costar™ Flat Bottom Cell Culture Plates (24 well) **Fisher Scientific Catalog #10732552**

⊗ Corning™ Costar™ Clear Polystyrene 96-Well Microplates 330µL With Lid **Fisher Scientific Catalog #10360691**

⊗ Formalin solution neutral buffered 10% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #HT501128-4L**

⊗ Crystal violet solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #HT90132**

⊗ VERO C1008 [Vero 76 clone E6 Vero E6] (ATCC® CRL-1586™) **ATCC Catalog #CRL-1586**

	Item	Reference	Storage conditions
	Agarose (Thermofisher)	16500500	RT
	DMEM (Thermofisher)	31966047	4°C
	Trypsin (Gibco)	11580626	4°C
	MEM (Thermofisher)	11095080	4°C
	Fetal Bovine Serum (Thermofisher)	A3160802	-20°C
	24 well plates (Thermofisher)	10732552	RT
	96 well plates round bottom (Thermofisher)	10360691	RT
	ETHANOL ABSOLUTE EXTRA PURE	24103-5L-R	RT
	FORMALIN SOLUTION, NEUTRAL BUFFERED, 10%	HT501128-4L	RT
	CRYSTAL VIOLET SOLUTION	HT90132-1L	RT

Safety warnings

 Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Before start

1. Wear gloves during the entirety of the procedure.
2. Use filtered tips only.
3. Change tips as many times as possible.
4. Decontaminate all tips and pipettes using a solution of diluted bleach.



Preparing Plaque Overlay (2x)

- 1 Mix  455 mL MEM with  20 mL FBS .



Note

The final concentration will be 4% in the 2x overlay.

- 2 Prewarm the **media** to  37 °C .

- 3 Dissolve  0.6 g Agarose in  30 mL H₂O .

Note

The final concentration will be 2% in the 2x overlay

- 4 Melt **Agarose** in the microwave until liquified.

- 5 Quickly add  25 mL melted Agarose to the prewarmed media from **Step 2**.



- 5.1 Close the bottle containing the **Agarose/media** mixture and shake vigorously for  00:00:30 , then again for 3-4 more time over the  00:10:00 time period.



- 6 Let the solution cool at  Room temperature .

Note

Solution can be stored at  4 °C and can be reused when needed.



Day 1

7 Plate 24 well plates with 7.5×10^4 **Vero E6** per well in 10% FBS/DMEM.

8 Let cells grow  Overnight at  37 °C .



Day 2 - Viral Dilutions

9 Set up a 96-well plate in order to dilute your viral solution (serial dilutions).

10 In each well, put  270 μ L MEM (serum-free) .



11 Add  30 μ L viral solution in the first row (A) and mix thoroughly.




Note

Steps 11 to 22 are to be performed in a BSL-3 laboratory/setting.

12 Discard your tips.

13 Transfer  30 μ L previous mix (solution in row A) to the second row (B).





- 14 Repeat that step from row C until row H (SARS-CoV-2 often reaches titers to 10^6).



Day 2 - Viral Infection

- 15 Transfer 250 μ L each dilution into each well .
- 16 Incubate at 37 $^{\circ}$ C for 01:00:00 . Every 00:15:00 , rock the plate a bit.
- 17 While cells are incubating with virus, prewarm 2x overlay media.
- 18 After 1 hour absorption time, add 250 μ L 2x overlay media on top of the 250uL of inoculum (final volume 500uL per well).
- 19 Incubate at 37 $^{\circ}$ C for ~ 65:00:00 .





Day 5

- 20 Prepare a solution of crystal violet diluted in EtOH as follows:  100 mL crystal violey
+  200 mL EtOH +  700 mL water . 
- 21 Put about  00500 μ L formalin 4% into each well on top of the overlay media. 
- 22 Let the well-plate sit for  00:20:00 to  00:30:00 at  Room temperature . 
- 23 Remove **media** from the well plate.
- 24 Add a solution of **0.25% Crystal Violet** and **20% Ethanol** in water to each well in order to stain viable cells.
- 25 Let the well-plate sit for  00:05:00 at  Room temperature . 
- 26 Plunge the well-plate into a first bucket filled with **1% bleach** diluted in water.
- 27 In order to properly clean the plate, plunge the well-plate into a second bucket filled with water.
- 28 Count plaques against a white background. 