

Dec 04, 2022

Whole blood T cell assay for NHPs, containment protocol

 Forked from [Immunophenotyping for NHPs, containment protocol](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1>

Jonathan Audet¹, Courtney Meilleur¹

¹Public Health Agency of Canada



Jonathan Audet

Public Health Agency of Canada

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1>

Protocol Citation: Jonathan Audet, Courtney Meilleur 2022. Whole blood T cell assay for NHPs, containment protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: December 02, 2022

Last Modified: December 04, 2022

Protocol Integer ID: 73500

Keywords: flow cytometry, nonhuman primate, whole blood, broncho-alveolar lavage, T cell assay, protocol for rhesus macaque, using whole blood, assay, concentration of the stimulant, antibody mix, cytokine, cynomolgus macaque, antibody, cells from bronchoalveolar lavage, tc1, stimulant, cell, containment lab

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This is a protocol used to perform T cell assays using whole blood (collected in EDTA tubes) or cells from bronchoalveolar lavage (BAL). We have successfully used this protocol for rhesus macaques, cynomolgus macaques, and African green monkeys (although antibody mix provided was titrated on macaques). It allows the characterization of T cells into Th1 and Th2 (or Tc1 and Tc2) phenotypes. It technically offers 64 different activation profiles (4 cytokines + CD107a + CD69).

The nature and concentration of the stimulant is left unspecified as we have used different ones depending on reagent availability.

This protocol was designed to deal with samples coming from containment labs (> CL2) at our facility. If you are using this protocol at a different facility please ensure that proper testing and approvals are in place. The protocol can be used for experiments completed entirely in CL2, simply go from step 17 directly to 31.

Guidelines

For BAL fluid, make sure to set FSC voltage much lower than for whole blood.

Materials

 Human TruStain FcX **BioLegend Catalog #422302**

 BD FACS Lysing Solution (10X) **Becton Dickinson (BD) Catalog #349202**

 BD Cytofix/Cytoperm **Becton Dickinson (BD) Catalog #554722**

 BD Perm/Wash buffer **Becton Dickinson (BD) Catalog #554723**

 Ghost Dye Red 780 **Tonbo Biosciences Catalog #13-0865-T500**

Equipment

FACS Tube	NAME
Tube	TYPE
Falcon	BRAND
14-959-2A	SKU

Equipment

2 ml screw-cap tubes	NAME
Microtubes	TYPE
Sarstedt	BRAND
72.694.006	SKU

 PBS



Protocol materials

✕ Ghost Dye Red 780 **Tonbo Biosciences Catalog #13-0865-T500**

✕ PBS

✕ Human TruStain FcX **BioLegend Catalog #422302**

✕ BD FACS Lysing Solution (10X) **Becton Dickinson (BD) Catalog #349202**

✕ BD Cytotfix/Cytoperm **Becton Dickinson (BD) Catalog #554722**

✕ BD Perm/Wash buffer **Becton Dickinson (BD) Catalog #554723**

Safety warnings



-

Before start

The antibody mix described in the methods was tested on Cynomolgus macaque whole blood and BAL fluid. All antibodies should cross-react with rhesus macaques and African green monkeys but the panel might need to be re-titrated.

Stimulation

- 1 Aliquot  100 µL whole blood into a tube or well.
- 2 Add  200 µL cRPMI containing 1.5 X the desired concentration of stimulant.
(Stimulant(s) and concentration depend on the experiment and can include anti-CD28 and anti-CD49d)
- 3 Mix and incubate at  37 °C 5% CO₂  Overnight .
- 4 Add monensin and brefeldin A, assuming 1e6 cells.
Anti-CD107a BV421 can be added here as well ( 5 µL)
- 5 Mix and incubate at  37 °C 5% CO₂  06:00:00 .

6h

Preparations

1h 10m

- 6 Prepare the staining mix. (for 1 sample:)

	A	B	C	D	E
	Supplier	Antibody	Clone	Channel	Volume per test
	BD Biosciences	CD3	SP34-2	Alexa Fluor 700	5
	BD Biosciences	CD8	RPA-T8	BV786	5
	BD Biosciences	CD45RA	5H9	PE-CF594	2.5
	BD Biosciences	CD4	L200	PerCP-Cy5.5	20
	BD Biosciences	CCR7	G043H7	BV711	5
	BD Biosciences	CD69	FN50	BV650	5



	A	B	C	D	E
		Total			42.5
		BD Brilliant Stain			57.5

- 7 Prepare the 1X FACS Lysing solution by diluting the 10X stock with Milli-Q water. You will need 2 ml of 1X solution for each sample.

 BD FACS Lysing Solution (10X) **Becton Dickinson (BD) Catalog #349202**

Surface Staining

30m

- 8 Put 5 μ l of TruStain FcX in each tube/well.

 Human TruStain FcX **BioLegend Catalog #422302**

- 9 Incubate 10 min at Room Temperature (RT).  00:10:00  Room temperature

10m

- 10 Add 5 μ l of the  Ghost Dye Red 780 **Tonbo Biosciences Catalog #13-0865-T500**
If staining BAL: Add 5 μ l of 1:100 diluted Ghost Dye.

- 11 Incubate 20 min at RT in the dark.  Room temperature

20m

 00:20:00

- 12 Add stain mix. Mix.  100 μ L

- 13 Incubate 20 min at RT in the dark.  Room temperature

20m

 00:20:00

RBC Lysis



- 14 Add 2 ml of 1X FACS Lysing solution. Vortex immediately, but gently.  2 mL
- 15 Incubate no more than 12 min at RT in the dark.  00:10:00 No more than 12 minutes 10m
- 16 Spin at 300 x g for 5 min.  300 x g, 20°C, 00:05:00
- 17 Decant supernatant.
- 18 Add 2 ml of PBS. Vortex.  2 mL
 PBS
- 19 Spin at 300 x g for 5 min.  300 x g, 20°C, 00:05:00
- 20 Decant supernatant.

Sample Inactivation

30m

- 21 Resuspend in Cytofix/Cytoperm.
 BD Cytofix/Cytoperm **Becton Dickinson (BD) Catalog #554722**
- 21.1 100 µl per 5×10^5 cells.
- 21.2 Use at least 400 µl for easy decanting.
- 22 Incubate at least 30 min at RT in the dark.  00:30:00 30m
 Room temperature In the dark.
- 23 Spin at 500 x g for 8-10 min.  500 x g, 20°C, 00:08:00



24 On a clean bench, decant supernatant.

25 Use same volume of Cytofix/Cytoperm as before to resuspend the cells.

 BD Cytofix/Cytoperm **Becton Dickinson (BD) Catalog #554722**

26 Transfer in a 2 ml screwcap tube. Shake the tube to cover all surfaces with Cytofix/Cytoperm.

Equipment

2 ml screw-cap tubes

NAME

Microtubes

TYPE

Sarstedt

BRAND

72.694.006

SKU

27 Transfer tubes from containment space to CL2 space according to the facility's approved protocols/SOPs.

28 Tubes can be opened in CL2 (in a BSC) no less than 30 min after the resuspension (step 17).

(Samples are generally processed the next day; keep at 4 C overnight, in the dark)

Intracellular staining

30m

29 Give tubes a quick spin in a tabletop centrifuge.

30 Ensure all tubes have at least a few hundred microliters of Cytofix/Cytoperm.

31 Transfer the samples into 1 ml of BD  BD Perm/Wash buffer **Becton Dickinson (BD) Catalog #554723** in FACS tubes.

32 Spin at 500 x g for 8-10 min.  500 x g, 20°C, 00:08:00

33 Decant supernatant.

34 Repeat steps 31-33.

35 Prepare the staining mix. (for 1 sample:)
(uses  BD Perm/Wash buffer **Becton Dickinson (BD) Catalog #554723**)

	A	B	C	D	E
	Supplier	Antibody	Clone	Channel	Volume per test
	BioLegend	IFNgamma	B27	APC	2.5
	BioLegend	TNF	MAb11	PE-Cy7	5
	BioLegend	IL-2	MQ1-17H12	Alexa Fluor 488	0.625
	BioLegend	IL-4	8D4-8	PE	2.5
		Total			10.625
		PermWash			89.375

36 Resuspend in 100 µl of intracellular staining mix.

37 Incubate  4 °C  00:30:00

30m



38 Add  1 mL  BD Perm/Wash buffer **Becton Dickinson (BD) Catalog #554723** .

39  500 x g, 4°C, 00:08:00

8m

40 Decant supernatant and blot tubes/plate.

41 Repeat steps 38-40.

42 Resuspend in PBS 1% formaldehyde.

43 Run on flow cytometer.

44 Gating strategy for whole blood: