

Nov 12, 2024

Version 2

Whole-genome CRISPR screening of stably expressing Cas9 cancer organoid lines V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4q6zrvo5/v2>

Tessa Fowler¹, Jade Smith¹, Adam Jackson¹, Agnieszka Andres¹, Emily Souster¹, Hazel Rogers¹, Alexandra Beck¹, Charlotte Beaver¹, Mathew Garnett¹

¹Wellcome Sanger institute



Tessa Fowler

Wellcome Sanger institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4q6zrvo5/v2>

External link: <https://www.sanger.ac.uk/group/garnett-group/>

Protocol Citation: Tessa Fowler, Jade Smith, Adam Jackson, Agnieszka Andres, Emily Souster, Hazel Rogers, Alexandra Beck, Charlotte Beaver, Mathew Garnett 2024. Whole-genome CRISPR screening of stably expressing Cas9 cancer organoid lines. protocols.io <https://dx.doi.org/10.17504/protocols.io.3byl4q6zrvo5/v2> Version created by **Tessa Fowler**



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 17, 2024

Last Modified: November 12, 2024

Protocol Integer ID: 100040

Keywords: Cancer Organoids, CRISPR/Cas9, gRNA Library titration, gRNA Library transduction, Antibiotic titration, Whole-genome CRISPR Screening, organoid, organoids, genome crispr screening, wide human crispr cas9 library, cas9 cancer organoid lines this protocol, expressing cas9 cancer organoid line, cas9 positive cancer organoid line, cas9 cancer organoid line, cas9 cancer organoid, cas9 cancer organoid lines in triplicate, efficient scale up of the organoid culture, crispr, grna copy per cell, positive cancer organoid line, organoid culture, suitable puromycin concentration for the selection, grna library titration, cancer organoid, grna library virus, puromycin antibiotic titration, lentiviral transduction, genome, other grna libraries if the following, antibiotic titration, suitable puromycin concentration, grna delivery, other grna library, available minimal genome, grna library, host organoid, grna copy, size of the grna library

Abstract

This protocol is for whole-genome CRISPR screening of stably expressing Cas9 cancer organoid lines in triplicate using the commercially available minimal genome-wide human CRISPR Cas9 library. The protocol uses lentiviral transduction as a method for gRNA delivery. This method can be adapted for other gRNA libraries if the following is known;

- The number of days required for screening
- The size of the gRNA library
- The required coverage of the library


To allow for efficient scale up of the organoid culture for screening, a 5% suspension culture method is used. This is not essential to be able to perform whole-genome CRISPR screening of stably expressing Cas9 cancer organoid lines, but provides a more scalable, ergonomic and cost efficient culturing method.

Prior to commencing the screen, a puromycin antibiotic titration is used to identify the most suitable puromycin concentration for the selection of Cas9 positive cancer organoid lines transduced with gRNA library virus.

A gRNA library titration is then performed to determine the volume of library virus required to transduce Cas9 cancer organoids at 30% transduction efficiency which is calculated using FACS analysis. To avoid host organoids taking up more than one gRNA copy per cell, an MOI of 30% should be aimed for.

This protocol uses a 21-23 day screen process.

Process Diagram



Process diagram outlining key steps for CRISPR-Cas9 screening in organoids

Guidelines

It is important that all stages are carried out using the stably expressing Cas9 version of each cancer organoid line.

Avoid maintaining organoids in a single-cell state for extended periods, and refrain from freeze/thaw cycles of virus stocks. Always use freshly prepared virus stocks to ensure optimal viral integrity and confirm that the virus concentration is appropriate for the experimental setup.

Adequate dissociation is crucial for effective viral exposure and gene delivery. Be aware that some organoid lines may exhibit lower uptake of viral particles, which could render them unsuitable for this process.

Puromycin titration;

- Thoroughly mix the organoid suspension at each step to achieve an even single-cell distribution before plating, utilising reservoirs and multi-channel pipettes to reduce ergonomic strain and ensure homogenous solutions.
- If required, CellTiter-Glo 2.0 reagent is light-sensitive, so try to avoid exposure when using it.

gRNA library titration;

The volumes of library titrated may differ depending on the batch of lentivirus or library used. To achieve a 30% MOI, there may therefore be a requirement to dilute the virus prior to titrating. For this reason, an increased number of titration points may also be required. This may require some R&D to determine optimal volumes to titrate.

gRNA library transduction and screen;

The gRNA library transduction is carried out in triplicate with 12.5×10^6 cells seeded per replicate. Expansion of the line to a minimum of 40×10^6 cells is required for this process.

Following selection at Day 6, the organoid suspension should remain in puromycin selection media for the remainder of the screen.

Unless otherwise stated, all steps should be performed under sterile conditions in a microbiological safety cabinet.



Materials

- ✕ TrypLE[®] Express Enzyme (1X), no phenol red **Thermo Fisher Catalog #12604021**
- ✕ 1X Dulbecco's Phosphate Buffered Saline (DPBS) **Thermo Fisher Scientific Catalog #14190094**
- ✕ 10mg/ml Puromycin **InvivoGen Catalog #ant-pr-1**
- ✕ CellTiter-Glo(R) 2.0 Assay **Promega Catalog #G9243**
- ✕ Black walled 96 well plate **Fisher Scientific Catalog #10419822**
- ✕ 37% Formaldehyde **AppliChem Catalog #A0823**
- ✕ eBioscience[®] Fixable Viability Dye eFluor[®] 780 **Thermo Fisher Catalog #65-0865-18**
- ✕ 10mg/ml Polybrene **Merck Millipore (EMD Millipore) Catalog #TR-1003-G**
- ✕ Costar[®] 6-well Clear TC-treated Multiple Well Plates Individually Wrapped Sterile **Corning Catalog #3516**
- ✕ Falcon Round bottomed 5ml tube with cell strainer lid **VWR International (Avantor) Catalog #734-0001**
- ✕ Corning[®] 50ml mini bioreactor **Corning Catalog #431720**
- ✕ Trypsin-EDTA (0.25%) phenol red **Thermo Fisher Scientific Catalog #25200072**
- ✕ Falcon[™] 15mL Conical Centrifuge Tubes **Fisher Scientific Catalog #14-959-53A**
- ✕ Y-27632 ROCK Inhibitor **Selleckchem Catalog #S1049**
- ✕ Corning[®] 125 mL Polycarbonate Erlenmeyer Flask with Vent Cap **Corning Catalog #431143**
- ✕ Corning[®] 250 mL Polycarbonate Erlenmeyer Flask with Vent Cap **Corning Catalog #431144**
- ✕ Microcentrifuge tube Safe-Lock write-on 1.5mL Eppendorf Tube **Eppendorf Catalog # 0030 120.086**
- ✕ Ultra-low attachment 6-well plate **Corning Catalog #CLS3471**
- ✕ Corning[®] Ultra-Low Attachment 75cm² U-Flask Canted Neck Cell Culture Flask with Vent Cap **Corning Catalog #3814**
- ✕ Falcon 50mL Conical Polypropylene Centrifuge Tube with Flat Screw Cap, Sterile **Scientific Laboratory Supplies Ltd Catalog #352070**
- ✕ Cultrex Reduced Growth Factor Basement Membrane Extract, Type 2 **Bio-Techne Catalog #3533-005-02**
- ✕ MinLibCas9 Library⁺ **addgene Catalog #164896**
- ✕ Corning[®] Cell Recovery Solution **Corning Catalog #CLS354253**

Organoid specific culture media- made to specification

Transduction media - organoid specific culture media + Y-27632 ROCK Inhibitor (final concentration 2.5 micromolar (μM))

Equipment

🔥 37 °C 5% CO₂ Incubator



Light microscope
Microbiological safety cabinet (CLASS II)
Centrifuge
Cell counter
Chemical fume hood
Plate reader
Pipette boy
Stripette
Pipettes and tips
Multichannel pipette and tips

Safety warnings



- Lentiviral vectors can infect human cells. However, they are not able to replicate, so the pathogenicity is considered negligible and the risk is reduced by ensuring the correct use of PPE i.e. correct gloves, lab coat and eye protection.
- It is recommended to use centrifuge buckets safety caps and to ensure the caps are only be opened in a microbiological safety cabinet to prevent any biohazard aerosol exposure post-centrifugation.
- Many steps involve repetitive movements. Ensure the use of multichannel and electronic pipettes where possible to reduce ergonomic stress.
- All lentiviral waste should be inactivated and disposed of via recommended local waste routes.

Before start

- If required, ensure that all media is pre-warmed to room temperature before use.
- Where required, prepare an aliquot of 1 mg/ml puromycin (working concentration) by diluting a 10 mg/ml stock 1:10 with sterile water.
- On Day 4 of the puromycin titration process, bring CellTiter-Glo 2.0 reagent to room temperature. (It is light-sensitive and so it is advisable to keep the reagent covered at all times to avoid exposure to light when using).
- Prior to gRNA library transduction, thaw an aliquot of (10 mg/ml) polybrene
- Where required, thaw the required volume of lentivirus for each process stage.
- When fixing organoids, dilute 37% formaldehyde solution 1:10 with DPBS and store at  4 °C



Puromycin titration

3d 0h 22m

1 Day 1: Titration plate set up

Note

This assay is set up using previously expanded organoids which are stably expressing Cas9.

- 1.1 Pre-warm organoid specific culture media to room temperature.
- 1.2 Aspirate media from each well of the organoid culture plate and add  2 mL TrypLE to each well.
- 1.3 Using a cell-scraper, detach BME2 drops containing the cancer organoids from the plate and transfer the organoid suspension to an appropriately sized tube.
- 1.4 Pipette the suspension up and down multiple times to dissociate organoids from the BME2.
- 1.5 Incubate at  37 °C 5% CO₂.
- 1.6 Check the organoid suspension under the microscope every  00:15:00 , to assess and monitor the dissociation of the organoids.

15m

Note

Mix the suspension thoroughly prior to each check to help dissociate the organoids. Stop the incubation once the organoids have broken down to single cells or a maximum time of

 01:00:00 has passed.

- 1.7 Centrifuge at  800 x g for  00:02:00 .

2m

1.8 Aspirate the supernatant and resuspend the pellet in 10 mL of organoid specific culture media.

1.9 Perform a cell count to calculate the total number of cells.

Note

A minimum of 12.5×10^6 cells are required at this stage to ensure enough cells remain for transfer to 5% BME2 suspension culture (step 1.17).

1.10 Resuspend 2.4×10^6 cells in 2.7 ml of organoid specific culture media + $300 \mu\text{L}$ BME2 this will give a final seeding density of 8×10^4 cells per well once plated, Rows B-G of Fig 1.

1.11 Prepare a control stock solution containing organoid specific culture media with 5% BME2 (Row A of Fig 1).

1.12 Set up the titration plate as detailed in Fig 1 below;




Fig 1: Puromycin titration plate set up
 Row A = 200 µl control stock solution per well (step 1.11)
 Rows B - G = 100 µl cell suspension (step 1.10)

Note

- Always seed 3 wells per row as the titration is carried out in triplicate.
- A 96-well plate can be used to titrate up to 4 organoid lines at a time.

1.13 Incubate the plate at  37 °C 5% CO₂ for  00:10:00 to allow the BME2 to polymerise.

10m

1.14 Using a 1 mg/ml stock, prepare puromycin antibiotic solutions at 2x concentrations in organoid specific culture media in tubes as outlined below (Table 1).

Safety information

Puromycin is toxic if swallowed and harmful in contact with skin.

Row (Fig 1)	2 x puromycin concentration (µg/ml)	Volume 1 mg/ml puromycin (µl)	Media volume (µl)	Total volume (µl)	Final plated puromycin concentration (µg/ml)
B	0	0	2500	2500	0
C	2	5	2495	2500	1
D	4	10	2490	2500	2
E	6	15	2485	2500	3
F	8	20	2480	2500	4
G	10	25	2475	2500	5

Table 1: 2 x puromycin concentrations using 1 mg/ml puromycin stock.

Note

- Prepare a minimum of 2.5 ml of each 2x antibiotic so that the volume is adequate for loading a multi-channel pipette without bubbles.
- Antibiotic dilutions should be prepared fresh on the day and kept at  4 °C until required.

1.15 Remove the plate from the incubator and pipette  100 µL of the relevant 2x puromycin stock into each well (Rows B-G) to achieve the final required puromycin concentration according to the plate layout (Fig 1).



1.16 Incubate the plate for  72:00:00 at  37 °C 5% CO₂.

1.17 Plate any remaining cells from step 1.10 into a 5% BME2 suspension culture by seeding cells into a solution of  19 mL organoid specific culture media +  1 mL BME2 using an ultra-low attachment T75 flask.

Note

A minimum of 10×10^6 cells are required at this stage to ensure survival in 5% suspension culture.

1.18 Incubate the ultra-low attachment T75 flask at  37 °C 5% CO₂.

Note

The seeded ultra-low attachment T75 flask is to be used for gRNA library titration set up (see section 3).

2 Day 4: Assessing viability using CellTiter-Glo 2.0 assay

2.1 Thaw CellTiter-Glo 2.0 reagent and equilibrate to room temperature prior to use.

Safety information

CellTiter-Glo is harmful to aquatic life with long lasting effects.

Note

- CellTiter-Glo reagent can be stored at -20 °C and is stable for up to 4 freeze-thaws; thawed reagent can be kept at 4 °C for up to 5 months.
- CellTiter-Glo is light-sensitive so should be kept covered, and used in a cell culture hood with the light off where possible.

2.2 Run a CellTiter-Glo 2.0 viability assay following the manufacturer's instructions.

 celltiterglo-2-0-assay-protocol.pdf

Note

It has been found that adding the CellTiter reagent at 1:5 rather than the manufacturer recommended 1:2 dilution works efficiently.
It is recommended by the manufacturer to use white plates. However, the luciferase signal was found to be too strong so using black plates is recommended.

2.3 Using the luminescence data plot a kill curve to ascertain the lowest concentration of puromycin which results in approximately 100% cell death after 72 hours. (Fig 2).

Note

To create the kill curve;

- Average the triplicate luminescence values to get a single value for each condition.
- Subtract the average background luminescence (Row A, media only) from the other averaged conditions (Rows B-G).
- Divide the average luminescence minus background for each concentration 1, 2, 3, 4 and 5 $\mu\text{g/ml}$ by the 0 $\mu\text{g/ml}$ average to obtain a relative percentage viability.
- Plot these percentage viabilities on a graph against the puromycin concentration.

Expected result




Fig 2: Example Kill curve showing a 'kill concentration' of 3 µg/ml

Guide RNA library titration of Cas9 expressing organoid lines

4m

3 Day 1: gRNA library titration set up

- 3.1 Pre-warm organoid specific culture media to room temperature and thaw a 10 mg/ml aliquot of polybrene.
- 3.2 Prepare transduction media by adding 12.5 µL of [M] 10 millimolar (mM) Y-27632 (Rock inhibitor) to 50 mL organoid specific culture media ([M] 2.5 micromolar (µM) final concentration).

**Note**

Organoids need to remain in transduction media during the setting up and plating of the guide RNA library titration stages within this protocol (steps 3.9 to 4.4).

- 3.3 Collect the organoid suspension culture (from step 1.18) in a 50 ml tube and centrifuge at  800 x g for  00:02:00 .

2m

- 3.4 Aspirate the supernatant and resuspend the pellet in  15 mL to  20 mL of TrypLE.

- 3.5 Pipette the suspension up and down multiple times to dissociate organoids from the BME2.

- 3.6 Incubate at  37 °C 5% CO₂.

- 3.7 Check the suspension under a microscope every 15 minutes to assess and monitor the dissociation of the organoids. Mix the suspension thoroughly prior to each check to help dissociate the organoids.

Note

Generally the suspension becomes cloudy once the majority of organoids are dissociated to small clumps of cells. Stop the incubation once the organoids have broken down to single cells or a maximum time of  01:00:00 has passed.

- 3.8 Centrifuge at  800 x g for  00:02:00 .

2m

- 3.9 Aspirate the supernatant and resuspend the pellet in an appropriate volume of transduction media.

Note

15-20 ml is advised.

- 3.10 Perform a cell count to calculate the total number of cells.
- 3.11 Prepare a preparation mix using the organoid suspension and transduction media to achieve a final concentration of 5.6×10^6 cells in a total volume of 11.2 mL plus 14 μL 10 mg/ml polybrene.

Note

This protocol provides the volumes required to set up a 6 point titration, plus 1 well of dead volume, seeding 8×10^5 cells/well with a final polybrene concentration of 10 $\mu\text{g/ml}$. If more or less points are required, adjust volumes accordingly.

- 3.12 Prepare 6x bioreactor tubes as outlined below in table 2 to achieve 8×10^5 cells per tube, adding appropriate transduction media and library virus volumes.

Bioreactor tube	Preparation mix (μL)	Library virus (μL) (Neat)	Transduction Media volume (μL)	Total volume added (μL)
1	1600	0	400	2000
2	1600	5	395	2000
3	1600	10	390	2000
4	1600	20	380	2000
5	1600	60	340	2000
6	1600	120	280	2000

Table 2: Library titration volumes.

- 3.13 Thaw an aliquot of the gRNA library to room temperature. If the gRNA library is not being used neat, dilute the required volume of the gRNA library in organoid complete culture media.

Note

- Any dilution made to the library at this stage will also need to be made when carrying out the gRNA library transduction at screen. [go to step #6.10](#)
- Once defrosted, the gRNA library should be used within 1 hour.
- Avoid freeze/thaw cycles of the gRNA library.

**Safety information**

Lentiviral vectors can infect human cells. Ensure correct use of PPE and utilise recommended waste routes to reduce the risk.

- 3.14 Add the appropriate volume of gRNA library virus (Table 2) to each bioreactor tube.
- 3.15 Mix well by pipetting and transfer the bioreactor tubes to the incubator at  37 °C 5% CO₂ for overnight incubation.
- 3.16 Passage all remaining organoids in 5% BME2 organoid specific suspension culture media.

Note

The passaged organoids are to be used for gRNA library transduction set up (section 6).

4 Day 2: Plating

- 4.1 Transfer the 6x bioreactor tubes to the centrifuge

Safety information

Lentiviral vectors can infect human cells. Ensure correct use of PPE and utilise recommended waste routes to reduce the risk.

To reduce the risk of aerosols, it is advised where possible that centrifuge buckets are sealed using safety caps and only opened in a microbiological safety cabinet.

- 4.2 Centrifuge at  800 x g for  00:02:00 .

2m

- 4.3 Aspirate the supernatant from the tubes and resuspend each organoid pellet in  1.9 mL of transduction media and add  100 µL of 100% BME2 to make a 5% BME2 solution. Mix well and transfer all 2 ml into each well of a ultra-low attachment 6 well plate (6wp). (Fig 3).




Fig 3: Example plate layout showing virus volumes per well. Virus volumes can be adjusted.

4.4 Incubate at 37 °C 5% CO₂ until Day 6.

Formaldehyde fixation of organoids

4m

5 Day 6: Fixing and staining organoids for flow cytometry analysis.

5.1 Prepare Live/Dead stain solution or antibodies.

Note

Other antibodies can be used but this protocol uses an e780 viability dye. For this reagent prepare a 1:10,000 dilution of e780 dye in PBS. Mix well and store at 4 °C (The solution can be used for 1 week from the time it was prepared).



5.2 Collect the suspension culture from each well of the ultra-low attachment 6wp into individual 2 ml tubes.

5.3 Centrifuge at  800 x g for  00:02:00 .

2m

5.4 Aspirate the supernatant from each tube and resuspend in  1 mL of Trypsin-EDTA (0.25%).

Safety information

Trypsin-EDTA may cause allergic reactions, asthma symptoms, or breathing difficulties if inhaled.

5.5 Incubate for  00:15:00 , until organoids have broken down to single cells.

15m

Note

Mix the solution every few minutes during the incubation. Some lines take longer to dissociate. Do not leave any longer than  00:30:00 .

5.6 Once organoids have broken down to single cells stop the reaction by adding  1 mL (diluting 1:1) of media containing serum.

5.7 Centrifuge at  800 x g for  00:02:00 .

2m

5.8 Aspirate supernatant and resuspend pellets in  200 μ L e780 dye solution.

5.9 Incubate at room temperature for  00:05:00 .

5m

5.10 Add  1 mL PBS to each tube



5.11 Centrifuge at  800 x g for  00:02:00 .

2m

5.12 Aspirate the supernatant and resuspend in  500 μ L of 3.7% formaldehyde. Mix well by pipetting to ensure cells are fixed as single cells.

Safety information

3.7% formaldehyde must be prepared and used only in a chemical fume hood, using chemical resistant gloves. Waste must be kept in the fume hood and disposed of via the recommended route.

5.13 Incubate at  4 °C for  00:10:00 .

10m

5.14 Centrifuge at  800 x g for  00:02:00 .

2m

5.15 Carefully aspirate the supernatant (in chemical fume hood).

Note

Organoid pellets may become transparent and therefore difficult to see. It may also be sticky so can easily stick to pipette tips.

5.16 Resuspend the pellet in  500 μ L (dependent on pellet size) PBS or alternative FACs buffer, and store at  4 °C until ready for analysis by flow cytometry.

5.17 Using collected flow cytometry data create a plot showing side scatter area vs BFP+ expression, gating the positively expressed BFP+ organoids.

5.18 For each titration point create an overlay plot (Fig 4) using your un-transduced control (well 1 of Fig 3).

Expected result




Fig 4: Example flow cytometry analysis of percentage of BFP (%BFP) expression. Blue = tested virus concentration. Red = e780 un-transduced control.

- 5.19 Using the %BFP expression for each titration point, plot a titration curve to calculate the volume of virus required to obtain 30% BFP expression. (Fig 5).

Note

A 30% transduction efficiency is aimed for to try and ensure that each host cell has only taken up one gRNA.

Expected result




Fig 5: Example gRNA library titration curve. For 30% BFP expression a volume of 41 µl of virus per 2 ml would be required.

Guide RNA library transduction and screen of Cas9 expressing organoid lines

2m

6 Day 1: gRNA library transduction

- 6.1 Pre-warm organoid specific culture media to room temperature. Thaw a 10 mg/ml aliquot of polybrene and a [M] 10 millimolar (mM) aliquot of Y-27632 (Rock inhibitor).
- 6.2 Collect the organoid suspension culture from step 3.16 [go to step #3.16](#) in a 50 ml tube and centrifuge at 800 x g for 00:02:00 . 2m
- 6.3 Aspirate the supernatant and resuspend the pellet in 30 mL of TrypLE.
- 6.4 Pipette the suspension up and down multiple times to dissociate organoids from the BME2.
- 6.5 Incubate at 37 °C 5% CO₂.
- 6.6 Check the organoid suspension under the microscope every 15 minutes to assess and monitor the dissociation of the organoids. Mix the suspension thoroughly prior to each



check to help dissociate the organoids.

Note

Generally the suspension becomes cloudy once the majority of organoids are dissociated to small clumps of cells. Stop the incubation once the organoids have broken down to single cells or a maximum time of  01:00:00 has passed.

6.7 Centrifuge at  800 x g for  00:02:00

2m

6.8 Aspirate the supernatant and resuspend the pellet in an appropriate volume of organoid specific culture media.

Note

25 -35 ml is advised

6.9 Perform a cell count to calculate the total number of cells.

Note

- For this protocol, using the minimal genome-wide human CRISPR Cas9 library requires a minimum of 40×10^6 cells at this stage.
- The number of cells required may vary for transduction using alternative gRNA libraries as the number of cells required at transduction is dependent on library size and required coverage.

6.10 Thaw an aliquot of the gRNA library to room temperature. If the gRNA library is not being used neat, dilute the required volume of the gRNA library in organoid complete culture media.

Note

- Any dilution made to the library at this stage must match the dilution carried out when setting up the gRNA library titration.  [go to step #3.13](#)
- Once defrosted, the gRNA library should be used within 1 hour.
- Avoid freeze/thaw cycles of the gRNA library.

Safety information

Lentiviral vectors can infect human cells. Ensure correct use of PPE and utilise recommended waste routes to reduce the risk.

- 6.11 Carry out each library transduction in triplicate. Prepare each replicate in an appropriate size erlenmeyer flask as per Table 3. Mix gently by pipetting.

This protocol uses the minimal genome-wide human CRISPR Cas9 library requires 12.5×10^6 cells per replicate at transduction.

	Number of cells	Volume of polybrene (10 µg/ml)	Volume of ROCKi (Y-27632) (10 mM)	Volume of gRNA library virus	Volume of cell suspension	Total volume
Transduction set up per replicate	12.5×10^6	34 µl	8.5 µl	Library titration result scaled for transduction labware size	34 ml minus volume of gRNA library virus	34 ml

Table 3: Replicate transduction mixture set up for minimal genome-wide human CRISPR Cas9 library. The volume of library virus required to obtain 30% transduction efficiency will be previously established following library titration and scaled to account for transduction labware size at screen.

- 6.12 Prepare an un-transduced control using the appropriate volumes of reagents listed in Table 4. Approximately 1.6×10^6 cells should be seeded into a 50 ml bioreactor tube.

Number of cells	Volume of polybrene (10 µg/ml)	Volume of ROCKi (Y-27632) (10 mM)	Volume of cell suspension	Total Volume
1.6×10^6	4 µl	1 µl	4 ml	4 ml

Table 4: Un-transduced control transduction mixture set up.

- 6.13 Incubate the control bioreactor and replicate flasks prepared in step 6.11 and 6.12 overnight at  37 °C 5% CO₂.



7 Day 2: Seeding

7.1

2m

Collect and transfer the replicate solutions into labelled 50 ml tubes and centrifuge alongside the control bioreactor tube at 800 x g for 00:02:00 .

Safety information

Lentiviral vectors can infect human cells. Ensure correct use of PPE and utilise recommended waste routes to reduce the risk.

To reduce the risk of aerosols, it is advised where possible that centrifuge buckets are sealed using safety caps and only opened in a microbiological safety cabinet.

7.2

For the replicates, aspirate the supernatant and resuspend each replicate pellet in 19 mL organoid specific culture media with 5 µL of 10 millimolar (mM) Y-27632 (Rock inhibitor) and 1 mL of BME2 to make a 5% BME2 solution.

7.3

Mix the suspension well by gentle pipetting and transfer each replicate into 1x ultra low-attachment T75 flask and incubate at 37 °C 5% CO₂.

7.4

For the un-transduced control sample, aspirate the supernatant from the tube and resuspend the pellet in 1.9 ml of organoid specific culture media containing 0.5 µL 10 millimolar (mM) Y-27632 (Rock inhibitor) and add 100 µL of BME2 to make a 5% BME2 solution.

7.5

Mix well by pipetting and transfer into 1 well of an ultra-low attachment 6wp and incubate at 37 °C 5% CO₂.

8 Day 6: Puromycin selection and flow cytometry analysis.

8.1

Using a pipette, collect 0.5 - 1 ml of organoid suspension from each replicate and transfer to individual 2 ml tubes. Collect all of the control and transfer to another 2 ml tube, this is done to test transduction efficiency using flow cytometry.

- 8.2 Fix each collected organoid sample and analyse by flow cytometry following steps 5.1 to 5.17. [↩ 5.1](#)

Note

The un-transduced control once fixed should be kept at [🌡 4 °C](#) for re-analysis at Day 21-23 (Assessment of final selection efficiency).

- 8.3 The pass rate for transduction efficiency is 30% (previously calculated from step 5.19) [↩ go to step #5.19](#) . However, a value between 15-50% is acceptable. This is to try and ensure only one gRNA has been transduced per cell. If a value outside of this range is obtained at this stage the screen should be failed (Fig 6). If the pass rate has been reached, continue to step 8.4.

Expected result




Fig 6: Example flow cytometry analysis of transduction efficiency for screen replicates showing %BFP expression. Red = transduced replicate. Blue = e780 un-transduced control.

- 8.4 For each replicate; collect the suspension culture in a 50 ml tube.
- 8.5 Mix the suspension well by pipetting to break down any larger clumps of BME2/aggregates

before centrifugation.

8.6 Centrifuge at  800 x g for  00:02:00 .

2m

Safety information

Lentiviral vectors can infect human cells. Ensure correct use of PPE and utilise recommended waste routes to reduce the risk.

To reduce the risk of aerosols, it is advised where possible that centrifuge buckets are sealed using safety caps and only opened in a microbiological safety cabinet.

8.7 Aspirate the supernatant and resuspend the pellet in  10 mL organoid specific culture media.

8.8 Prepare an ultra-low attachment T75 flask with the appropriate remaining volume of media, BME2 and puromycin (Table 5 and 6) (Puromycin should be added at a concentration previously calculated in  [go to step #2.3](#) Fig 2).

Note

Ultra-low attachment T75 flasks have a working volume of between 20 - 40 ml. This is dependent on the confluency of the organoid suspension culture.

Safety information

Puromycin is toxic if swallowed, harmful in contact with skin

Total volume per flask (ml)	Cell suspension (ml)	Volume of media added to flask (ml)	Volume of BME2 (ml)
20	10	9	1
30	10	18.5	1.5
40	10	28	2

Table 5: Volumes for preparing 5% BME media solution



Final puromycin concentration (µg/ml)	Flask volume 20 ml	Flask volume 30 ml	Flask volume 40 ml
1	20 µl	30 µl	40 µl
2	40 µl	60 µl	80 µl
3	60 µl	90 µl	120 µl
4	80 µl	120 µl	160 µl
5	100 µl	150 µl	200 µl

Table 6: Puromycin selection volumes to add per specific total flask volumes.

8.9 Mix the organoid suspension well by pipetting, then transfer the organoid suspension to the individual ultra-low attachment T75 flasks.

8.10 Incubate at  37 °C 5% CO₂.

9 Day 9-20: Maintenance of screens.

Note

Media change, or passage 1:2 during week 3 as appropriate.
Do not discard any organoids at this point and re-seed all at passage. Screens are maintained in puromycin until the point of pellet.

9.1 Collect suspension for each replicate in 50 ml tubes.

9.2 Centrifuge at  800 x g for  00:02:00 .

2m

Safety information

To reduce the risk of aerosols, it is advised where possible that centrifuge buckets are sealed using safety caps and only opened in a microbiological safety cabinet.

9.3 Aspirate the supernatant and resuspend the pellet in an appropriate volume of TrypLE.

Note

Small pellets (less than 2 ml in size) can be resuspended in up to 10 ml, whilst larger pellets may need to be resuspended in up to a maximum of 40 ml per tube.



9.4 Mix well by pipetting and incubate at 37 °C 5% CO₂ for 00:10:00 .

10m

Note

Check organoid suspension under the microscope after 5 minutes to assess and monitor the dissociation of the organoids. Pipette the organoid suspension up and down to help dissociate the organoids. Generally the suspension becomes cloudy once the majority of organoids are dissociated to smaller clumps of cells. Do not leave any longer than

00:30:00 .

9.5 Centrifuge at 800 x g for 00:02:00 .

2m

Safety information

To reduce the risk of aerosols, it is advised where possible that centrifuge buckets are sealed using safety caps and only opened in a microbiological safety cabinet.

9.6 Aspirate the supernatant and resuspend each pellet in 20 mL organoid specific culture media.

9.7 For a 1:2 split, prepare 2x ultra-low attachment T75 flasks with the appropriate remaining volume of media, BME2 and puromycin per replicate (Table 5 and 6). 8.8

9.8 Mix the organoid suspension well by pipetting, then transfer half the organoid suspension to each new ultra-low attachment T75 flask.

9.9 Incubate at 37 °C 5% CO₂.

10 Day 21-23: Assessment of final selection efficiency

Note

All steps in section 10 should be completed prior to, but can be completed on the same day as section 11.

- 10.1 Collect 0.5-1 ml of organoid suspension from each replicate and transfer to separate 2 ml tubes.
- 10.2 Fix each collected organoid sample and analyse by flow cytometry following steps 5.1 to 5.17. [➡ go to step #5.1](#)

Note

The un-transduced control, fixed at Day 6 [➡ go to step #8.2](#) and kept at 4 °C should be re-analysed alongside the fixed replicates at Day 21-23.

- 10.3 The pass rate for selection efficiency is between 60-100%. This ensures efficient selection of organoids that have been transduced with the viral library. If a value below this range is obtained at this stage the screen should be failed. (Fig 7). If the pass rate has been reached continue to step 11.1.

Expected result




Fig 7: Example flow cytometry analysis of screen replicates at assessment of final selection efficiency showing %BFP expression. Blue = transduced replicate. Red = e780 un-transduced control.

11 Day 21-23: Pelleting using cell recovery solution

- 11.1 Collect the organoid suspension for each replicate in 50 ml tubes.



11.2 Centrifuge at  800 x g for  00:02:00 .

2m

Safety information

To reduce the risk of aerosols, it is advised where possible that centrifuge buckets are sealed using safety caps and only opened in a microbiological safety cabinet.

11.3 Aspirate the supernatant and resuspend each pellet in an appropriate amount of cell recovery solution (up to  30 mL per replicate) to remove BME. Mix well by pipetting.

11.4 Incubate on ice for  01:00:00 .

1h

11.5 Once the sample tubes have been on ice for  00:30:00 ; mix the suspensions well by pipetting and take aliquots from each replicate to perform a cell count.

30m

11.6 Following  01:00:00 on ice, based on the required number of cells per pellet, transfer the required volume of suspension per pellet into 15 ml tubes.

1h

Note

For the minimal library used in this example 25×10^6 cells are recommended per pellet. If more than 15 ml of suspension is required, multiple rounds of centrifugation are advised.

11.7 Centrifuge at  800 x g for  00:02:00 .

2m

11.8 Aspirate the supernatant and resuspend in 10 ml of  4 °C PBS to wash and remove the cell recovery solution from the pellets.

Note

It is advised to dispense the PBS quickly onto the pellet and to avoid mixing the suspension as the pellet will be prone to sticking to the stripette resulting in a loss of organoids.



11.9 Centrifuge at  800 x g for  00:02:00 .

2m

11.10 Aspirate the PBS and store the pellets at  -80 °C ready for downstream processing.

Protocol references

Verity Goodwin, Emily Souster, Charlotte Beaver, Adam Jackson, Rizwan Ansari, Mathew Garnett , Fiona Behan 2020. Guide RNA Library Transduction of Cas9 Cancer Cell Lines. **protocols.io** [dx.doi.org/10.17504/protocols.io.bg2njyde](https://doi.org/10.17504/protocols.io.bg2njyde)

Gonçalves, E., Thomas, M., Behan, F. M., Picco, G., Pacini, C., Allen, F., Vinceti, A., Sharma, M., Jackson, D. A., Price, S., Beaver, C. M., Dovey, O., Parry-Smith, D., Iorio, F., Parts, L., Yusa, K., & Garnett, M. J. (2021). Minimal genome-wide human CRISPR-Cas9 library. *Genome biology*, 22(1), 40. <https://doi.org/10.1186/s13059-021-02268-4>

Price, S., Bhosle, S., Gonçalves, E. et al. A suspension technique for efficient large-scale cancer organoid culturing and perturbation screens. *Sci Rep* 12, 5571 (2022). <https://doi.org/10.1038/s41598-022-09508-y>

Emily Souster, Verity Goodwin, Adam Jackson, Charlotte Beaver, Rizwan Ansari, Fiona Behan, Mathew Garnett 2020. Puromycin titration of cancer cell lines. **protocols.io** [dx.doi.org/10.17504/protocols.io.bg2gjybw](https://doi.org/10.17504/protocols.io.bg2gjybw)

Verity Goodwin, Emily Souster, Charlotte Beaver, Adam Jackson, Rizwan Ansari, Mathew Garnett , Fiona Behan 2020. Guide RNA Library Titration of Cas9 Cell Lines. **protocols.io** [dx.doi.org/10.17504/protocols.io.bgxujxnw](https://doi.org/10.17504/protocols.io.bgxujxnw)