


Oct 05, 2023

🌐 **XTT Assay for Detection of Bacterial Metabolic Activity in water-based Polyester Polyurethane**

📖 [PLOS One](#)

✓ Peer-reviewed method



DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l27zbg1y/v1>

Nallely Magaña-Montiel¹, Luis Muriel-Millán¹, Liliana Pardo-López¹, Liliana Pardo-López²

¹Departamento de Microbiología Molecular, Instituto de Biotecnología, UNAM, Av. Universidad #2001, Col. Chamilpa, 62210 Cuernavaca, Morelos, Mexico;

²UNAM

Nallely Magaña-Montiel: nallely.magana@ibt.unam.mx;

Luis Muriel-Millán: luis.muriel@ibt.unam.mx

Liliana Pardo-López: liliana.pardo@ibt.unam.mx

IBT - protocols workshop



Nallely Magaña-Montiel

Instituto de Biotecnología, UNAM

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l27zbjg1y/v1>

Protocol Citation: Nallely Magaña-Montiel, Luis Muriel-Millán, Liliana Pardo-López, Liliana Pardo-López 2023. XTT Assay for Detection of Bacterial Metabolic Activity in water-based Polyester Polyurethane . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l27zbjg1y/v1>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 23, 2023

Last Modified: October 03, 2024

Protocol Integer ID: 79267

Keywords: XTT, plastic biodegradation, marine bacteria, Gulf of Mexico, polyurethane, microbial biodegradation assay, detection of bacterial growth, bacterial metabolic activity in water, measuring bacterial metabolic activity, bacterial growth, pseudomonas putida kt2440, bacterial metabolic activity, bacterial strain, metabolic activity of the strain, different bacterial strain, based polyester polyurethane, xtt assay for detection, polyester polyurethane, strain, escherichia coli, based plastic dispersion, soluble orange formazan, plastic dispersion

Funders Acknowledgements:

Conacyt PhD Scollarship

Grant ID: 417065

PAPIIT-DGAPA

Grant ID: IG200223

Disclaimer

Spotlight Video

The video below is a supplement with extra context and tips, as part of the protocols.io Spotlight series, featuring conversations with protocol authors.

https://www.youtube.com/embed/TWhfImLyr_s?si=6guVQ8pc30lf7EGY

Abstract

In microbial biodegradation assays, the detection of bacterial growth in water-based plastic dispersions can be difficult to measure using traditional methods because of the turbidity of culture media and the formation of flocculi. Here, we present a protocol for the detection of bacterial growth in Impranil®DLN, a polyester polyurethane (PU) water-based dispersion. By measuring bacterial metabolic activity, as an indicator of cell viability with the water-soluble 2, 3-bis [2-methyloxy-4-nitro-5-sulfophenyl]-2H-tetrazolium-5-carboxanilide (XTT) salt. Viable growing cells, i.e., those cells that can utilize PU as a carbon source, will reduce the yellow-colored XTT to a water-soluble orange formazan by the action of dehydrogenase enzymes of the respiratory chain. For the standardization of the protocol, we used *Pseudomonas putida* KT2440 and *Escherichia coli* BL21 strains as positive and negative controls, respectively. We determined the metabolic activity of the strains grown with citrate or both citrate and impranil as carbon sources. *P. putida* KT2440 showed higher XTT-detected metabolic activity in the presence of PU than when it was grown only with citrate, indicating that the strain also used PU as a carbon source. In contrast, the negative control did not show differences in metabolic activity between the growth conditions. Our protocol can be adapted to different bacterial strains and culture media.

Image Attribution

Laboratory of Marine Biotechnology, Institute of Biotechnology (IBT) UNAM

Guidelines

Practical considerations

- The assay can be adapted for other bacterial strains.
- Also, you can use an automated microplate spectrophotometer to directly incubate and measure the transformation of XTT to formazan continuously (i.e., every hour).
- Glassware preparation: before use in this procedure, clean glassware is washed with the solvents to extract contaminants and avoid plastic residues.



Materials

LEGEND

¹ for preparation of bacterial culture

² for prior preparation of glass material

³ for XTT assay

Reusable materials:

1. Erlenmeyer 50 mL flasks¹
2. Erlenmeyer 250 mL flasks¹
3. 50 mL measuring cylinders¹
4. Glass container for storing solution 1:1 Methanol:Chloroform²
5. Bacteriological loops¹
6. 10 mL pyrex tubes¹
7. 50 mL graduated cylinder¹
8. Ice buckets¹

Disposable materials:

1. Sterile 96 well culture microplates³
2. Sterile polystyrene petri dishes¹
3. Sterile 50 mL polypropylene centrifuge tubes¹
4. Neoprene & Nitrile Chemical resistant Gloves²
5. Pipette tips^{1,3}
6. Nitrile gloves^{1,3}

Equipment:

1. Multichannel pipette 20–200 µL with 12 channels³
2. Fume hood²
3. Digital clock^{1,3}
4. Vortex^{1,3}
5. Epoch™ 2 Microplate Spectrophotometer Biotek®^{1,3}
6. Centrifuge 5415 R Eppendorf¹
7. Incubated shaker Lab Companion¹

Chemicals:

1. XTT sodium salt powder (Sigma-Aldrich) Catalog #X4626³
2. Basal mineral medium (BM) (Composition in step 8)¹

3. Impranil®DLN^{1,3}
4. Luria-Bertani (LB) agar (NaCl 10 g • L⁻¹, peptone 10 g • L⁻¹ and yeast extract 5 g • L⁻¹, ACS grade)¹
5. Chloroform ACS grade²
6. Methanol ACS grade²
7. MgSO₄ ACS grade¹
8. Sodium citrate ACS grade^{1,3}

Strains:

1. *Pseudomonas putida* KT2440 (positive control)^{1,3}
2. *Escherichia coli* BL21 (negative control)^{1,3}

Software:

1. Microsoft excel³
2. Gen5 Data Analysis Software from Epoch™ 2³

Protocol materials

 XTT sodium salt Merck MilliporeSigma (Sigma-Aldrich) Catalog #X4626

Safety warnings

! This procedure involves the use of hazardous chemicals (chloroform and methanol) for the glassware preparation.

XTT contains no substances with occupational exposure limit values.

1. Read the corresponding safety data sheets for each chemical in the procedure.
2. Use personal protective equipment throughout the procedure: nitrile gloves, wear a lab coat and safety glasses.
3. Dispose of all chemical waste in appropriately labeled containers.

Before start

■ Washing glassware

All glassware should be washed twice using 2 mL of methanol-chloroform mixture (1:1) and allowed to dry in a fume hood.



Preparation of bacterial strains

1 Scrape some of the frozen surface of the glycerol stock using a sterile loop and streak the bacteria onto a Luria-Bertani (LB) agar plate.

2 Incubate the culture at  30 °C during  24:00:00

1d

Note

Use adequate incubation temperature for your strain.



3 Take one isolated colony and inoculate it into  5 mL LB broth.

4 Incubate  Overnight at  30 °C with shaking at  180 rpm .



4.1 Measure the optical density (OD) at 600 nm of cultures using a spectrophotometer.

Note

If needed, prepare a dilution of the culture (e.g., 10^{-1}) using fresh LB broth and measure its OD₆₀₀. The OD₆₀₀ of the original culture is calculated by multiplying the obtained OD by the dilution factor.

5 Inoculate 250-mL flasks containing  50 mL of Basal Medium standardize supplemented with Instant Ocean Sea Salt ($0.06 \text{ g} \cdot \text{L}^{-1}$) (BM), peptone ($10 \text{ g} \cdot \text{L}^{-1}$) and yeast extract ($5 \text{ g} \cdot \text{L}^{-1}$) (BM-PY broth) with an aliquot of the previous culture to obtain an OD₆₀₀ of approximately 0.1.



Note

Our research group uses marine salts to standardize the protocol and screen marine bacteria for subsequent toxicity assays of bacterial culture in zebrafish embryos. For your specific requirements, feel free to use an appropriate minimal medium and carbon source. Composition of Basal Medium in $\text{g}\cdot\text{L}^{-1}$: 0.8 K_2HPO_4 , 0.2 KH_2PO_4 , 0.3 NH_4Cl , 0.19 Na_2SO_4 , 0.07 CaCl_2 , 0.005 $\text{FeSO}_4\cdot 7\text{H}_2\text{O}$, 0.16 MgCl_2 , and 0.0002 Na_2MoO_4 (all J.T. Baker® ACS grade).

Citation

Muriel-Millán LF, Rodríguez-Mejía JL, Godoy-Lozano EE, Rivera-Gómez N, Gutierrez-Rios RM, Morales-Guzmán D, Trejo-Hernández MR, Estradas-Romero A, Pardo-López L. (1970) . Functional and Genomic Characterization of a *Pseudomonas aeruginosa* Strain Isolated From the Southwestern Gulf of Mexico Reveals an Enhanced Adaptation for Long-Chain Alkane Degradation. *Front. Mar. Sci.* 6:572.

<https://doi.org/10.3389/fmars.2019.00572>

LINK

- 6 Incubate the new culture at  180 rpm, 30°C for the time needed to reach the exponential growth phase: approximately  04:00:00 (see the note below).

4h



Note

This 4 hour timeframe applies to our experimental and cultural conditions. It is important to be knowledgeable of the exponential growth rate of your bacterial strains.

- 7 Collect the cells by centrifugation at  6000 rpm, 4°C, 00:20:00 and discard the supernatant.

20m





- 7.1 Wash twice the cellular pellet by suspension in  20 mL of sterile 10mM MgSO_4 followed by centrifugation at  6000 rpm, 4°C, 00:10:00 to remove all traces of the old culture medium.
- 7.2 Resuspend the cells in  5 mL of fresh sterile 10mM MgSO_4 and reserve to be used as inoculum for the next steps. Keep the washed cells on ice to facilitate their handling and preparation.

10m



Culture media preparation

- 8 The strains' ability to grow using the commercial PU coating Impranil® DLN as a carbon source is evaluated in Basal Medium, which is always supplemented with Instant Ocean Sea Salt ($0.06 \text{ g} \cdot \text{L}^{-1}$).
- 9 **XTT solution**
Prepare a solution of  XTT sodium salt Merck MilliporeSigma (Sigma-Aldrich) Catalog #X4626 at 2 $\text{mg} \cdot \text{mL}^{-1}$ in BM.
- 10 **BM-citrate solution**
Prepare a solution of 20 mM sodium citrate in BM.

Note

The presence of easily metabolizable carbon sources such as citrate has been shown to promote the degradation of xenobiotic compounds ([Fonseca 2011](#), Johnsen et al. 2002).

Citation

Johnsen Anders R., Bendixen Karen, Karlson Ulrich (1970)
. Detection of Microbial Growth on Polycyclic Aromatic Hydrocarbons in Microtiter Plates by Using the Respiration Indicator WST-1.
Applied and Environmental Microbiology.

doi: [10.1128/AEM.68.6.2683-2689.2002](https://doi.org/10.1128/AEM.68.6.2683-2689.2002)

[LINK](#)

- 11 **BM-citrate-Impranil solution**
Prepare a BM-citrate solution added with Impranil®DLN ($1 \text{ mg} \cdot \text{mL}^{-1}$)



XTT experiment for detection of bacterial growth

12 XTT experiment standardization



Add  150 μL of the BM solution and  50 μL of the XTT solution into a 96-well microplate and measure the UV-Vis spectrum in a range of 300 to 700 nm with a microplate spectrophotometer. Also, obtain the spectrum of BM-citrate and BM-citrate-Impranil in the same UV-Vis range. Perform in triplicate.

Note

Perform experiment standardization for each condition, i.e., for each medium. The signal intensity is subject to several parameters, including concentration of formazan salts, incubation time, number of viable cells, and cellular metabolic activity. Optimization of test conditions requires consideration of all these parameters (Ghasemi et al. 2021, Riss et al. 2023).

Citation

Ghasemi M, Turnbull T, Sebastian S, Kempson I (1970)
. The MTT Assay: Utility, Limitations, Pitfalls, and Interpretation in Bulk and Single-Cell Analysis.

<https://doi.org/10.3390/ijms222312827>

LINK

Citation

Terry L Riss, PhD, Richard A Moravec, BS, Andrew L Niles, MS, Sarah Duellman, Hélène A Benink, Tracy J Worzella, and Lisa Minor
(1970). Cell Viability Assays. Assay Guidance Manual .

<https://www.ncbi.nlm.nih.gov/books/NBK144065/>

LINK

13 To obtain the maximum absorbance range of XTT with viable cells:

Prepare *Pseudomonas putida* KT2440 (positive control) cultures in BM-citrate, BM-citrate-Impranil, and BM without any carbon source (biotic control), by inoculating  20 mL of each medium on 50-mL flasks up to obtain an OD₆₀₀ approximately of 0.3 ($\sim 1 \times 10^8$ cells) in triplicate.


To obtain the maximum absorbance range, add  150 μ L of each culture medium of *Pseudomonas putida* KT2440 and  50 μ L of the XTT solution into a 96-well microplate and measure UV-Vis spectrum in a range of 300 to 700 nm with a microplate spectrophotometer, immediately after the addition of XTT and again after 1-3 hours of incubation at  180 rpm, 30°C in dark.

Note

The assay can be adapted for other bacterial strains.
Also, you can use an automated microplate spectrophotometer to directly incubate and measure the transformation of XTT to formazan continuously (i.e. every hour).

Expected result

Cells with metabolic activity reduce XTT (yellow) to formazan (orange), which has a wavelength range of maximum absorbance between 450-500 nm (Cell Counting & Health Analysis, 2023; Maldonado et al. 2017).



Citation

Maldonado Cubas, Juan & Casañas Pimentel, Rocio & Merlin, Iván & San Martín-Martínez, Eduardo.
(1970)
. La espectroscopia UV-Vis en la evaluación de la viabilidad de células de cáncer de mama.
Lat. Am. J. Phys. Educ..

https://www.researchgate.net/publication/340915352_La_espectroscopia_UV-Vis_en_la_evaluacion_de_la_viabilidad_de_celulas_de_canc
LINK

Citation

Cell Counting & Health Analysis, Merck guides, online 2023 (1970)
 . Cell Viability and Proliferation XTT Assay Protocol Guide.

<https://www.sigmaaldrich.com/MX/es/technical-documents/protocol/cell-culture-and-cell-culture-analysis/cell-counting-and-health-LINK>

Detection of Bacterial Growth in Polyester Polyurethane

- 14 To evaluate microbial growth, measure the absorbance of each treatment every hour with an automated microplate spectrophotometer at Abs_{max} (range of 450nm - 500nm) with a reference wavelength of 630nm immediately after the addition of XTT and again every hour.



Note

The use of a reference wavelength considerably reduces the noise of the particles in the medium (background subtraction at 630–690 nm). Also, subtract the absorbance at time zero (Abs_i) from subsequent readings to obtain the change in absorbance.

$$\text{Absorbance} = [\text{Abs}_{\text{max}} - \text{Abs}_{630}] - \text{Abs}_i$$

Culture media:

	A	B	C	D	E
		Abiotic control	Pos. control	Neg. control	Test microorganism
BM	x	x	x	x	x



	A	B	C	D	E
	BM-citrate	x	x	x	x
	BM-citrate- Impranil	x	x	x	x

"x" indicates that must be included.

Note

Note that BM in positive control, negative control, and test microorganism correspond to biotic controls without any carbon source.

- 15 Add  150 μL of the different cultures into a 96-well microplate containing  50 μL of XTT in each well.

We tested the metabolic activity of *P. putida* KT2440 (Franklin et al. 1981) and *E. coli* BL21 as positive and negative controls, respectively.

Citation

F. C. Franklin, M. Bagdasarian, M. M. Bagdasarian, and K. N. Timmis (1981)
. Molecular and functional analysis of the TOL plasmid pWWO from *Pseudomonas putida*
and cloning of genes for the entire regulated aromatic ring meta cleavage pathway.
Proc.NatLAcad.Sci.USA.

<https://doi.org/10.1073/pnas.78.12.7458>

LINK

**Note**

The assay can be adapted for other bacterial strains.

16 Immediately after XTT addition, incubate the plate  24:00:00 at  30. °C

1d

16.1

To evaluate microbial growth, measure absorbance every hour with an automated microplate spectrophotometer reader EPOCH2 (BioTek Instruments Inc.) both at 470nm and 630nm (reference wavelength for background subtraction).



Expected result

	1	2	3	4	5	6	7	8	9	10	11	12
A	BM-1	BM-1	BM-2	BM-2	BM-3	BM-3	BM-C1	BM-C1	BM-C2	BM-C2	BM-C3	BM-C3
B	BM-IC1	BM-IC1	BM-IC2	BM-IC2	BM-IC3	BM-IC3						
C	BM-1 C-	BM-1 C-	BM-2 C-	BM-2 C-	BM-3 C-	BM-3 C-	BM-C1 C-	BM-C1 C-	BM-C2 C-	BM-C2 C-	BM-C3 C-	BM-C3 C-
D	BM-IC1 C-	BM-IC1 C-	BM-IC2 C-	BM-IC2 C-	BM-IC3 C-	BM-IC3 C-						
E	BM-1 C+	BM-1 C+	BM-2 C+	BM-2 C+	BM-3 C+	BM-3 C+	BM-C1 C+	BM-C1 C+	BM-C2 C+	BM-C2 C+	BM-C3 C+	BM-C3 C+
F	BM-IC1 C+	BM-IC1 C+	BM-IC2 C+	BM-IC2 C+	BM-IC3 C+	BM-IC3 C+						
G												
H												

Expected results for XTT assay. Viable cells transform the yellow XTT to formazan (orange). The signal intensity depends both on the number of viable cells and the cellular metabolic activity. BM-C: BM-citrate, BM-IC: BM-citrate-Impranil, C-: negative control, C+: positive control. The numbers indicate biological replicates each one performed in duplicate (technical replicates).

Step 10

Johnsen Anders R., Bendixen Karen, Karlson Ulrich. Detection of Microbial Growth on Polycyclic Aromatic Hydrocarbons in Microtiter Plates by Using the Respiration Indicator WST-1

doi: 10.1128/AEM.68.6.2683-2689.2002

Terry L Riss, PhD, Richard A Moravec, BS, Andrew L Niles, MS, Sarah Duellman, H  l  ne A Benink, Tracy J Worzella, and Lisa Minor. Cell Viability Assays

<https://www.ncbi.nlm.nih.gov/books/NBK144065/>

Maldonado Cubas, Juan & Casañas Pimentel, Rocio & Merlin, Iván & San Martín-Martínez, Eduardo.. La espectroscopia UV-Vis en la evaluación de la viabilidad de células de cáncer de mama

https://www.researchgate.net/publication/340915352_La_espectroscopia_UV-Vis_en_la_evaluacion_de_la_viabilidad_de_celulas_de_canc

Cell Counting & Health Analysis, Merck guides, online 2023. Cell Viability and Proliferation XTT Assay Protocol Guide

<https://www.sigmaaldrich.com/MX/es/technical-documents/protocol/cell-culture-and-cell-culture-analysis/cell-counting-and-health->

F. C. Franklin, M. Bagdasarian, M. M. Bagdasarian, and K. N. Timmis. Molecular and functional analysis of the TOL plasmid pWWO from *Pseudomonas putida* and cloning of genes for the entire regulated aromatic ring meta cleavage pathway.

<https://doi.org/10.1073/pnas.78.12.7458>

Muriel-Millán LF, Rodríguez-Mejía JL, Godoy-Lozano EE, Rivera-Gómez N, Gutierrez-Rios RM, Morales-Guzmán D, Trejo-Hernández MR, Estradas-Romero A, Pardo-López L. Functional and Genomic Characterization of a *Pseudomonas aeruginosa* Strain Isolated From the Southwestern Gulf of Mexico Reveals an Enhanced Adaptation for Long-Chain Alkane Degradation.

<https://doi.org/10.3389/fmars.2019.00572>