

Oct 23, 2024

Version 1

Zika Virus NS2B-NS3 protease fusion protein his10-tagged protein expression and purification: large scale 1 L cultures V.1

 Forked from [Parallel rapid expression and purification of proteins for crystallography \(PREPX\): large scale 1 L cultures](#)



DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld8o17v5b/v1>

Michael Fairhead¹

¹university of oxford

ASAP Discovery



Michael Fairhead

University of Oxford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld8o17v5b/v1>

External link: <https://orcid.org/0000-0001-5361-3933>

Protocol Citation: Michael Fairhead 2024. Zika Virus NS2B-NS3 protease fusion protein his10-tagged protein expression and purification: large scale 1 L cultures. protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ld8o17v5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 21, 2024

Last Modified: October 23, 2024

Protocol Integer ID: 110435

Keywords: parallel protein purification, Recombinant protein, Escherichia coli, ZIKV, Zika, NS2B-NS3 protease, purification of zika virus ns2b, zika virus ns2b, ns3 protease fusion protein, zika virus, ns3 protease fusion, recombinant protein, tagged protein expression, ns3 protease, protein expression, protein, purification, protease, tag for biophysical assay, escherichia coli, tag cleavage, fusion protein

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

This protocol details the expression and purification of Zika Virus NS2B-NS3 protease fusion protein with a his10-tag for biophysical assay at a 1 L culture scale. Recombinant proteins are expressed in *Escherichia coli* using the autoinduction method and then purified in using a IMAC, desalt, tag cleavage, reverse IMAC and gel filtration work flow.

Attachments



[nyikb5id7.docx](#)

24KB



[PAGE24-00720.pdf](#)

1.5MB

Guidelines

Method overview

Standard workflow is expression via autoinduction followed by purification using IMAC/PD-10/revIMAC and serial gel filtration

Materials

Zika Virus NS2B-NS3 protease fusion protein

Vector: pNIC (Kanamycin resistance)

Cell line: E. coli strain BL21(DE3)-RR

Tags and additions: N-terminal His10

Construct protein sequence

- MHHHHHHHHHHHGSMSGKSVDYIERAGDITWEKDAEVTGNSPRLDVALDESGDFSLVEDDGPPMREGGGGSGGGGG SGALWDVPAPKEVKKGETTDGVYRVMTRLLGSTQVGVGVMQEGVFHTMWHVTKGSALRSGEGRLDPYWGDVKQDL VSYCGPWKLDAAWDGHSEVQLLAVPPGERARNIQTLPGIFKTKDGDIGAVDYPAGTSGSPILDKCGRVIGLYGNGVVIK NGSYVSAITQGRREEETPVERK
- MW = 27.116 kDa
- Extinction co-efficient = $43430 \text{ M}^{-1} \text{ cm}^{-1}$
- pI = 5.67

 PD-10 desalting columns packed with Sephadex G-25 resin **Cytiva Catalog #17085101**

 PD-10 Spin Adapter **Cytiva Catalog #28923245**

 LabMate PD-10 Buffer Reservoir **Cytiva Catalog #18321603**

 His GraviTrap **Cytiva Catalog #11003399**

 Super Broth **Formedium Catalog #SPB0102**

 AirOtop Enhanced Flask seals **Generon Catalog #899425**

 Nalgene®; Unwire®; Test Tube Racks: Resmer®; Manufacturing Technology, for 30mm tubes, white **Thermo Fisher Catalog #5970-0030**

 AIM – Terrific Broth Base including Trace elements **Formedium Catalog #AIMTB0210**

 Ultra Yield 2.5L Flask, Sterile **Generon Catalog #931136-B**

Optional but useful

 BENCHMIXER™ XL MULTI-TUBE VORTEXER **Benchmark Scientific Catalog #BV1010**

Materials (1 L cultures) for Expression:

- Plates with LB-agar+antibiotics
-  1 L of autoclaved autoinduction TB + 20 g/L glycerol + antibiotics
-  1 mL of 10 % Antifoam 204 (Sigma) in ethanol
-  2.5 L Ultra Yield flasks (fitted with loose foil cover**)

Materials (1 L cultures) for Purification:

■ 1L of Base Buffer

	A	B
	HEPES	10 mM
	Glycerol	5%
	NaCl	500 mM
	TCEP, pH 7.5	0.5 mM

-  100 mL of [M] 3 Mass Percent imidazole pH 7.5.
-  100 mL of 10 % Triton X-100 in water.
-  10 mL Lysozyme solution (100 x).
-  1 mL homemade benzonase (1000x). Maybe substituted for 10 mg/mL of commercial DNase I
- 2 x His GraviTrap column per litre of culture to be purified (Cytiva) fitted with LabMate extender (Cytiva) and PD-10 spin adapter (Cytiva) in 24 place Nalgene rack.
- 2 x PD-10 desalting column per litre of culture to be purified fitted with LabMate extender and PD-10 spin adapter in 24 place Nalgene rack.
- 2 × 50 mL centrifuge tubes per litre of culture to be purified in a 24 place Nalgene rack.

Safety warnings

- ! Triton x-100 is currently restricted for use in the EU and cannot be used without an exemption certificate REACH Annex XIV (Jan 2021). It can be readily substituted with IGEPAL CA-630 (which is likely to be subject to the same restrictions in the near future). Alternatives that also maybe used are Tergitol 15-S-9 or Tween-20 or octyl glucoside.



Expression

12h 20m

- 1 Transform *BL21 (DE3)* with the appropriate plasmid and spread onto LB-agar plate + 50 µg/mL kanamycin and incubate Overnight 37 °C *.
- 2 Use several colonies to inoculate 10 mL of superbrotth + 1 % glucose + 50 µg/mL kanamycin in a 50 mL tube and 250 rpm, 37°C, 16:00:00
- 3 Use 10 mL to inoculate 1 L AIM-TB + 50 µg/mL kanamycin + 0.01% Antifoam 204 in a 2.5 L Ultra Yield baffled flask (Thomson) fitted with an AirOtop enhanced flask seal (Thomson)
- 4 Grow 250 rpm, 37°C, 04:00:00 shaking.
- 5 Grow 250 rpm, 18°C, 24:00:00 shaking.
- 6 Harvest at 4000 x g, 4°C, 00:20:00 .
- 7 Scrape out pellet and place in plastic polygrip bag and place in -80 °C freezer. Final wet cell weight is 45 g/L of culture

16h

Cell lysis

3h 30m

- 8 Place the polygrip bag on a flat surface and smash cell pellet into small pieces and transfer the cell pellet into 500 mL beaker.
- 9 Add 3 mL Base Buffer/g cell pellet (10 millimolar (mM) HEPES, 500 millimolar (mM) NaCl, 5 % Glycerol, 0.5 millimolar (mM) TCEP, 7.5) + 0.5 mg/mL Lysozyme, 1 µg/ml Benzonase or 10 µg/ml DNase I, 1 % Triton X-100***, 20 millimolar (mM) Imidazole.



**Note**

***Triton X-100 is currently restricted for use in the EU and cannot be used without an exemption certificate REACH Annex XIV (Jan 2021). It can be readily substituted with IGEPAL CA-630 (which is likely to be subject to the same restrictions in the near future). Alternatives that could be used are Tergitol 15-S-9 or Tween-20 or octyl glucoside.

10 Use stripette to dissolve pellet and transfer  45 mL (maximum) to a 50 mL tube (4 tubes in total).

11 Leave  00:30:00  Room temperature .

30m

12 Place tubes in  -80 °C freezer for 2h or overnight if preferred.



13 Thaw in  Room temperature water bath for  01:00:00 and mix.

1h



14 Centrifuge at  4000 x g, 4°C, 01:00:00 .

1h

**Purification**

3h 30m

15 Apply supernatant from 2 × 50 mL tubes to 1 mL His GraviTrap column (Cytiva) fitted with a LabMate PD-10 Buffer Reservoir (Cytiva). His GraviTrap columns are pre-equilibrated in Base Buffer +  20 millimolar (mM) Imidazole.

16 Wash with  10 mL Base Buffer +  20 millimolar (mM) Imidazole.



17 Slot His GraviTrap column into PD-10 desalting column (Cytiva) fitted with a LabMate PD-10 Buffer Reservoir. PD-10 desalting columns are pre-equilibrated in Base Buffer w/o Imidazole.

18 Elute protein, with  2.5 mL of Base Buffer +  500 millimolar (mM) Imidazole, directly onto PD-10 column.




- 19 Remove His GraviTrap column.
- 20 Place PD-10 desalting column into 50 mL falcon tube and add  3.5 mL Base Buffer w/o Imidazole and collect the flow through from the column. 
- 21 Measure the absorbance at 280nm (A280)
- 22 Concentrate to 20 mg/mL using a 10,000 MWCO centrifugal filter (Amicon, Merck-Millipore) and inject the sample (2 mL) onto 125 mL superose 12 pg column pre-equilibrated in Base Buffer. Collect 2 mL fractions using a flow rate at 1 mL/minute. A superdex 75 pg column or equivalent may also be used.
- 23 Pool the peak fraction(s) (e.g. 19-24 in the chromatogram below) and concentrate to 10 mg/mL using a 10,000 MWCO centrifugal filter (Amicon, Merck-Millipore) Snap freeze protein using liquid nitrogen as 0.1 mL aliquots and store the sample at -80°C until use

Purification


3h 30m

- 23.1 chromatogram using Base Buffer as the mobile phase



ZIKV NS2B-NS3 SEC chromatogram profile on Superose 12 pg 125 mL using Base Buffer as the mobile phase

23.2



SDS-PAGE ZIKV NS2B-NS3 after SEC, Samples run on NuPAGE Bis-Tris 4-12% gels (Invitrogen), 200V 40 minutes using MES running buffer (Invitrogen)

Column regeneration: PD-10



24 Wash PD-10 columns with  50 mL -  100 mL of Milli-Q water.



25 Store all columns in water at  4 °C . For long term storage, use 20 % Ethanol

Column regeneration: His GraviTrap

26 Wash the IMAC columns with  40 mL Milli-Q.



27 Wash IMAC columns  10 mL 20 % Ethanol +  0.1 Mass Percent EDTA*.



*PUT NICKEL WASTE IN APPROPRIATE CONTAINER FOR DISPOSAL!

28 Wash IMAC columns  40 mL Milli-Q.



29 Wash IMAC columns  10 mL  1 Mass Percent NaOH.



30 Wash IMAC columns  40 mL Milli-Q.



31 Wash IMAC columns  10 mL  1 Mass Percent Acetic Acid + 1 % Triton X-100.



32 Wash IMAC columns  40 mL Milli-Q.

33 Wash IMAC columns with  0.5 mL  100 millimolar (mM) Nickel Sulfate +  20 millimolar (mM) Tris.HCl pH 8*.



Note

*PUT NICKEL WASTE IN APPROPRIATE CONTAINER FOR DISPOSAL!



34 Wash IMAC columns with  40 mL Milli-Q.

35 Store all columns in water at  4 °C . For long term storage use 20 % Ethanol