




Nov 23, 2023

DNA extraction protocol for animal blood samples using the E.Z.N.A blood mini kit

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov141xpvr2/v1>



Louis-Stéphane Le Clercq^{1,2}, Désirée ee Dalton³, Antoinette Kotzé^{1,2}, Paul Grobler¹

¹University of the Free State; ²South African National Biodiversity Institute; ³Teesside University

Biological clock measures the association between the circadian and epigenetic clock as predictors of migration and age
Tech. support email: leclercq.l.s@gmail.com



Louis-Stéphane Le Clercq

University of the Free State, South African National Biodive...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov141xpvr2/v1>

External link: <https://sites.google.com/view/lsleclercq/projects/phd-project>

Protocol Citation: Louis-Stéphane Le Clercq, Désirée ee Dalton, Antoinette Kotzé, Paul Grobler 2023. DNA extraction protocol for animal blood samples using the E.Z.N.A blood mini kit. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov141xpvr2/v1>

Manuscript citation:

Le Clercq, L.S., 2023. *Biological clock measures: Assessing the association between the circadian and epigenetic clock as predictors of migration phenology and biological aging in wildlife* (Doctoral thesis, University of the Free State).

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 25, 2021

Last Modified: November 23, 2023

Protocol Integer ID: 51045

Keywords: DNA Extraction, Omega, Manual, Blood, Columns, Animals, dna extraction protocol for animal blood sample, dna extraction protocol, dna extraction, animal blood sample, preparation of genomic dna, purified dna, dna, genomic dna for consistent pcr, genomic dna, anticoagulated whole blood, such as pcr, blood, extraction, whole blood, rapid method for the isolation

Funders Acknowledgements:

National Research Foundation (RSA)

Grant ID: 112062

Abstract

The E.Z.N.A. Blood DNA Mini Kit provides an easy and rapid method for the isolation of genomic DNA for consistent PCR and Southern analysis. Up to 250 µL fresh, frozen, or anticoagulated whole blood can be readily processed at one time. The E.Z.N.A. Blood DNA Mini Kit can also be used for the preparation of genomic DNA from buffy coat, serum, plasma, saliva, buccal swabs, and other body fluids. The E.Z.N.A. Blood DNA Kit allows for single or multiple simultaneous processing of multiple samples. There is no need for phenol/chloroform extractions, and time-consuming steps are eliminated (e.g. precipitation using isopropanol or ethanol). Purified DNA obtained with the E.Z.N.A. Blood DNA Kit is ready for applications such as PCR, restriction digestion, and Southern blotting.

Image Attribution

<https://www.omegabiotek.com/product/e-z-n-a-blood-dna-mini-kit/>

Guidelines

Read and review the included product manual before starting.

Materials

-  E.Z.N.A. Blood DNA kit **Whitehead Scientific Catalog #D3392-01**
-  RNase I - 25,000 units **New England Biolabs Catalog #M0243L**
-  100% Ethanol
-  Isopropanol
- Tabletop microcentrifuge capable of at least 13,000 x g
- Nuclease-free 2 mL microcentrifuge tubes
- Water bath, incubator, or heat block capable of 65°C
- Vortexer

Protocol materials

-  Isopropanol
-  E.Z.N.A. Blood DNA kit **Whitehead Scientific Catalog #D3392-01**
-  RNase I - 25,000 units **New England Biolabs Catalog #M0243L**
-  100% Ethanol

Safety warnings

 None

Ethics statement

Protocol approval for the present study was obtained from the protocol committee of the Department of Genetics, University of the Free State (approval number: Res18/2020). Ethics approvals were obtained from the University of the Free State (approval number: UFS-AED2020/0015/1709) as well as the South African National Biodiversity Institute (approval number: SANBI/RES/P2020/30). Appropriate research permits were also obtained from South African regulatory authorities including the Department of Agriculture, Land Reform, and Rural Development (Section 20 permit: 12/11/1/1/18(1824JD)) and the Department of Environmental Affairs (Threatened Or Protected Species (TOPS) permit: O-52903).

Before start

- Prepare HBC Buffer and DNA Wash Buffer according to the directions.
- Set water bath, incubator, or heat block to 65°C.
- Heat the Elution Buffer to 65°C.

Lyse

1 Transfer the Sample into a sterile microcentrifuge tube and bring the volume up to 250 μ L with 10mM Tris-HCl, PBS, or Elution Buffer (provided).



2 Add 25 μ L Proteinase K Solution (provided), 250 μ L BL Buffer (provided), and 5 μ L RNase I - 25,000 units **New England Biolabs Catalog #M0243L** (50 mg/mL).



3 Vortex at maximum speed for 00:00:15 seconds.

15s



4 Incubate at 65 $^{\circ}$ C for 00:10:00 minutes. Vortex briefly once during incubation.

10m



Heating block set to 65 degrees Celsius.

Adjust binding conditions



- 5 Add  260 μ L  100% Ethanol . Vortex at maximum speed for  00:00:20 seconds.

20s



- 6 Centrifuge briefly to collect any drops from the inside of the lid.



Bind

- 7 Insert a HiBind DNA Mini Column (provided) into a 2 mL Collection Tube.



Green HiBind columns to be placed in clear collection tubes. Can label samples on cap.

- 8 Transfer the entire sample to the column.



- 9 Centrifuge at  11000 x g, 21°C, 00:01:00

1m



10 Discard the filtrate and the Collection Tube. Insert the HiBind DNA Mini Column into a new 2 mL Collection Tube.

11 Add  500 μ L HBC Buffer.



Note

HBC Buffer must be diluted with 100% isopropanol before use.

12 Centrifuge at  11000 x g, 21°C, 00:01:00

1m



13 Discard the filtrate and reuse Collection Tube.

Wash

14 Add  700 μ L DNA Wash Buffer (provided).



Note

DNA Wash Buffer must be diluted with 100% ethanol before use.

15 Centrifuge at  10000 x g, 21°C, 00:01:00

1m



16 Discard the filtrate and reuse the Collection Tube.

17  go to step #14 Repeat once

Dry



18 Centrifuge the empty HiBind DNA Mini Column at  13000 x g, 21°C, 00:02:00 .

2m

Note

It is important to dry the column membrane before elution. Residual ethanol may interfere with downstream applications.

Elute

7m

19 Transfer the HiBind DNA Mini Column into a nuclease-free 2 mL microcentrifuge tube.

20 Add  50 µL Elution Buffer heated to  65 °C .



21 Incubate the HiBind DNA Mini Column at  65 °C for  00:05:00 .

5m



22 Centrifuge at  13000 x g, 21°C, 00:02:00 .

2m



23 **Optional:** Apply filtrate to column and repeat centrifugation.



Store

24 Store eluted DNA at  -20 °C .

Samples extracted using this protocol were submitted to NCBI BioSample and linked to BioProject [PRJNA737185](https://dx.doi.org/10.17504/protocols.io.ewov141xpvr2/v1).

Expected result

Nanodrop One results for samples:

	Sample	A260	A280	A260/A230	A260/A280	ng/uL
	1	0.185	0.096	0.359	1.934	9.26
	2	0.125	0.059	0.732	2.131	6.26
	3	0.159	0.102	0.528	1.558	7.94
	4	0.082	0.047	0.405	1.723	4.08
	5	0.107	0.058	0.815	1.838	5.33
	6	0.110	0.063	1.020	1.742	5.51
	7	0.143	0.070	0.518	2.033	7.14
	8	0.077	0.029	1.443	2.615	3.83
	9	0.060	0.039	1.156	1.534	2.99
	10	0.059	0.032	0.398	1.844	2.95
	11	0.082	0.044	1.039	1.855	4.10
	12	0.119	0.073	1.025	1.637	5.96
	13	0.072	0.032	0.938	2.237	3.60
	14	0.098	0.051	1.051	1.918	4.88
	15	0.092	0.049	1.475	1.869	4.59

The extractions yielded pure DNA (A260/A280 approximately 1.6 to 1.8) with sufficiently high concentrations of around 4 ng/uL.