

Feb 08, 2016

## Isolation of cyanophages by liquid enrichment assay

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.dp65rd](https://dx.doi.org/10.17504/protocols.io.dp65rd)

Mathias Middelboe<sup>1</sup>, Amy M. Chan<sup>1</sup>, and Sif K. Bertelsen<sup>1</sup>

<sup>1</sup>Manual of Aquatic Viral Ecology

VERVE Net

Suttle Laboratory of Mar...



Amy Chan

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dp65rd>

External link: [http://www.aslo.org/books/mave/MAVE\\_118.pdf](http://www.aslo.org/books/mave/MAVE_118.pdf)

**Protocol Citation:** Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen 2016. Isolation of cyanophages by liquid enrichment assay. protocols.io <https://dx.doi.org/10.17504/protocols.io.dp65rd>



**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** September 01, 2015

**Last Modified:** March 28, 2018

**Protocol Integer ID:** 1502

**Keywords:** liquid enrichment assay, isolation of cyanophage, liquid enrichment culture, enrichment cultures for bacteriophage, liquid enrichment, cyanophage, bacteriophage, water sample, use of enrichment culture, larger volumes of water sample, enrichment, assay, enrichment culture, host cell, virus, larger sample volume

## Abstract

Another convenient way to screen larger sample volumes is to perform liquid enrichment cultures (Suttle 1993). By this approach, larger volumes of water samples can be screened for cyanophages, thus enabling detection of "rare" viruses. Similar to the use of enrichment cultures for bacteriophages, the disadvantages include the fact that lysis of the host cells are not always obvious, especially if the initial titer is low or if the host culture is not clonal or unialgal.

## Guidelines

If low titers are expected, the viruses in the sample can be concentrated via TFF to make a virus concentrate (VC) (Suttle et al. 1991; Wommack et al. 2010, this volume). To increase the detection limit, several different VCs can be combined and added to the same culture. Another convenient way to screen larger sample volumes is to perform liquid enrichment cultures (Suttle 1993). By this approach, larger volumes of water samples can be screened for cyanophages, thus enabling detection of rare viruses. Similar to the use of enrichment cultures for bacteriophages, the disadvantages include the fact that lysis of the host cells are not always obvious, especially if the initial titer is low or if the host culture is not clonal or unialgal. Also, more steps are required to dilute out nonreplicating viruses to obtain pure clonal isolates. Because it is an end-point dilution assay, only the most abundant phages will be isolated. The principle is the same as for the isolation of bacteriophages mentioned earlier except that the sample volumes screened are in the order of liters instead of milliliters.

## Troubleshooting

- 1 Prefilter at least 3 L water sample through a glass fiber filter.

**Note**

(e.g., Advantec type GC50, Whatman type GF/C, or Gelman type A/E)

- 2 Filter sample through 0.2  $\mu\text{m}$  or 0.45  $\mu\text{m}$  low protein binding PVDF filter.

**Note**

(e.g., Millipore Durapore filter)

- 3 Dispense the filtered water samples (e.g., 0.5 L or more) into culture vessels, e.g., 1 L or larger Erlenmeyer flasks.

- 4 Add nutrients to the filtered water to support growth of the target cells.

**Note**

For example, F/2 nutrients for marine cyanobacteria or BG-11 nutrients (Rippka et al. 1979) for freshwater cyanobacteria.

- 5 Seed the filtered water with 1% to 10% v/v of host culture.

**Note**

Target cells must be in exponential growth to avoid loss of potential viral infection (e.g., viruses adsorbed to dead or dying cells will not cause infection and subsequent production of progeny virus).

- 6 As a control, replace the filtered environmental sample with virus-free (0.02  $\mu\text{m}$  filtered or heat-killed) water sample.

**Note**


The volume of this culture need not be as large as the experimental flasks. This control is to make sure that there is not anything in the water sample that would inhibit growth of the target cells (e.g., chemical inhibition).

- 7 Incubate the flasks at the temperature and light conditions appropriate for the cyanobacteria and look for signs of lysis.

**Note**

This could take 2 to 3 weeks, depending on the growth rate of the host as well as the initial titer of cyanophages. It is recommended that the in vivo chlorophyll fluorescence be monitored regularly. A small decrease in relative fluorescence could indicate the presence of a lytic virus.

- 8 Remove an aliquot of the enrichment culture, and pellet remaining cells by centrifugation (e.g., 20 minutes at 6000g).

 00:20:00

- 9 Filter the supernatant through a 0.22  $\mu\text{m}$  or 0.45  $\mu\text{m}$  PVDF filter.

- 10 Store the lysate at 4°C until further analysis.

- 11 Verify the presence of lytic phages by liquid assay (or plaque assay).

- 12 To propagate/amplify the lytic agent, the liquid bioassay is repeated using the putative lytic agent as the test sample.

**Protocol**

NAME

**Propagation/Amplification of Lytic Agent**

CREATED BY

Amy Chan

[Preview](#)

- 12.1 Set up bioassay using 5 mL or larger culture tubes in triplicates.



12.2 Add between 5 to 50  $\mu\text{L}$  of each sample below to target cells in log phase.

Note

- a. Whole lysate (unfiltered).
- b. Filtered lysate (0.22 or 0.45  $\mu\text{m}$ )
- c. Negative control (no addition, or use filtered media)

12.3 Monitor in vivo chlorophyll fluorescence for about 1 week, look for decrease in relative fluorescence compared with control cultures.

12.4 If the cultures lyse, then the lytic agent is most probably a virus.

12.5 Propagate the lytic agent several times to dilute out nonreplicating viruses.

12.6 Filter the lysate and use it to obtain pure clonal stocks.