

Feb 09, 2016

Prophage induction in marine *Synechococcus*

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.dpc5iv>

John H. Paul and Markus Weinbauer¹


¹Manual of Aquatic Viral Ecology

VERVE Net



Verve Team

University of Arizona



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dpc5iv>

External link: http://www.aslo.org/books/mave/MAVE_030.pdf

Protocol Citation: John H. Paul and Markus Weinbauer 2016. Prophage induction in marine *Synechococcus*. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dpc5iv>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: August 25, 2015

Last Modified: November 09, 2017

Protocol Integer ID: 1476

Keywords: prophage induction in marine synechococcus, marine synechococcus, manual of aquatic viral ecology, aquatic viral ecology, detection of lysogeny, lysogeny, marine environment, method for prophage induction, prophage induction

Abstract

This protocol gives a method for prophage induction in marine *Synechococcus*.

Paul, J. H., and M. Weinbauer. 2010. Detection of lysogeny in marine environments, p. 30–33. In S. W. Wilhelm, M. G. Weinbauer, and C. A. Suttle [eds.], *Manual of Aquatic Viral Ecology*. ASLO.

Guidelines

Materials

Sterile 96 well microtiter plates

Indicator host culture (i.e., *Synechococcus* WH7803)

References

Hurley, M. A., and M. E. Roscoe. 1983. Automated statistical analysis of microbial enumeration by dilution series. *J. Appl. Bacteriol.* 55:159–164.

Suttle, C. A., and A. M. Chan. 1994. Dynamics and distribution of cyanophages and their effect on marine *Synechococcus* spp. *Appl. Environ. Microbiol.* 60:3167–3174.

Weinbauer, M. G., and C. A. Suttle. 1996. Potential significance of lysogeny to bacteriophage production and bacterial mortality in coastal waters of the Gulf of Mexico. *Appl. Environ. Microbiol.* 62:4374–4380.

Before start

The samples for prophage induction are pretreated by the technique of viral reduction (Weinbauer and Suttle 1996).

Induction

- 1 Filter each sample through a 0.2- μ m filter to a volume of approximately 5 mL to remove most of the ambient viruses.
- 2 Add virus-free (0.02- μ m filtered) water prepared from the same sample and reduce the volume a second time.
- 3 Return the retentate to its original volume by adding virus-free seawater, dividing into aliquots, and incubating with and without inducing agent.
- 4 Amend treated samples with the inducing agent mitomycin C at a concentration of 1 μ g/mL or with the inducing agent of choice.

Enumeration

- 5 Employ the most probable number (MPN) method to enumerate the cyanophage population (Suttle and Chan 1994).
- 6 Prepare a one- to five-dilution series of the environmental or prophage induction treatment sample using 96-well microtiter plates (Costar, Corning Inc.).
- 7 Freshly dilute a susceptible *Synechococcus* host 1:10 and place in each well.
- 8 Prepare control plates similarly using sterile SN media in the first column of wells.
- 9 Prepare three replicate treatment and control plates from each site.
- 10 Incubate the plates until good growth of the host organism is evident (10–14 days).
- 11 Score wells as positive for virus if lysis of the host organism is evident as a well clearing.
- 12 Calculate viral abundance for each plate using an MPN program (Hurley and Roscoe 1983).



Data analysis

- 13 Evaluate treatment and control cyanophage and *Synechococcus* counts by paired t test between samples using Minitab statistical software.
- 14 Perform comparison of induction results and environmental parameters using linear regression and X^2 analysis, also using Minitab.