

Apr 28, 2018

Ampicillin solution 100 mg/mL

DOI

<https://dx.doi.org/10.17504/protocols.io.puadnse>

Diego Bonatto¹

¹Centro de Biotecnologia da UFRGS



Diego Bonatto

Centro de Biotecnologia da UFRGS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.puadnse>

Protocol Citation: Diego Bonatto 2018. Ampicillin solution 100 mg/mL. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.puadnse>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 28, 2018



Last Modified: April 28, 2018

Protocol Integer ID: 11874

Keywords: Escherichia coli, Ampicillin, Culture medium, ampicillin solution, ml ampicillin solution, ampicillin solution 100 mg, ampicillin, plasmid

Abstract

Ampicillin solution for E. coli strains transformed with plasmid containing AmpR gene

Materials

MATERIALS

⊗ Ethanol 70% [Note: freshly prepared]

⊗ Ampicillin, sodium salt USP **Bio Basic Inc. Catalog #AB0028.SIZE.25g**



Preparation

- 1 Dissolve 1 g of ampicillin in 10 mL of EtOH 70% (v/v). Store solution @ -20 °C.