


May 13, 2020

Cell-free lysate (E. coli) preparation with sonication

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.8nthven>



Nadanai Laohakunakorn¹

¹University of Edinburgh



Nadanai Laohakunakorn

University of Edinburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8nthven>

Protocol Citation: Nadanai Laohakunakorn 2020. Cell-free lysate (E. coli) preparation with sonication. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8nthven>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 23, 2019

Last Modified: May 13, 2020

Protocol Integer ID: 29107

Keywords: cell-free protein synthesis, synthetic biology, in vitro transcription translation, preparation with sonication production, sonication production, free lysate, optional induction of t7 rnap, t7 rnap, optional dialysis, constant energy sonication s12 centrifugation run, cell

Abstract

Production of cell-free lysate from E. coli BL21 Star (DE3) with optional induction of T7 RNAP. Adapted from Kwon and Jewett 2015. Features:

- Variable starting culture sizes (10 mL - 1 L)
- Constant energy sonication
- S12 centrifugation
- Run-off and optional dialysis

Citation

Kwon YC, Jewett MC (2015)

. High-throughput preparation methods of crude extract for robust cell-free protein synthesis. Scientific reports.

<https://doi.org/10.1038/srep08663>

LINK

Protocol successfully used at the University of Edinburgh by Nandanai Laohakunakorn, and at EPFL by Barbora Lavickova and the 2017 iGEM team.



Materials

- Big centrifuge: Eppendorf 5810R with A-4-62 swing-bucket rotor (holds 50 mL tubes)
- Small centrifuge: Eppendorf 5424R with FA-45-24-11 rotor (holds 2 mL tubes)
- Vibra Cell 75186 sonicator / Qsonica Q125 + CL-18 probe
- spectrophotometer
- incubator with shaker
- autoclave
- magnetic stirrer

- LB/2xYTP/2xYTPG medium autoclaved
- E. coli strain of interest

- Tris base (Sigma T1503-100G)
- magnesium glutamate (L-glutamic acid hemimagnesium salt tetrahydrate) (Sigma T1503-100G)
- potassium glutamate (L-glutamic acid potassium salt monohydrate) (Sigma 49601-500G)

- DTT (1,4-dithiothreitol) (Sigma 10708984001)
- acetic acid
- deionized water

- autoclaved tips and petri dishes
- 500 mL Erlenmeyer flasks sterile
- culture tubes sterile
- 50 mL falcon tubes
- 2 mL eppendorf tubes
- 1 L beakers
- 10k MWCO dialysis cassettes (Slide-A-Lyzer, 3mL, Life Technologies)
- magnetic stir bar
- spectrophotometer cuvettes
- liquid nitrogen, dewar, -80 storage



Bacterial growth

1 Prepare all materials for bacterial culture

This protocol will make  200 mL of culture to yield around  1 mL of lysate with a protein concentration of around 40 - 80 mg/mL.

1.1 Reconstitute LB media from premix as per instructions on box.

Note

Protocol has also been successfully carried out with 2x YTP and 2x YTPG (34°C, original). Cells grow slower on YTP media (~5-6h to OD 1.5-1.8) compared to ~4h in LB.

1.2 Autoclave LB media for 00:20:00 at 121 °C along with tips and dishes

2 Grow overnight mini-culture

2.1 Add 5 mL of LB medium to a culture tube and inoculate with small amount of bacteria from glycerol stock.

Note

Strains used successfully with this protocol are

- BL21 (DE3) and BL21 star (DE3) (lacking RNaseE), T7 and E. coli RNAP
- Rosetta (DE3) for oscillators (contains rare tRNA encoding plasmid, good for eukaryotic protein expression)
- Top10, Top10-GamS (GamS inhibits RecBCD degradation of linear DNA), E. coli RNAP, alpha-complementation
- M15 (E. coli RNAP, alpha-complementation)
- DH5alpha (E. coli RNAP, alpha-complementation)

2.2 Grow overnight culture in an incubator at 37 °C and 200 RPM. Make sure cap is in loose position to allow for air exchange

3 Grow lysis culture (the next morning)



- 3.1 **Measure and record** OD600 of overnight culture at **10x dilution**: pipette  900 μ L LB and  100 μ L overnight culture into a cuvette and take absorbance readings at 600 nm.

Expected result

OD600 should be around 4 (0.4 at 10x dilution)

- 3.2 Add  1 mL of overnight culture to  200 mL of LB medium in a 500 mL Erlenmeyer flask.


- 3.3 Incubate the culture for  04:00:00 ,  37 $^{\circ}$ C , and 200 RPM.

4h



Expected result

At the end of the step, the OD₆₀₀ should be ~1.5-2. While Kwon and Jewett grow to OD₆₀₀ = 3, we grow instead for a fixed time. They observe varying robustness of yield to final OD which is strain dependent. **In the end we should probably grow to fixed OD.**



Activity of cells harvested at different densities, grown in 1L 2xYTPG and sonicated at 1.5 mL/556 J. From Kwon and Jewett 2015, Figure 3, CC-BY license.

- 3.4 **(If required)** induce after 2 hours with 0.4 mM IPTG (800 μL of 100 mM stock in 200 mL culture) to express T7 polymerase in the BL21 (DE3) strains.
- 4 A few minutes before incubation step finishes, prepare centrifuge and tubes
 - 4.1 Cool down big centrifuge (Eppendorf 5810R) to 4 °C with fast temp mode; this takes around 10 mins.
 - 4.2 Weigh one 50 mL Falcon tube along with its cap, and record. Clearly label the tube as well as its cap; this tube will be used to spin the final pellet.



- 5 As soon as the incubation finishes, put the Erlenmeyer flask On ice to arrest growth.

Centrifugation and cleaning

- 6 Put 50 mL of culture into each of four 50 mL Falcon tubes (one of which has been labelled), and spin at 4000 rpm , = 3220g 4 °C 00:20:00 . Tubes in swing-bucket rotor A-4-62. 20m
- 7 Put tubes back On ice ; keep bacterial pellet cold as much as possible.
- 8 Carefully discard supernatant using a pipette, then add 10 mL of **Buffer A (no DTT)** to each pellet. Resuspend by carefully pipetting up and down. Finally, transfer all four parts into the single labelled tube.
- 9 Balance centrifuge with a second tube containing 40 mL water. Then spin at 4000 rpm , = 3220g 4 °C 00:10:00 . 10m
- 10 Carefully discard supernatant using a pipette, then add 10 mL of Buffer A to the pellet. Then spin at 4000 rpm , = 3220 g 4 °C 00:10:00 . 10m
- 11 go to step #10 and repeat spin once more.
- 12 Finally, carefully discard as much supernatant as possible.
- 13 Weigh the tube containing the pellet (make sure cap is on), and record:

--	--

Final weight of tube + pellet



The **wet mass** of cells is the difference between this weight and the dry weight measured earlier. This is typically around ~1 g.

- 14 Flash freeze the cells using liquid nitrogen, respecting all safety procedures (wear protective glasses as there is danger of tube explosion etc.). Store the cells at

-80 °C .



Note

Flash-freezing is optional (the protocol will work without) but it serves two purposes:

- Convenient pause point
- Frozen cells, once rethawed, lyse slightly more easily than fresh cells. So this step must be kept consistent to ensure reproducible protocol.

Sonication

- 15 Keeping the pellet On ice , add **Buffer A with DTT**: 1 mL buffer per 1 g **wet mass** as determined earlier. **This ratio is critical.** The DTT must be prepared fresh as it will degrade otherwise.



- 16 Resuspend by vortexing, and put back On ice

- 17 Put 1 mL **exactly** of resuspended cells into a new 2 mL Eppendorf tube. This volume is critical.






Note

The larger the volume the more robust the lysis:

BL21 Star



200 400 600 800
Energy

Yield as a function of volume and sonication energy. From Kwon and Jewett Figure 4, CC-BY license.

So if possible, scale up the production to achieve (ideally) at least 5 mL lysis volumes (= 5 g wet mass).

- 18 Place the tube in an **ice bath** - ice mixed with water - and place the sonicator tip inside the tube so that it is immersed as much as possible but **no part of the tip touches the tube surface**. This is one of the most variable parts of the protocol. A suggested technique is to drop the tip to the bottom of the tube, and then raise it slightly. Holding the tube using a cool rack (e.g. Fisher 15592801) is recommended.
- 19 Sonicate until **400J** has been achieved, using 50% amplitude and pulses of 10s : 10s (i.e. energy for 10s followed by pause for 10s). This typically takes around  00:01:24 but the time is variable. **The total sonication energy is a critical parameter.**



Separation and Clarification

- 20 Cool down the small centrifuge (Eppendorf 5424R) to  4 °C using the fast temp mode: this takes around 20 minutes.
- 21 Balance and centrifuge the lysate at  12000 x g , = 11304 rpm  4 °C  00:10:00
- 22 Transfer the supernatant to a new tube. To prevent any transfer of bacterial debris, it is critical that you **do not take all the lysate**.
- 23 Place the lysate in an incubator at  37 °C , 200 RPM for  01:30:00 , to degrade remaining DNA and RNA. This is called the run-off reaction.

10m



1h 30m



**Note**

Even though Kwon and Jewett observe deterioration in BL21 production yield with run-off, we do not see the same behaviour. If E. coli RNAP will be used, do dialysis as well as run-off.

24 Making sure the centrifuge is still cool, centrifuge the lysate at

12000 x g , = 11304 rpm 4 °C 00:10:00

10m



25 Transfer the supernatant to a new tube. Again, **do not take all the lysate**. Aliquot around

3 µL into separate tube for a Bradford assay.


26 **(IF REQUIRED) carry out dialysis**

26.1 Prepare **Buffer A + DTT** and add 900 mL to a 1 L beaker at 4 °C . Add a magnetic stir bar.

26.2 Rehydrate the required number of 10k MCWO cassettes in **Buffer A+DTT** in a beaker at 4 °C for 00:02:00

26.3 Load cassettes with 2.5 mL max (partial loading possible) of extract and dialyze for 03:00:00 4 °C with stirring.

3h



- 26.4 Finally load extract into centrifuge and centrifuge the lysate at  12000 x g , = 11304 rpm  4 °C  00:10:00
- 26.5 Transfer the supernatant to a new tube. Again, **do not take all the lysate**. Proceed to step 27.
- 27 Keeping everything  On ice , aliquot the lysate into small tubes as required (around  25 µL per tube is recommended).
- 28 Flash-freeze the remaining tubes in liquid nitrogen and store at  -80 °C
- 29 To perform the Bradford assay, dilute  1 µL of the lysate in  99 µL of Buffer A. Mix  5 µL of this diluted lysate with  250 µL of Bradford reagent. Incubate for



00:05:00

before measuring the absorbance using a Nanodrop. The final protein concentration is typically 40 - 80 mg/mL.

Buffer A and B preparation

30 Preparation of 1M stock solutions of tris acetate, magnesium glutamate, and potassium glutamate

- 30.1 Weigh  12.11 g of tris base (121.14 g/mol), add it to  100 mL of deionized water, and adjust to pH 8.2 with acetic acid, to make 1M tris acetate stock.
- 30.2 Weigh  38.86 g of magnesium glutamate (388.61 g/mol), and add it to  100 mL of deionized water, to make 1M magnesium glutamate stock.
- 30.3 Weigh  20.32 g of potassium glutamate (203.23 g/mol), and add it to  100 mL of deionized water, to make 1M potassium glutamate stock.

31 Preparation of 1L of buffer A

- 31.1 Add  10 mL of 1M tris acetate pH 8.2,  14 mL of 1M magnesium glutamate, and  60 mL of 1M potassium glutamate, and fill to  1 L with deionized water to make buffer A (10 mM tris acetate, 14 mM magnesium glutamate, 60 mM potassium glutamate).
- 31.2 Store at  4 °C
- 31.3 Add  2 µL of 1M DTT to  1000 µL Buffer A (final concentration 2 mM) for the lysis solution, or  1.8 mL of 1M DTT to  900 mL of Buffer A for the dialysis solution. This should be made fresh every time but if short-term storage is required, put at  -20 °C



Citations

Kwon YC, Jewett MC. High-throughput preparation methods of crude extract for robust cell-free protein synthesis.

<https://doi.org/10.1038/srep08663>