

Feb 09, 2023

HuBMAP | GE/Vanderbilt MALDI IMS and Cell DIVE™ Modality Overview

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg39k41g25/v1>

Nathan Heath Patterson¹, Elizabeth Neumann², Christine Surrette³, Soumya Ghose³, Liz McDonough³,
Jamie Allen¹, Jeff Spraggins¹, Fiona Ginty³

¹Vanderbilt University; ²UC Davis; ³GE Research

Human BioMolecular Atl...

GE Research



Liz McDonough

GE HealthCare

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg39k41g25/v1>

Protocol Citation: Nathan Heath Patterson, Elizabeth Neumann, Christine Surrette, Soumya Ghose, Liz McDonough, Jamie Allen, Jeff Spraggins, Fiona Ginty 2023. HuBMAP | GE/Vanderbilt MALDI IMS and Cell DIVE™ Modality Overview. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg39k41g25/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 21, 2022

Last Modified: February 09, 2023

Protocol Integer ID: 74344

Keywords: vanderbilt university cell dive collaboration, human biomolecular atlas program, hubmap, vanderbilt maldi im, cell

Funders Acknowledgements:

HuBMAP (Vanderbilt University)

Grant ID: U54DK134302

HuBMAP (GE Research)

Grant ID: 5UH3CA246594-02

Abstract

This is an overview of all protocols currently in use for the GE/Vanderbilt University Cell DIVE collaboration for the Human BioMolecular Atlas Program (HuBMAP). It includes links to each of the individual protocols that make up this project workflow.



MALDI IMS

- 1 Collection of post-surgical tissue.

Protocol

NAME

Collection and Post-Surgical Excision of Human Kidney Tissue through the Cooperative Human Tissue Network

CREATED BY

Jamie Allen

Preview

- 2 Stabilize and freeze tissues.

Protocol

NAME

Freezing Fresh Tissue

CREATED BY

Jamie Allen

Preview

- 3 Cryosection tissues into micrometer thick sections, alternating between thaw mounting onto indium tin-oxide and positively charged glass slides (proceed to step 4), or collecting several tissue sections within an microcentrifuge tube for proteomics analysis.



Protocol

NAME

Cryostat Sectioning of Tissues for 3D Multimodal Molecular Imaging

CREATED BY

Jamie Allen

Preview

MALDI IMS

4 Perform autofluorescence microscopy on all tissue sections

Protocol

NAME

Autofluorescence Microscopy Data Acquisition

CREATED BY

Jamie Allen

Preview

5 Perform Matrix Application

Protocol

NAME

Automatic Deposition of CHCA Matrix for MALDI Analysis of Lipids

CREATED BY

Jamie Allen

Preview



- 6 Perform high resolution IMS analysis of matrix coated tissue sections.

Protocol

NAME

High Resolution Imaging Mass Spectrometry Analysis using Bruker Daltonics Platforms

CREATED BY

Jamie Allen

Preview

MALDI IMS

- 7 Obtain autofluorescence microscopy images of tissues after IMS analysis

Protocol

NAME

Post-IMS Autofluorescence Microscopy

CREATED BY

Jamie Allen

Preview

Preparing Sample for MxIF

- 8 Perform Matrix Removal & Tissue Fixation



Protocol

NAME

Fixation Protocol for Fresh Frozen Tissue Samples (post-MALDI)

CREATED BY

Jamie Allen

Preview

Cell DIVE

- 9 Characterize antibodies (primary/secondary, direct conjugates, and zenon labelled antibodies) and determine any antigen effects from the Cell DIVE dye inactivation process.

[Cell DIVE™ Platform | Antibody Characterization for Multiplexing](#)

[Cell DIVE™ Platform | Antibody Staining & Imaging](#)

- 10 Prepare direct conjugates for study.

[Cell DIVE™ Platform | Antibody Purification Chemistry](#)

[Cell DIVE™ Platform | Ab Conjugation: Initial Conjugation & Scale up Conjugation](#)

- 11 Perform Cell DIVE™ multiplexed data acquisition on the final cohort.

Note

Staining is done manually using a humidity chamber and images are acquired on the Leica Cell DIVE imager utilizing a coverslipless imaging approach