

Jun 21, 2023

LB Broth

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm375nl3p/v1>

sgoodwin^{1,2}

¹University of Sheffield; ²Sheffield iGEM



sgoodwin

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm375nl3p/v1>

Protocol Citation: sgoodwin 2023. LB Broth. protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvm375nl3p/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: June 21, 2023

Last Modified: June 21, 2023



Protocol Integer ID: 83812

Keywords: broth, lb, growth, process

Abstract

Process for producing LB broth growth medium



Making the LB broth

3h 10m

1 Dilute  20 g of LB powder in  800 mL distilled water

5m

2 Stir until fully dissolved with magnetic stirrer

5m

3 Transfer to a Duran flask, and sterilise by autoclaving

3h