



Sep 06, 2022

Lentivirus production

 [Nature Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr4xn3gmk/v1>

Itika Saha¹, F. Ulrich Hartl¹, Mark S. Hipp^{2,3}

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry, Am Klopferspitz 18, 82152 Martinsried, Germany;

²Department of Biomedical Sciences of Cells and Systems, University Medical Center Groningen, University of Groningen, Antonius Deusinglaan, 1, 9713 AV Groningen, The Netherlands;

³School of Medicine and Health Sciences, Carl von Ossietzky University Oldenburg, 26129 Oldenburg, Germany



Felix Kraus

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr4xn3gmk/v1>

External link: <https://www.biorxiv.org/content/10.1101/2022.02.18.481043v1.full>



Protocol Citation: Itika Saha, F. Ulrich Hartl, Mark S. Hipp 2022. Lentivirus production. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr4xn3gmk/v1>

Manuscript citation:

The AAA+ chaperone VCP disaggregates Tau fibrils and generates aggregate seeds Itika Saha, Patricia Yuste-Checa, Miguel Da Silva Padilha, Qiang Guo, Roman Körner, Hauke Holthusen, Victoria A. Trinkaus, Irina Dudanova, Rubén Fernández-Busnadiego, Wolfgang Baumeister, David W. Sanders, Saurabh Gautam, Marc I. Diamond, F. Ulrich Hartl, Mark S. Hipp bioRxiv 2022.02.18.481043; doi: <https://doi.org/10.1101/2022.02.18.481043>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 05, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 69582

Keywords: ASAPCRN, lentivirus production, production of lentivirus, lentivirus production this protocol, lentivirus

Abstract

This protocol describes the production of lentiviruses to transduce HEK293T cells and has to be performed in a biosafety level 2 laboratory

Safety warnings

! Has to be performed in a biosafety level 2 laboratory.



Lentivirus production

1d 1h 10m

- 1 Plate $\sim 3.6 \times 10^6$ Lenti-X HEK293T cells (Takara) in 10 cm dish in 10 mL standard DMEM. Cells should be $\sim 80\%$ confluent at the time of transfection.

Note

NOTE: Only low passage cells should be used.

- 2 Next day, remove 5 mL medium and replenish with fresh medium.
- 3 Warm up reduced serum medium e.g. Opti-MEM (Gibco) and transfection reagent to room temperature (RT). This protocol was performed with Lipofectamine 3000 transfection reagent (Thermo).
- 4 Add 24 μL Lipofectamine 3000 to 600 μL Opti-MEM, mix by vortexing and incubate 5 min at RT. 5m
- 5 In another tube, mix 6 μg plasmid containing gene of interest, 5 μg packaging plasmid psPAX2 (RRID:Addgene_12260), 1 μg envelope plasmid pMD2.G (RRID:Addgene_12259) and 24 μL P3000 reagent (provided by manufacturer along with Lipofectamine 3000 reagent) in 600 μL Opti-MEM, mix by vortexing and incubate 5 min at RT.
- 6 Mix contents of both tubes and incubate for 15 min at RT. 15m
- 7 Add DNA-lipid complex to cells dropwise.
- 8 2 days later, collect virus-containing medium and centrifuge for 5 min at 1,000 x g. 2d
- 9 Collect supernatant in a fresh tube and proceed with concentration.

Concentration

1d 1h 10m



- 10 Add Lenti-X concentrator (Takara) to clarified virus-containing medium at 1:4 dilution and mix well by gently inverting tube.
- 11 Incubate overnight or 2 h at 4 °C.
- 12 Next day, centrifuge for 45 min at 1,500 x g at 4 °C followed by gently aspirating supernatant.
- 13 Resuspend viral pellet in 100 - 1000 µL PBS, aliquot and store at -80 °C until use.

1d

45m