

Jan 12, 2026

Nuclease-MDA for Enriching Monkeypox Virus Genomes for NGS

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l22524l1y/v1>

Masayasu Misu¹, Takeshi Kurosu¹, Tomoki Yoshikawa¹, Madoka Kawahara¹, Kohei Oishi¹, Masayuki Shimojima¹, Hideki Ebihara¹

¹Department of Virology I, National Institute of Infectious Diseases



Tomoki Yoshikawa

National Institute of Infectious Diseases

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l22524l1y/v1>

Protocol Citation: Masayasu Misu, Takeshi Kurosu, Tomoki Yoshikawa, Madoka Kawahara, Kohei Oishi, Masayuki Shimojima, Hideki Ebihara 2026. Nuclease-MDA for Enriching Monkeypox Virus Genomes for NGS. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4r3l22524l1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 20, 2023

Last Modified: January 12, 2026

Protocol Integer ID: 83667

Keywords: Oxford Nanopore Technology, RNA virus, Sequence method, MinION, Nanopore sequencing, RCA-NGS, mda for enriching monkeypox virus, enrichment of the monkeypox virus, enriching monkeypox virus, viral genome, genome sequencing of mpvx, monkeypox virus, intact viral genome, idt for illumina dna, unwanted host dna via micrococcal nuclease treatment, illumina dna prep, high proportion of viral read, illumina dna, genome sequencing, viral read, genome amplification, cowpox virus, virus concentration step, micrococcal nuclease treatment, unwanted host dna, other poxvirus, rna, employing phi29 polymerase, genome, phi29 polymerase, nuclease, sequencing, oxford nanopore technology, vaccinia virus, iseq100 system

Funders Acknowledgements:

Ministry of Health, Labor and Welfare of Japan

Grant ID: 20HA2005 and 23HA2004

Japan Agency for Medical Research and Development

Grant ID: 21fk0108426j0001, 25fk0108660h0003, 22fk0108502j0101 and 23fk0108581j0701

Japan Society for the Promotion of Science

Grant ID: JP22K085

Abstract

This protocol describes the enrichment of the monkeypox virus (MPXV) genome using the Nuclease-MDA method, optimized for Oxford Nanopore Technologies (ONT) sequencing with V14 chemistry.

This protocol can be applied using the ONT Native Barcoding Kit(SQK-NBD114.24, SQK-NBD114.96), Rapid Barcoding Kit V14 (SQK-RBK114.24 or SQK-RBK114.96) or Illumina DNA Prep, (M) Tagmentation (24 Samples, IPB) (Cat#20060060), IDT for Illumina DNA/RNA UD Indexes Set A (Cat#20027213), and the iSeq100 system.

This approach eliminates unwanted host DNA via micrococcal nuclease treatment, preserving the intact viral genome. The viral genome is then amplified using multiple displacement amplification (MDA), an isothermal whole-genome amplification (WGA) technique employing Phi29 polymerase. This method provides a high proportion of viral reads (~90%) without requiring specific PCR or virus concentration steps.

This method enables high-depth whole-genome sequencing of MPXV and can be readily adapted to other poxviruses, such as vaccinia virus and cowpox virus, with appropriate primer modifications.



Materials

-  Micrococcal Nuclease - 320,000 gel units **New England Biolabs Catalog #M0247S**
-  High Pure Viral Nucleic Acid Kit **Roche Catalog #11858874001**
- **Dr.GentLE Precipitation Carrier - Takara Catalog #9094**
-  Agencourt AMPure XP **Beckman Coulter Catalog #A63880**
- **GenomiPhi V3 Ready-To-Go DNA Amplification Kit - Cytiva Catalog #25-6601-24**
-  T7 Endonuclease I - 250 units **New England Biolabs Catalog #M0302S**
-  NEBNext FFPE DNA Repair Mix - 24 rxns **New England Biolabs Catalog #M6630S**
-  NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns **New England Biolabs Catalog #E7546S**
-  Blunt/TA Ligase Master Mix - 50 rxns **New England Biolabs Catalog #M0367S**
-  NEBNext Quick Ligation Module - 20 rxns **New England Biolabs Catalog #E6056S**
- **Native Barcoding Kit 24 V14 - Oxford Nanopore Technologies Catalog #SQK-NBD114.24**
-  Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
-  Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230**
-  DNA LoBind Tube 1.5ml **Eppendorf Catalog #022431021**
-  0.2 ml PCR Tube strips **Eppendorf Catalog #0030124359**
- **100 % ethanol**
- **70 % ethanol**
- **TE(pH8.0)**
- **nuclease-free H₂O**



Protocol materials

- ⊗ Micrococcal Nuclease - 320,000 gel units **New England Biolabs Catalog #M0247S**
- ⊗ T4 RNA Ligase 2, truncated KQ - 2,000 units **New England Biolabs Catalog #M0373S**
- ⊗ Superscript IV **Thermo Fisher Scientific Catalog #18090050**
- ⊗ 0.2 ml PCR Tube strips **Eppendorf Catalog #0030124359**
- ⊗ NEBNext FFPE DNA Repair Mix - 24 rxns **New England Biolabs Catalog #M6630S**
- ⊗ NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns **New England Biolabs Catalog #E7546S**
- ⊗ Agencourt AMPure XP **Beckman Coulter Catalog #A63880**
- ⊗ High Pure Viral Nucleic Acid Kit **Roche Catalog #11858874001**
- ⊗ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ⊗ Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230**
- ⊗ DNA LoBind Tube 1.5ml **Eppendorf Catalog #022431021**
- ⊗ T7 Endonuclease I - 250 units **New England Biolabs Catalog #M0302S**
- ⊗ Blunt/TA Ligase Master Mix - 50 rxns **New England Biolabs Catalog #M0367S**
- ⊗ NEBNext Quick Ligation Module - 20 rxns **New England Biolabs Catalog #E6056S**
- ⊗ Qubit dsDNA BR (Broad Range) assay **Thermo Fisher Scientific Catalog #Q32850**

Safety warnings

- ! Follow your facility's regulations and biosafety practices.

Before start

It is recommended to check for the absence of bacterial contamination (e.g., *Mycoplasma* spp.).



Preparation for virus supernatant

- 1 Centrifuge the working stock virus to remove debris.

6000 x g, Room temperature, 00:10:00

10m

- 2 Transfer 180 μ L virus supernatant to a 1.5ml screw cap tube.

- 3 Unwanted DNA and RNA mainly originating from the virus-infected cells are digested using



Micrococcal Nuclease - 320,000 gel units **New England Biolabs Catalog #M0247S**

- 3.1 Total 200 μ l reaction

1h

- 179 μ L virus supernatant
- 20 μ L 10X Micrococcal Nuclease Reaction Buffer
- 1 μ L Micrococcal nuclease

Mix by pipetting and spin down.

37 °C water bath 01:00:00

The viral genomic DNA extraction

- 4 The viral genomic DNA extraction is performed using



High Pure Viral Nucleic Acid Kit **Roche Catalog #11858874001**

- 4.1 Add 200 μ L of binding buffer, 4 μ L PolyA carrier RNA and 50 μ L Proteinase K.

10m

Mix well by pipetting and inverting the tube thoroughly, and spin down.

72 °C 00:10:00

- 4.2 Add 100 μ L of binding buffer and mix well by pipetting or inverting the tube thoroughly, and spin down.



4.3 Transfer the whole sample to a High Pure Filter Tube.

8000 x g, Room temperature, 00:01:00

1m

Discard the flow-through liquid and Collection Tube, and insert the Filter Tube into a new Collection Tube.

4.4 Add 500 μL of inhibitor removal buffer.

1m

8000 x g, Room temperature, 00:01:00

Discard the flow-through liquid and Collection Tube, and insert the Filter Tube into a new Collection Tube.

4.5 Add 450 μL of wash buffer.

1m

8000 x g, Room temperature, 00:01:00

Discard the flow-through liquid and Collection Tube, and insert the Filter Tube into a new Collection Tube.

4.6 Add 450 μL of wash buffer.

1m

13000 x g, Room temperature, 00:01:00 and discard the flow-through liquid.

Discard the Collection Tube and insert the Filter Tube into a 1.5 ml tube -

DNA LoBind Tube 1.5ml **Eppendorf Catalog #022431021** .

4.7 Add 50 μL elution Buffer.

1m

13000 x g, Room temperature, 00:01:00

Note

The purified 50 μL viral genomic DNA can be stored at -80°C .

Amplification of DNA by multiple displacement amplification (MDA)

5 DNA is amplified by MDA using **GenomiPhi V3 Ready-To-Go DNA Amplification Kit - Cytiva Catalog #25-6601-24**.



5.1 Total 20 µL reaction

3m

- 10 µL purified viral genomic DNA
- 10 µL 2X denaturation buffer

Mix by pipetting and spin down.

95 °C 00:03:00

4 °C on ice

5.2 Add 20 µL denatured sample to Ready to go GenomiPhi cake.

4h 10m

30 °C 04:00:00

65 °C 00:10:00

Note

The amplified 20 µL viral genomic DNA can be stored at -80°C.

DNA purification by AMPure XP

6 The DNA is purified using

Agencourt AMPure XP Beckman Coulter Catalog #A63880

6.1 Add 10 µL (X0.5 volume) AMPure beads

6.2 Incubate on rotor mixer.

5m

00:05:00 Room temperature

6.3 Spin down and pellet on a magnet.

1m

Wait for 00:01:00 and pipette off the supernatant.

6.4 Wash twice by 100 µL 70 % ethanol, remove the ethanol using a pipette, and discard.



6.5 Spin down and pipette off any residual ethanol.

6.6 Resuspend pellet in  40 μL nuclease-free H_2O .

5m

Incubate on a rotor mixer.

 00:05:00 Room temperature

6.7 Spin down and pellet the beads on the magnet until the elute is clear and colorless.

6.8 Remove and retain  40 μL elute into a new tube.

7 DNA concentration is measured using a Qubit 4 Fluorometer with

2m

 Qubit dsDNA BR (Broad Range) assay **Thermo Fisher Scientific Catalog #Q32850**

-  199 μL 1X working solution

-  1 μL DNA

Mix by vortexing.

Incubate  00:02:00  Room temperature and measure.

Note

Confirm the total amplified DNA to be over 1500 ng.

The DNA can be stored at  -30 $^{\circ}\text{C}$ for at least two months.

**Note**

At this step, whole-genome sequencing using an Illumina next-generation sequencer is also possible. In this case, the T7 endonuclease treatment and any subsequent steps in this protocol are not required. Instead, you should follow Illumina's library preparation protocol. We have confirmed that this works using the Illumina DNA Prep, (M) Tagmentation (Cat#20060060) and the IDT for Illumina DNA/RNA UD Indexes Set A (Cat#20027213).

Similarly, the Oxford Nanopore Rapid Barcoding Kit 24 V14 (SQK-RBK114.24) can also be used. Here, the T7 endonuclease treatment is also not required. From this point, you should proceed by following the Rapid Barcoding Kit protocol, which significantly reduces library preparation time.

T7 endonuclease treatment

- 8 The amplified DNA by MDA is digested using

 T7 Endonuclease I - 250 units **New England Biolabs Catalog #M0302S** .

The following protocol is modified based on the Native barcoding amplicons (with EXP-NBD104, EXPNBD114, and SQK-LSK109) protocol (NBA_9093_v109_revA_12Nov2019) provided by the Oxford Nanopore Technologies website.

- 8.1 Total 30 µl reaction

30m

-  x µL (1.0 µg) DNA
-  3 µL NEBuffer 2
-  1.5 µL T7 endonuclease I
-  25.5-x µL nuclease-free H₂O

Mix by pipetting and spin down.

 37 °C  00:30:00

- 9 The DNA is purified using

 Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

Add  15 µL (X0.5 volume) AMPure beads

 [go to step #6.3](#)

Resuspend pellet in  13 µL nuclease-free H₂O.

DNA repair and end-prep



10 The purified DNA is end-prepped using

NEBNext FFPE DNA Repair Mix - 24 rxns **New England Biolabs Catalog #M6630S**

and

NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns **New England Biolabs Catalog #E7546S**

10.1 Total 15 µl reaction

35m

- 12 µL DNA
- 0.875 µL NEB Next FFPE DNA repair buffer
- 0.5 µL NEB Next FFPE DNA repair Mix
- 0.875 µL Ultra II end-prep reaction buffer
- 0.75 µL Ultra II end-prep reaction Mix

Mix by pipetting and spin down.

20 °C 00:30:00

65 °C 00:05:00

11 The DNA is purified using

Agencourt AMPure XP **Beckman Coulter Catalog #A63880** .

Add 15 µL (X1.0 volume) AMPure beads



Resuspend pellet in 10 µL nuclease-free H₂O.

Note

The DNA can be stored at 4°C overnight.

Native barcode ligation

12 The end-prepped DNA is ligated with native barcode using **Native Barcoding Kit V14 - Oxford Nanopore Technologies Catalog #SQK-NBD114.24** or **#SQK-NBD114.24** with

Blunt/TA Ligase Master Mix - 50 rxns **New England Biolabs Catalog #M0367S** .

12.1 Total 20 µl reaction

20m

- 8.5 µL DNA

- 1.5 µL native barcode
- 10 µL Blunt/TA ligase master mix

Mix by pipetting and spin down.

Room temperature 00:20:00

- 12.2 Add the following volume of EDTA to each well and mix thoroughly by pipetting and spin down briefly.

Note

Ensure you follow the instructions for the cap colour of your EDTA tube.

	EDTA cap colour	Volume per well
	For clear cap EDTA	2 µl
	For blue cap EDTA	4 µ

EDTA is added at this step to stop the reaction.

- 13 Pool all the barcoded samples in a 1.5 ml Eppendorf DNA LoBind tube.

- 14 The DNA is purified using

Agencourt AMPure XP Beckman Coulter Catalog #A63880

Add 8.8 or 9.6 µL (X0.4 volume) AMPure beads multiply by the original number of tubes.

- 14.1 Incubate on rotor mixer.

00:05:00 Room temperature

5m

- 14.2 Spin down and pellet on a magnet.

Wait for 00:01:00 and pipette off the supernatant.

1m



14.3 Wash twice by  700 µL 70 % ethanol, remove the ethanol using a pipette, and discard.

14.4 Spin down and pipette off any residual ethanol.

14.5 Resuspend the pellet in 15 µL of nuclease-free H₂O.

If the pooled sample volume exceeds 100 µL (i.e., more than 5 samples), resuspend the pellet in x0.2 volume of the pooled sample (e.g., 20 µL when 5 samples are pooled) of nuclease-free H₂O.

Incubate on a rotor mixer.

 00:10:00

 37 °C

Every 2 minutes, agitate the sample by gently flicking for 10 seconds to encourage DNA elution.

14.6 Spin down and pellet the beads on the magnet until the elute is clear and colorless.

14.7 Remove and retain the elute into a new tube.

10m

Adaptor ligation and clean-up

20m

15 Adaptor Ligation with pooled samples is performed using **Ligation Sequencing Kit - Oxford Nanopore Technologies Catalog #SQK-NBD114.24 or #SQK-NBD114.96 with**

 NEBNext Quick Ligation Module - 20 rxns **New England Biolabs Catalog #E6056S**

15.1 Total 20 µl reaction

20m

-  12 µL DNA
-  2 µL Native Adapter (NA)
-  4 µL NEB Next Quick Ligation Reaction Buffer(5X)
-  2 µL Quick T4 DNA ligase

Mix gently and incubate.

 Room temperature

 00:20:00

16 The adaptor-ligated DNA is purified using

 Agencourt AMPure XP **Beckman Coulter Catalog #A63880** .



Prepare AMPure XP beads for use; resuspend by vortexing.
Transfer amplified DNA sample to 1.5ml low binding tube.

- 16.1 Add  10 μ L (X0.5 volume) AMPure XP reagent and mix by pipetting.

Incubate on a rotor mixer.

 00:05:00

 Room temperature

- 16.2 Spin down and pellet on a magnet. Wait for  00:01:00 and pipette off the supernatant.

- 16.3
 - Wash twice by  100 μ L Short Fragment Buffer(SFB) and remove the SFB using a pipette and discard.

- 16.4 Spin down and pipette off any residual SFB.

- 16.5
 - Resuspend pellet in  15 μ L Elution Buffer (EB)

5m

 00:05:00

 Room temperature and tapping occasionally.

- 16.6 Spin down and pellet the beads on the magnet until the elute is clear and colourless.

- 16.7 Remove retain  15 μ L elute into a new tube.

- 17 DNA concentration is measured using a Qubit 4 Fluorometer with

2m

 Qubit dsDNA BR (Broad Range) assay **Thermo Fisher Scientific Catalog #Q32850**

.

-  199 μ L 1X working solution

-  1 μ L DNA

Mix by vortexing.

Incubate  00:02:00  Room temperature and measure.



Note

The molar concentration of the DNA sample can be converted based on the length of the major band confirmed by electrophoresis after T7 endonuclease treatment. Typically, the fragment lengths are around 2000 base pairs.

- 18 Make up the library to  12 μL at 10-20 fmol

Sequencing by MinION

- 19 Sequencing according to the manufacturer's instructions.