

Oct 06, 2022

## Obtain LB broth medium

DOI

[dx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1)

An.Huang<sup>1</sup>

<sup>1</sup>XJTLU



An.Huang

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1>

**Protocol Citation:** An.Huang 2022. Obtain LB broth medium. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 05, 2022



**Last Modified:** October 06, 2022

**Protocol Integer ID:** 70864

**Keywords:** ordinary lb broth medium, lb broth medium lb broth medium, lb broth medium, basic culture media for microorganism, microorganism, basic culture media

## Disclaimer

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK



The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

LB broth medium is a basic culture media for microorganisms. This protocol helps make an ordinary LB broth medium.

## Troubleshooting



1 Dissolve  25 g LB broth powder in  1000 mL distilled water.

2 Autoclave the medium.

Safety information

Use autoclave safely.

3 Store the medium at  4 °C .