

Oct 24, 2025

🌐 Protein Extraction Protocol from the Salivary Glands of Hemipterans: Application in *Mahanarva spectabilis*

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnrr9ml5d/v1>



Angelo José Rinaldi¹, Monique da Silva Bonjour¹, Humberto Josué de Oliveira Ramos¹, Alexander Machado Auad²

¹UFV; ²Embrapa

Bioline

Metabolomics Protocols ...

1 more workspace



Angelo José Rinaldi

UFV

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnrr9ml5d/v1>

Protocol Citation: Angelo José Rinaldi, Monique da Silva Bonjour, Humberto Josué de Oliveira Ramos, Alexander Machado Auad 2025. Protein Extraction Protocol from the Salivary Glands of Hemipterans: Application in *Mahanarva spectabilis*.

protocols.io <https://dx.doi.org/10.17504/protocols.io.14egnrr9ml5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 23, 2025

Last Modified: October 24, 2025

Protocol Integer ID: 230603

Keywords: salivary glands, protein extraction, proteomics, TCA, acetone, urea, CHAPS, PMSF, *Mahanarva spectabilis*, Hemiptera, protein extraction protocol from the salivary gland, protein extraction process from the salivary gland, protein extraction process, extraction of insoluble protein, protein extraction protocol, extraction of soluble protein, resolubilization for proteomic analysis, proteomic analysis, extraction, salivary gland, insoluble protein, protein precipitation with tca, protein precipitation, soluble protein, spittlebug mahanarva, hemipteran, hemiptera, application in mahanarva

Abstract

This protocol describes the collection and protein extraction process from the salivary glands of the spittlebug *Mahanarva spectabilis* (Hemiptera: Cercopidae). The method combines dissection under a stereomicroscope, mechanical and chemical lysis, centrifugation, protein precipitation with TCA/acetone, and resolubilization for proteomic analysis. It includes two workflows:

- **Workflow A** – Extraction of soluble proteins (supernatant)
- **Workflow B** – Extraction of insoluble proteins (pellet)

Guidelines

- Always perform the entire procedure on **ice (0–4 °C)** or in a **refrigerated environment** to prevent protein degradation.
- Prepare all **buffers and reagents fresh** or store them no longer than one week at 4 °C (without PMSF).
- Add **protease inhibitors (PMSF or cocktail)** immediately before use to ensure full protection of proteins.
- **Label all tubes and reagents clearly** to avoid cross-contamination between samples.
- Calibrate and pre-cool all instruments (centrifuge, sonicator, etc.) prior to use.
- Keep the work area clean and organized — dedicate one ice tray and one set of tools per sample batch.
- Use **low-binding plasticware** (e.g., microtubes) for protein samples to minimize sample loss.
- For reproducibility, note the **number of glands** processed per sample and the **approximate weight** of tissue.
- Check pH of buffers before use; ensure that Tris and PBS buffers are within ± 0.1 pH unit of the target value.
- Perform all protein-handling steps using **sterile, chilled pipette tips** to avoid contamination.



Materials

- Stereomicroscope
- Fine forceps and scalpel
- Microtubes (1.5 mL, low protein-binding) and Falcon tubes (15/50 mL)
- Manual or mechanical tissue homogenizer
- Probe sonicator (with amplitude control)
- Refrigerated centrifuge
- Orbital shaker
- Ice bath and water bath (37 °C)
- Pipettes and sterile tips
- –80 °C ultrafreezer
- Spectrophotometer (Bradford or BCA assay)
- SDS-PAGE or LC-MS/MS system

Troubleshooting

Problem

Low protein yield, possible cause: PMSF degradation, insufficient tissue

Solution

Prepare PMSF fresh; increase number of glands

Problem

Protein degradation, possible cause: Delay or poor cooling

Solution

Keep on ice at all times; add PMSF promptly

Problem

Pellet difficult to dissolve, possible cause: Over-drying

Solution

Gradually rehydrate or increase buffer volume

Problem

Assay interference, possible cause: Detergent or urea incompatibility

Solution

Dilute sample or use compatible assay kits

Safety warnings

- ❗ **Trichloroacetic acid (TCA)** is *corrosive* — handle only under a **fume hood**, wearing gloves and eye protection.
- **Acetone** is *highly flammable* — keep away from open flames, sparks, and heating elements.
- **Phenylmethylsulfonyl fluoride (PMSF)** is *toxic and unstable in aqueous solution* — prepare fresh in **isopropanol** immediately before use, and dispose of properly after use.
- **Urea** decomposes slowly, producing cyanate — do not store buffers containing urea for more than 1 week.
- Use **face protection and hearing protection** during **probe sonication** to prevent aerosol formation and hearing damage.
- Always wear **gloves, lab coat, and protective eyewear** when handling biological material and reagents.
- Dispose of chemical and biological waste following your institution's **biosafety and hazardous waste protocols**.

Before start

Work at **4 °C** or on **ice** at all times.

- Use complete **personal protective equipment (PPE)**: lab coat, nitrile gloves, safety goggles, face shield, and ear protection for sonication.
- **Prepare all buffers and solutions fresh** or store them at 4 °C for no more than one week (without PMSF).
- Add **protease inhibitors (PMSF)** immediately before use.
- **Pre-chill** acetone and TCA/acetone solution to −20 °C for at least 1 h before use.
- Ensure the centrifuge is **pre-cooled to 4 °C**.

Reagents and Preparation

1 Acidified Water (pH $\approx$ 3.0)

Use: Initial storage of dissected glands.

Preparation (100 mL):

- Add **8.3 μ L of concentrated HCl (12 M)** to **100 mL Milli-Q water**.
- Mix well and check that pH $\approx$ 3.0.
- Store at 4 °C for up to one week.

2 Protein Extraction Buffer (Workflow A)

Final composition:

- 50 mM Tris-HCl (pH 7.5)
- 150 mM NaCl
- 1 mM EDTA
- 1 % Triton X-100
- 1 mM PMSF (add freshly before use)

Preparation (100 mL):

1. Dissolve **0.605 g Tris base** in 80 mL Milli-Q water.
2. Add **0.876 g NaCl** and **0.037 g EDTA·Na₂**.
3. Add **1 mL Triton X-100** (use wide-bore tip).
4. Adjust pH to **7.5** with 1 M HCl.
5. Bring to 100 mL with Milli-Q water.
6. Filter-sterilize (0.22 μ m).
7. Immediately before use, add **PMSF to 1 mM final** from a 100 mM isopropanol stock.
8. Keep on ice during use.

3 10 % (w/v) TCA in Acetone (Chilled)

Use: Protein precipitation.

Preparation (100 mL):

1. Dissolve **10 g trichloroacetic acid (TCA)** in **100 mL anhydrous acetone**.
2. Mix until fully dissolved (use fume hood).
3. Store at **-20 °C** in an amber bottle for up to one month.

4 100 % Acetone (Chilled)

Use: Washing of protein pellets.

- Keep acetone at **-20 °C** prior to use.
- Always use chilled acetone for washes.

5 Solubilization Buffer (Workflow B)

Final composition:

- 8 M urea
- 2 % CHAPS
- 50 mM Tris-HCl (pH 8.0)

- 1 mM EDTA
- 1 mM PMSF (fresh)

Preparation (100 mL):

1. Dissolve **48 g urea** in 60 mL Milli-Q water (gentle stirring).
2. Add **2 g CHAPS**, mix until dissolved.
3. Add **0.605 g Tris base** and **0.037 g EDTA·Na₂**.
4. Adjust pH to **8.0** with 1 M HCl.
5. Bring to 100 mL final volume.
6. Store at 4 °C for up to 7 days.
7. Add **PMSF to 1 mM** immediately before use.

6 PBS 1× (Phosphate-Buffered Saline)

Use: Washing step for pellets.

Final composition:

- 137 mM NaCl
- 2.7 mM KCl
- 10 mM Na₂HPO₄
- 1.8 mM KH₂PO₄
- pH 7.4

Preparation (1 L):

1. Dissolve **8.0 g NaCl**, **0.2 g KCl**, **1.44 g Na₂HPO₄**, and **0.24 g KH₂PO₄** in 800 mL Milli-Q water.
2. Adjust pH to **7.4** with 1 M HCl.
3. Bring to 1 L with water.
4. Store at 4 °C for up to one month.

7 Safety Notes

- **TCA** is **corrosive** → handle in fume hood with gloves and goggles.
- **Acetone** is **flammable** → keep away from ignition sources.
- **PMSF** is **toxic and unstable in water** → prepare only under a fume hood and discard properly.
- **Urea** decomposes to cyanate → prepare fresh when possible.
- Dispose of all waste according to institutional chemical-safety regulations.

**Workflow A — Soluble Protein Fraction (Supernatant)****8 Workflow A — Soluble Protein Fraction (Supernatant) Gland Collection**

- Dissect salivary glands from adult *M. spectabilis* under a stereomicroscope.
- Place immediately into microtubes containing acidified water (pH 3).

Keep on ice or store at −80 °C until extraction.



- 9 **Buffer Addition**
 - Transfer glands into 1.5 mL microcentrifuge tubes.
 - Add **500 µL of extraction buffer**. Keep on ice.
- 10 **Homogenization and Lysis**
 - Homogenize manually or mechanically until uniform suspension forms.
 - Sonicate (probe): **5 cycles of 30–60 s at 20–30 % amplitude**, with **1 min on ice between cycles**.
- 11 **Cold Incubation**
 - Incubate the homogenate on ice for **30 min** for optimal solubilization.
- 12 **Centrifugation (Clarification)** 
 - Spin at **12 000 × g for 10 min at 4 °C**.
 - Transfer the supernatant (soluble proteins) to a new tube.
 - Store pellet at **–80 °C** for Workflow B.
- 13 **Protein Precipitation**
 - Add **an equal volume** of chilled **10 % TCA/acetone (1:1)**.
 - Mix gently and incubate on ice for **30 min**.
- 14 **Post-Precipitation Centrifugation**
 - Centrifuge at **12 000 × g for 10 min at 4 °C**.
 - Carefully discard the supernatant.
- 15 **Pellet Washing**
 - Wash pellet **three times** with **1 mL chilled acetone (100 %)**.
 - Each wash: mix gently, spin 12 000 × g for 5 min at 4 °C.
 - Air-dry the pellet for 3–5 min (do not overdry).
- 16 **Resuspension and Storage**
 - Resuspend pellet in a suitable buffer (e.g., **8 M urea + 2 % CHAPS**).
 - Vortex or pipette until fully dissolved.
 - Store at **–20 °C** or **–80 °C**.
 - Quantify protein concentration (Bradford or BCA) before SDS-PAGE or MS analysis.

Workflow B — Insoluble Protein Fraction (Pellet)

- 17 **Pellet Washing**
 - Add **1 mL cold PBS 1×**, vortex gently, and spin **12 000 × g for 5 min at 4 °C**.
 - Discard supernatant.
- 18 **Pellet Solubilization**
 - Add **500 µL solubilization buffer** (8 M urea + 2 % CHAPS + PMSF 1 mM).
 - Homogenize by pipetting or gentle vortexing.

**19 Incubation for Solubilization**

- Shake at room temperature for **1–2 h** or incubate at **37 °C for 30 min** to assist solubilization.

20 Centrifugation

- Spin **12 000 × g for 10 min at 4 °C**.
- Transfer the supernatant (solubilized proteins) to a new tube.

21 Preparation for Analysis

- Quantify protein concentration (Bradford or BCA).
- For SDS-PAGE: add sample buffer containing SDS + β -mercaptoethanol and heat to **95 °C for 5 min**.
- For proteomics: perform enzymatic digestion (e.g., trypsin) followed by LC-MS/MS.

22 Storage

- Store samples at **–20 °C or –80 °C** in aliquots to avoid freeze–thaw cycles.

Protocol references

- Lomate & Bonning (2016). *Distinct properties of proteases and nucleases in the gut, salivary gland and saliva of Nezara viridula*. *Scientific Reports* 6:27587.
- Ramm et al. (2015). *Morphology and Proteome Characterization of the Salivary Glands of Halyomorpha halys*. *Journal of Insect Physiology* 76:29–37.
- Alhaddad et al. (2011). *Comparative Behavioral and Protein Study of Salivary Glands of an Insect*. *Annals of the Entomological Society of America* 104(3):543–551.
- He et al. (2024). *Comparative Transcriptomic Analysis of Salivary Glands of Hemiptera*. *BMC Genomics* 25(1):87.
- Ranganathan et al. (2019). *Insight into Salivary Gland Proteomes of Two Polyphagous Hemipteran Bugs*. *Proteomics* 19(10):e1800436.

Acknowledgements

Fundação de Amparo à Pesquisa do Estado de Minas Gerais (FAPEMIG)
Empresa Brasileira de Pesquisa Agropecuária (Embrapa) Gado de Leite JF/MG
Universidade Federal de Viçosa (UFV)