

Jul 10, 2024

Protein interaction modeling

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8215rl2w/v1>

Shiyi Wang¹

¹Duke University



Shiyi Wang

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8215rl2w/v1>

Protocol Citation: Shiyi Wang 2024. Protein interaction modeling. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8215rl2w/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 10, 2024

Last Modified: July 10, 2024

Protocol Integer ID: 103176

Keywords: ASAPCRN, protein interaction modeling, modeling protein interaction modeling, using alphafold, multimer, protein

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-020607

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Protein interaction modeling using AlphaFold 2.0 Multimer

- 1 ****Structure Prediction in AlphaFold 2.0 Multimer**** - Use full-length mouse Ezrin (UniProt ID: Q4KML7) and Atg7 (UniProt ID: Q9D906) for structure prediction. - Predict the structure of Ezrin in the open conformation: - Independently predict Ezrin N- and C-termini. - Thread them onto an open ERM hinge structure using PyMOL2 molecular visualization software.
- 2 ****Modeling Ezrin:Atg7 Interaction**** - Predict the structures of Ezrin and the Atg7 homodimer using AlphaFold 2.0 Multimer. - Dock the predicted structures to model the Ezrin:Atg7 interaction and associated conformational changes.
- 3 ****Assess Conformational Changes**** - Use CABS-flex 2.0 (REF: <https://doi.org/10.1093/nar/gky356>) and PyMOL2 to assess conformational changes in the Ezrin:Atg7 protein complex.
- 4 ****Select High-Confidence Structures**** - For all predicted structure models, select the highest-confidence structures based on the predicted local distance difference test (pLDDT).
- 5 ****Energy Minimization**** - Subject the highest-confidence structures to energy minimization via AMBER relaxation.
- 6 ****Import and Visualize Minimized Structures**** - Import the minimized structures into PyMOL2. - Use PyMOL2 to produce representative models of the structures.