

Jan 24, 2024

Protocol for preparing post-mortem tissue for intracellular neuromelanin quantification in H&E-stained sections.

 [Nature Communications](#)

 In 3 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7pk1qgwz/v1>

Marta Gonzalez-Sepulveda^{1,2}, Thais Cuadros^{1,2}, Miquel Vila^{1,3,4,5}

¹Neurodegenerative Diseases Research Group, Vall d'Hebron Research Institute (VHIR)-Network Center for Biomedical Research in Neurodegenerative Diseases (CIBERNED), 08035 Barcelona, Spain;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

⁴Institute of Neurosciences, Autonomous University of Barcelona (INc-UAB), Bellaterra, Barcelona, Spain;

⁵Catalan Institution for Research and Advanced Studies (ICREA), Barcelona, Spain

Vilalab Public

Nuria



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7pk1qgwz/v1>

External link: <https://doi.org/10.1038/s41467-024-53168-7>

Protocol Citation: Marta Gonzalez-Sepulveda, Thais Cuadros, Miquel Vila 2024. Protocol for preparing post-mortem tissue for intracellular neuromelanin quantification in H&E-stained sections.. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g7pk1qgwz/v1>

Manuscript citation:

Laguna, A., Peñuelas, N., Gonzalez-Sepulveda, M. *et al.* Modelling human neuronal catecholaminergic pigmentation in rodents recapitulates age-related neurodegenerative deficits. *Nat Commun* **15**, 8819 (2024). <https://doi.org/10.1038/s41467-024-53168-7>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 24, 2024

Last Modified: January 24, 2024

Protocol Integer ID: 94079

Keywords: intracellular neuromelanin quantification, stained section

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020505

Abstract

Protocol for preparing post-mortem tissue for intracellular neuromelanin quantification in H&E-stained sections

Tissue Processing

- 1 Formalin-fixed paraffin-embedded tissue blocks from the pontine and midbrain regions were sectioned at 5µm, collected onto adhesive microscope slides and allowed to dry at room temperature for 24 hours.
- 2 Slides were incubated in the oven at 60°C for 30 minutes to melt the paraffin.
- 3 To remove the paraffin, sections were submerged in Xylene (Panreac Applichem 211769.2714) for 3 × 3 minutes, followed by rehydration in decreasing ethanol concentrations (100% ethanol for 2 × 5 minutes, 95% ethanol for 2 × 5, 70% ethanol for 2×5 minutes) and distilled H₂O for 5 minutes.
- 4 Standard hematoxylin-eosin (H&E) staining.

Image Acquisition

- 5 Sections were scanned using 20x objective (NA=0.8) with pre-set focusing and exposure parameters for optimal neuromelanin signal quality with an automated Slide Scanner Olympus (SLIDEVIEW VS200, Tokyo, Japan).
- 6 Identical parameters were applied to each scanned section to ensure consistency in capturing neuromelanin.