



Nov 20, 2025

Recipe 3 for Avibacterium paragallinarum Strains from Yangzhou

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw4rmwlmk/v1>

wujunfeng¹

¹Nanjing University of Posts and Telecommunications



wujunfeng

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw4rmwlmk/v1>

Protocol Citation: wujunfeng 2025. Recipe 3 for Avibacterium paragallinarum Strains from Yangzhou. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvw4rmwlmk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 20, 2025



Last Modified: November 20, 2025

Protocol Integer ID: 233038

Keywords: avibacterium paragallinarum strain, deep16sculture ai framework for avibacterium paragallinarum strain, yangzhou this recipe, yangzhou, recipe, yangzhou university

Abstract

This recipe was generated by our Deep16SCulture AI framework for *Avibacterium paragallinarum* Strains from Yangzhou University.

Materials

- Bacto Agar: 15,000 mg
- Glucose: 5,000 mg
- L-Cysteine: 500 mg
- Tryptone: 30,000 mg
- Yeast extract: 20,000 mg
- Hemin solution: 25 ml
- NAD: 10 mg (added)
- Water: 1,000 ml

Optional: Glucose.



- 1 Add hemin/NAD post-autoclaving.