

Protocol status: Working

Created: November 01, 2017

Last Modified: June 22, 2018

Protocol Integer ID: 8513

Abstract

Thrombus formation was measured by the viscoelastic test: ROTEM (Tem International GmbH c/o Dutch Affiliate, Tilburg, The Netherlands). Standard assays were used according to the manufacturer's recommendations: EXTEM (ref.: 503-05, Tem International GmbH c/o Dutch Affiliate, Tilburg, The Netherlands), FIBTEM (ref.: 503-06, Tem International GmbH c/o Dutch Affiliate, Tilburg, The Netherlands), and HEPTEM (ref.: 503-09, Tem International GmbH c/o Dutch Affiliate, Tilburg, The Netherlands). All samples were measured within 1 h after blood collection. Furthermore, by means of EXTEM and FIBTEM results, the contribution of platelet count to the thrombus formation was calculated as the PLTEM parameter [22].

