



Jan 21, 2022

Version 1

Sanger sequencing of SARS-CoV-2 Spike protein V.1

 [PLOS One](#)

 Peer-reviewed method

DOI

<https://dx.doi.org/10.17504/protocols.io.bx6kprcw>

Tiago S. Salles*¹, Andrea Cony Cavalcanti*², Fabio Burack da Costa¹, Renata Campos Azevedo¹

¹Laboratório de Interação Vírus-Célula, Departamento de Virologia, Inst. de Microbiologia Paulo de Góes, Universidade Federal do Rio de Janeiro, Rio de Janeiro, Brazil;

²Laboratório Central de Saúde Pública Noel Nutels - LACEN-RJ, Rio de Janeiro, Brazil

PLOS ONE Lab Protocols

Low-cost, high-quality ...



fabio_burack

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bx6kprcw>

External link: <https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0262170>

Protocol Citation: Tiago S. Salles*, Andrea Cony Cavalcanti*, Fabio Burack da Costa, Renata Campos Azevedo 2022. Sanger sequencing of SARS-CoV-2 Spike protein. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bx6kprcw>

Manuscript citation:

Salles TS, Cavalcanti AC, da Costa FB, Dias VZ, de Souza LM, et al. (2022) Genomic surveillance of SARS-CoV-2 Spike gene by sanger sequencing. PLOS ONE 17(1): e0262170. <https://doi.org/10.1371/journal.pone.0262170>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: September 13, 2021

Last Modified: October 06, 2021

Protocol Integer ID: 53164

Keywords: Sanger sequencing, SARS-CoV-2, Spike protein, sanger sequencing of sar, sanger sequencing, urgent surveillance of sar, sequencing methodology, sar, spike protein, protocol for complete nucleotide, mutation, pandemic, vaccine immunity, complete nucleotide

Abstract

The SARS-CoV-2 responsible for the ongoing COVID pandemic reveals particular evolutionary dynamics and an extensive polymorphism, mainly in Spike protein. Monitoring the S protein mutations is crucial for successful controlling measures and detect variants that can evade vaccine immunity. Even after the costs reduction imposed by the pandemic, the new generation sequencing methodologies remain unavailable to many scientific groups. Therefore, to support the urgent surveillance of SARS-CoV-2 S protein, this work describes a protocol for complete nucleotide sequencing of the S protein using the Sanger technique. Thus, any laboratory with experience in sequencing can adopt this protocol.

Guidelines

This protocol works well with recent extracted RNA samples, ideally with ct values lower than 20.



Materials

- Thermal cycler
- PCR tubes 0.2mL
- Filter pipette tips: 1-10 μ L+ 10-100 μ L
- Micropipettes: 1-10 μ L+ 10-100 μ L
- Superscript III one-step RT-PCR kit (Invitrogen, Carlsbad, CA, USA)
- Primers
- Nuclease free water
- Extracted RNA from SARS-CoV-2 positive samples with low ct values
- Agarose
- TAE buffer (Tris-acetate-EDTA)
- LGC biotecnologia - Blue Green loading Dye I
- Horizontal Electrophoresis cube
- UV Transilluminator
- Nanodrop spectrophotometer

1 RT-PCR

- Program the thermal cycler before setting up the reaction. The thermal cycler should be preheated to 45–60°C.
- Keep all components, reaction mixes, and samples on ice. After preparation of the samples, transfer them to the preheated thermal cycler and immediately start the RT-PCR program.
- Reaction mix should be prepared according to table 1.

	A	B	C
	Component	Volume	Final concentration
	2X Reaction Mix	12.5 µL	1x
	Sense primer (10 µM)	1.75 µL	0.7 µM
	Anti-sense primer (10 µM)	1.75 µL	0.7 µM
	SuperScript [™] III RT/Platinum [™] Taq High Fidelity Enzyme Mix[1]	0.5 µL	
	Template RNA	5 µL	
	Nuclease free water	to 25 µL	

Table 1 - RT-PCR master mix

- **RT-PCR cycle is described below:**

60°C	1 min		
50°C	45 min		Reverse Transcription 1 cycle
94°C	2 min		



95 °C	15 s	Denaturation		
53 °C	30 s	Annealing		40 cycles
68 °C	60 s	Extension		

68 °C	7 min	Final extention		1 cycle
4 °C	∞			

■ **Primers used:**

			Position	Melting
	temperature	Size		
P1 forward		GTTTGTTTTTCTTGTTTTATT	(21551-21574)	43.5 °C
		923bp		
P1 reverse		ACAGTGAAGGATTTCAACGTACAC	(22450-22474)	55.3 °C
P2 forward		CGTGATCTCCCTCAGGGTTTT	(22190-22211)	56.8
	°C	620bp		
P2 reverse		TCAGCAATCTTTCCAGTTTGCC	(22810-22832)	56.1
	°C			
P3 forward		GTAATTAGAGGTGATGAAGTCAGA	(22751-22775)	51.8 °C
		892bp		
P3 reverse		ACATAGTGTAGGCAATGATGGA	(23621-23643)	53.6
	°C			
P4 forward		CTTGCGTGTTTATTCTACAG	(23445-23466)	51.4
	°C	979bp		
P4 reverse		GCTTGTGCATTTTGTTGACC	(24403-24424)	55.6
	°C			
P5 forward		AGACTCACTTTCTTCCACAGCA	(24355-24377)	56.1
	°C	342bp		
P5 reverse		AGATGATAGCCCTTTCCACA	(24699-24719)	53.3
	°C			
P6 forward		TTCTGCTAATCTTGCTGCTACT	(24610-24632)	54.0
	°C	766bp		
P6 reverse		GTTTATGTGTAATGTAATTTGACTCC	(25348-25372)	50.7
	°C			
P7 forward		TAGAGAAAACAACAGAGTT	(21492-21511)	45.6
	°C	712bp		



P7	reverse	TGAGGGAGATCACGCACTAA	(22184-22204)	55.4
°C				
P8	forward	TTCTGCTAATCTTGCTGCTACT	(24610-24632)	54.0
°C	825bp			
P8	reverse	CCTTGCTTCAAAGTTACAGTTCCA	(25409-25433)	55.6
°C				

2 Agarose gel Electrophoresis

- Dissolve agarose in 1.X TAE Buffer (40 mM Tris-acetate, 1 mM EDTA) to a final concentration of 1.5% agarose.
- Heat the solution in a microwave and let it cool at room temperature.
- Apply the agarose gel to the casting tray with the appropriate well comb and let it solidify.
- Mix 5 µL of PCR product with 5µL of loading buffer and 1µL Blue Green Loading dye I (LGC Biotecnologia), apply it to the gel, and run the electrophoresis at 120V.
- Visualize the gel with UV Transilluminator

3 Preparing Samples for sequencing

- Sequencing procedures are performed according to BigDye™ Terminator v3.1 Cycle Sequencing Kit (Applied Biosystems) guidelines.
- Measure DNA Concentration with a nanodrop spectrophotometer
- Dilute template to 200 ng/µl with nuclease-free water
- Dilute primers to 1 µM with nuclease-free water. Only one primer is used for each sequencing reaction, leading to two reactions per sample. Each reaction will need 1µl of diluted primer.
- Label and ship the samples and primers according to the sequencing service provider guidelines.

sequencing reaction:

- BigDye™ Terminator 3.1 Ready Reaction Mix - 8 µL
- primer (1 µM) - 6.4 µL
- Template (200ng) - 2 µL
- Deionized water (RNase/DNase-free) to - 20 µL

sequencing cycle:



96°C	1 min	Incubation	1 cycle
96 °C	10 s	Denaturation	
50 °C	05 s	Annealing	25 cycles
60 °C	4 min	Extension	
4 °C	∞		1 cycle

4 Sequence analysis

- Upon receiving the electropherograms, edit them with proper programs like chromas or Bioedit.
- For better editing accuracy, align the forward sequence with the reverse complement of the reverse sequence.
- Use BLAST search to confirm if the sequenced product corresponds to the desired target.
- Create contigs by aligning the planned overlaps between each target fragment and form one consensus sequence covering the full ORF of SARS-CoV-2 Spike protein
- Edited sequences can now be analyzed with the CoVsurver mutations app, provided by GISAID, to trace the mutation patterns of each sample and study their effects on protein structure.
- Sequences can also be deposited in the GISAID database of SARS-CoV-2.