

Apr 21, 2018

 SOB Media

 Forked from [SOB Media](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.pnfdmbn>

New England Biolabs¹

¹New England Biolabs



Diego Bonatto

Centro de Biotecnologia da UFRGS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.pnfdmbn>

External link: <https://www.neb.com/protocols/2012/06/21/making-your-own-electrocompetent-cells>

Protocol Citation: New England Biolabs 2018. SOB Media. protocols.io <https://dx.doi.org/10.17504/protocols.io.pnfdmbn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: April 21, 2018

Last Modified: April 21, 2018

Protocol Integer ID: 11687

Keywords: sob media



- 1 2% tryptone
- 2 0.5% yeast extract
- 3 10 mM NaCl
- 4 2.5 mM KCl
- 5 10 mM MgCl₂
- 6 10 mM MgSO₄