

Aug 29, 2019

Version 2

## Vandy - Mouse Myocardial Infarction V.2

DOI

[dx.doi.org/10.17504/protocols.io.6xghfjw](https://dx.doi.org/10.17504/protocols.io.6xghfjw)



Lin Zhong<sup>1</sup>, Jeffrey Rottman<sup>1</sup>, Chee Lim<sup>1</sup>

<sup>1</sup>Vanderbilt University

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6xghfjw>

External link: <https://mmpc.org/shared/document.aspx?id=227&docType=Protocol>



**Protocol Citation:** Lin Zhong, Jeffrey Rottman, Chee Lim 2019. Vandy - Mouse Myocardial Infarction. **protocols.io**

**<https://dx.doi.org/10.17504/protocols.io.6xghfjw>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 29, 2019

**Last Modified:** August 29, 2019

**Protocol Integer ID:** 27336

**Keywords:** cardiovascular mortality, myocardial infarction, mouse myocardial infarction summary, many mouse models relevant to human cardiovascular disease, surgical induction of myocardial injury, myocardial injury, corresponding process of heart injury recovery, heart injury recovery, reperfusion injury, outcome from myocardial infarction, human cardiovascular disease, cardiovascular disease, coronary artery, reperfusion, common cause of cardiovascular mortality, heart, ischemia

## Abstract

### Summary:

The most common cause of cardiovascular mortality in man is the outcome from myocardial infarction. Accordingly, it is necessary to study the corresponding process of heart injury recovery in many mouse models relevant to human cardiovascular disease. This protocol describes the surgical induction of myocardial injury via 1) permanent ligation of a coronary artery in the heart (myocardial infarction), 2) cryoinjury, and 3) transient occlusion of a coronary artery followed by reperfusion (ischemia-reperfusion injury).

## Materials

### Reagents and Materials:

- Pentobarbital
- Buprenorphine
- Betadine
- 70% alcohol
- PE-50 tube
- Cryo-probe
- 6-0 sutures
- 7-0 sutures
- 8-0 sutures

## Troubleshooting

- 1 Mice are anesthetized with pentobarbital (50 mg.kg, IP).
- 2 The ventral neck and left parasternal region is shaved and disinfected with Betadine followed by 70% alcohol.
- 3 The mouse is positioned supine on a heating pad and a small incision is made through the skin underlying the trachea.
- 4 The trachea is exposed, a small puncture is made in the trachea, and endotracheal intubation is performed using a PE-50 tube.
- 5 The endotracheal tube is connected to a small rodent ventilator (Harvard Apparatus) for mechanical ventilation of the mouse.
- 6 ECG electrode leads are placed subcutaneously to monitor the ECG during myocardial infarction(MI).
- 7 With the use of a surgical microscope, a left thoracotomy is performed and the fourth intercostal space is entered using scissors and blunt dissection.
- 8 An 8-0 silk suture is placed through the myocardium into the anterolateral left ventricular wall around the left anterior descending (LAD) coronary artery.
- 9 **MI:** The suture is tied resulting in permanent ligation of the LAD and subsequent MI is confirmed via the use of the ECG monitor. Sham undergo the same procedure with the exception that the suture is not tied around the LAD, thus no ligation occurs.

**Ischemia-Reperfusion Injury:** A sterile PE-50 tube is placed against the LAD and the suture is tied against the tube resulting in ligation of the LAD. During the ischemic period, the surgical area is covered with sterile gauze soaked in warm saline. Following the ischemic period, the tube is removed and the LAD ligature cut resulting in restoration of blood flow and reperfusion.

**Cryoinfarction:** A small cryo-probe (3 mm diameter) is cooled in liquid nitrogen for 2 minutes and placed in a circumscribed region of the free LV wall (30 seconds) in order to achieve a reproducible transmural lesion.



- 10 The chest is closed in layers with 7-0 sutures.
- 11 The mouse is gradually weaned from the ventilator to resolve any possible pneumothorax.
- 12 Once spontaneous respiration resumes, the endotracheal tube is removed and the trachea incision is closed with 8-0 suture. The skin is then closed with 6-0 suture.
- 13 The mouse is maintained on the heating pad until fully recovered from anesthesia.
- 14 Buprenorphine is administered SC immediately following surgery and every 8-12 hr for 72 hr.