

Feb 16, 2026

Whole-Brain Axonal Projection Quantification Using ABBA and QuPath

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwn2o9vmk/v1>

Cristian González-Cabrera¹

¹LIN



Cristian González-Cabrera

Leibniz Institute für Neurobiologie

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwn2o9vmk/v1>

Protocol Citation: Cristian González-Cabrera 2026. Whole-Brain Axonal Projection Quantification Using ABBA and QuPath. protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvwn2o9vmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 16, 2026

Last Modified: February 16, 2026

Protocol Integer ID: 243362

Keywords: brain axonal projection quantification, brain quantification of axonal projection, allen mouse brain atlas, brain quantification, axonal projection, brain section, standardized anatomical region boundary, fluorescence intensity measurement in qupath, registration with abba, using abba, qupath, fluorescence intensity measurement, fluorescence intensity, qupath this protocol, abba

Abstract

This protocol describes whole-brain quantification of axonal projections using atlas-based registration with ABBA [51] followed by fluorescence intensity measurement in QuPath [52]. Brain sections are first aligned to the Allen Mouse Brain Atlas using ABBA to obtain standardized anatomical region boundaries. Fluorescence intensity is then quantified in QuPath and background-corrected. Data are aggregated per animal to avoid pseudoreplication.

Materials

- ABBA (Allen Brain Atlas Alignment) [51].
- Allen Mouse Brain Atlas reference.
- QuPath (fluorescence analysis mode) [52].
- Spreadsheet software for aggregation.



Procedure

- 1 Import section images into ABBA.
- 2 Align sections to the Allen Mouse Brain Atlas using the ABBA workflow.
- 3 Manually verify alignment quality.
- 4 Export atlas-defined region boundaries or masks for downstream quantification.
- 5 Open registered section images in QuPath.
- 6 Confirm correct pixel calibration.
- 7 Import atlas-derived ROI boundaries if exported from ABBA.
- 8 Use atlas-aligned boundaries to define anatomically consistent ROIs.
- 9 Treat hemispheres separately if hemispheric analysis is required.
- 10 Apply ROI boundaries consistently across animals.
- 11 Define a local background ROI on the same section.
- 12 Select an area without specific axonal labeling.
- 13 Maintain consistent background selection criteria across sections.



- 14 Extract mean fluorescence intensity for each ROI.
- 15 Extract mean intensity for corresponding background ROI.
- 16 Compute background-corrected intensity as: ROI mean - background mean.
- 17 Average corrected values per region per animal.
- 18 Retain hemisphere-specific values if required before final averaging.
- 19 Ensure each animal contributes one value per region.

Protocol references

[51] Chiaruttini, N., Castoldi, C., Reque, L., et al. ABBA+BraiAn, an integrated suite for whole-brain mapping, reveals brain-wide differences in immediate-early genes induction upon learning. *Cell Reports* (2025). <https://doi.org/10.1016/j.celrep.2025.115876>

[52] Bankhead, P., Loughrey, M. B., Fernández, J. A., et al. QuPath: Open source software for digital pathology image analysis. *Scientific Reports* 7, 16878 (2017). <https://doi.org/10.1038/s41598-017-17204-5>



Acknowledgements

Outputs

- Background-corrected fluorescence intensity per region per animal.
- Hemisphere-specific values when applicable.
- Exportable dataset for statistical analysis.

Critical Steps

- Verify ABBA atlas alignment accuracy before quantification.
- Use identical imaging settings within datasets.
- Maintain consistent ROI and background selection criteria.
- Avoid treating multiple sections from the same animal as independent samples.