

Jan 06, 2020

Version 2

4-methylumbelliferyl heptanoate (MUH) - Cell Viability Assay V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.bazjif4n>


Peter Vangheluwe¹, Shaun Martin¹, Mujahid Azfar¹

¹KU Leuven



Peter Vangheluwe

KU Leuven



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bazjif4n>

Protocol Citation: Peter Vangheluwe, Shaun Martin, Mujahid Azfar 2020. 4-methylumbelliferyl heptanoate (MUH) - Cell Viability Assay. [protocols.io](https://dx.doi.org/10.17504/protocols.io.bazjif4n) <https://dx.doi.org/10.17504/protocols.io.bazjif4n>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 06, 2020

Last Modified: May 31, 2024

Protocol Integer ID: 31499

Keywords: Cell Viability Assay, MUH assay, 4-methylumbelliferyl heptanoate, cell-mediated cytotoxicity, cell proliferation, ASAPCRN, cell viability assay, methylumbelliferyl heptanoate, mediated cytotoxicity, principle behind this assay, assay, fluorimetric assay, hydrolyse the fluorochrome, viable cell, fluorochrome, muh reagent, activity of esterase, cell, esterase

Abstract

The 4-methylumbelliferyl heptanoate (*MUH*) assay is a A fluorimetric assay used to study cell-mediated cytotoxicity. The principle behind this assay employs the activity of esterases found in viable cells. These esterases then hydrolyse the fluorochrome in the MUH reagent which can then be read using a plate reader.

Attachments



[Cell Viability - MUH...](#)

361KB

