

Jan 29, 2016

# Estimation of viral-induced phytoplankton mortality using the modified dilution method

DOI

dx.doi.org/10.17504/protocols.io.dpf5jm

Susan A. Kimmance and Corina P. D. Brussaard<sup>1</sup>

<sup>1</sup>Manual of Aquatic Viral Ecology

**VERVE Net** 



**Anna Noordeloos** 

# Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: https://dx.doi.org/10.17504/protocols.io.dpf5jm

External link: http://www.aslo.org/books/mave/MAVE\_065.pdf

Protocol Citation: Susan A. Kimmance and Corina P. D. Brussaard 2016. Estimation of viral-induced phytoplankton mortality using the modified dilution method. protocols.io <a href="https://dx.doi.org/10.17504/protocols.io.dpf5jm">https://dx.doi.org/10.17504/protocols.io.dpf5jm</a>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: August 26, 2015

Last Modified: February 27, 2018

Protocol Integer ID: 1479

**Keywords:** estimation of viralinduced phytoplankton mortality, viral lysis rates of natural phytoplankton, viralinduced phytoplankton mortality, partition phytoplankton mortality, phytoplankton mortality, induced phytoplankton mortality, phytoplankton growth rate, subsequent impact on phytoplankton growth rate, manual of aquatic viral ecology, viral lysis by dilution, aquatic viral ecology, viral lysis rate, phytoplankton, natural phytoplankton, dilution method, dilution, dilution assay, modified dilution method, viral lysi

### **Abstract**


The modified dilution assay aims to partition phytoplankton mortality into virus- versus grazing-induced fractions and has previously been applied to several different environments to determine viral lysis rates of natural phytoplankton. The method involves creating a gradient of both grazing and viral lysis by dilution with different proportions of grazer- and virus-free filtrate, and assessing the subsequent impact on phytoplankton growth rates.

Kimmance, S. A., and C. P. D. Brussaard. 2010. Estimation of viralinduced phytoplankton mortality using the modified dilution method, p. 65–73. In S. W. Wilhelm, M. G. Weinbauer, and C. A. Suttle [eds.], Manual of Aquatic Viral Ecology. ASLO.



### Guidelines

Experiments are set up based on the original protocol of the modified dilution method of Evans et al. (2003), and the following method is appropriate for both natural and cultured phytoplankton-virus systems (Fig. 1). The approach is based around combining mesoplankton-free seawater (whole water) with either grazer-free or virusfree diluents. To prevent contamination of seawater samples by handling, vinyl gloves are worn throughout the water collection and experimental setup. To create mesoplankton-free whole water for experiments, seawater is gently siphoned through 200-µm mesh (or reverse sieved, Baudoux et al. 2006) using silicone tubing into clean, polypropylene carboys. Additional (reverse) prefiltration may be required, e.g., during sampling of colonial phytoplankton (Baudoux et al. 2006) or in productive, coastal, or sediment-filled waters, to prevent clogging during the later fine-scale filtration. Gentle filtration and avoidance of air bubbling are crucial to prevent destruction of the grazer, virus, and phytoplankton populations. To protect the phytoplankton communities within the seawater sample from excessive light (light-shock), experiments using natural phytoplankton populations are best set up predawn in dimmed-light conditions using light-proofed carboys. It is also recommended that experiments be performed at the same time of day because of the synchronicity of phytoplankton cell division and potential diel effects on viral infection processes. Filtration and handling should be conducted at in situ temperature. Handling time once the seawater has been collected should be kept to a minimum before the start of the experimental incubation. Ideally experimental setup should be within an hour of water collection if possible.



**Figure 1**: Modified dilution assay experimental design. Mesoplankton-free whole water (1) is combined with either <0.45-μm filtrate (2) or kDa filtrate (3) in the correct proportions to create the parallel t0 dilution series: <0.45μm series, with reduced grazing mortality (4) kDa series, with reduced grazing and viral mortality (5). Replicate



sample bottles from the <0.45- $\mu$ m and kDa dilution series are then created (6 and 7) and incubated under experimental conditions.

#### References

Baudoux, AC., A. A. M. Noordeloos, M. J. W. Veldhuis, and C. P. D. Brussaard. 2006. Virally induce Phaeocystis globosa during two spring blooms in temperate coastal waters. Aquat. Microb. Ecol. 44	-
——, M. J. W. Veldhuis, H. J. Witte, and C. P. D. Brussaard. 2007. Viruses as mortality agents of picophytoplankton in the deep chlorophyll maximum layer during IRONAGES III. Limnol. Oceanogr. 5	2:2519-2529.
——, ——, A. A. M. Noordeloos, G. van Noort, and C. P. D. Brussaard. 2008. Estimates of virus-induced mortality of picophytoplankton in the North Sea during summer. Aquat. Microb. Ecol. 52:69	
Evans, C., S. D. Archer, S. Jacquet, and W. H. Wilson. 2003. Direct estimates of the contribution of v microzooplankton grazing to the decline of a Micromonas spp. population. Aquat. Microb. Ecol. 30:2	-
Kimmance, S. A., W. H. Wilson, and S. D. Archer. 2007. Use of the modified dilution technique to esti versus grazing mortality of phytoplankton: limitations associated with method sensitivity in natural v Microb. Ecol. 49:207-222.	

## **Troubleshooting**

### Before start

Prior to each experiment 0.2– $0.45~\mu m$  filters, silicon tubing, carboys, and bottles should be acid-washed in 5% HCL and rinsed thoroughly with Milli-Q water and then twice with the sample. The 0.2–0.45- $\mu m$  and kDa filters should be flushed prior to use with freshly prepared Milli-Q and the first few liters of both filtrates discarded. Ensure that all air bubbles are removed from the filter capsule and tubing during filtration by opening the filter valves. To speed up this process you may use two filters in parallel.



To create the grazer-free diluent whole seawater is gravity filtered through acid-washed 0.2–0.45 μm filters (e.g., PALL Acropak™ Supor® membrane capsules) into a clean carboy.

#### Note

A new filter should be used for each dilution experiment.

After gravity filtration half of the grazer-free diluent is passed through a tangential flow filtration system with kilodalton (kDa) pore size to create the virus-free diluent.

#### Note

The designated pore size is variable but typically between 10–100 kDa (Evans et al. 2003; Baudoux et al. 2006; Kimmance et al. 2007). This pore size range is sufficient for removal of both bacteriophage and larger algal viruses such as *Emiliania huxleyi* viruses. The effectiveness of the modified dilution technique is dependent on the efficiency of the 0.2–0.45-µm filtration step to create a gradient of grazing pressure and the kDa filtration step to remove viruses from the 0.2–0.45-µm filtrate, thus creating an additional viral gradient. Assessing the differences in virus abundance between the diluents in comparison with natural samples determines how effective the additional filtration step was at removing viruses and thus producing a gradient of viral pressure.

- The grazer-free and virus-free diluents are added to 10-L polycarbonate bottles in the correct proportions to create the parallel  $t_0$  dilution series, e.g., 20%, 40%, 70%, and 100% whole water.
- 4 The mesoplankton-free whole water is then gently added by siphoning.
- From each of these t<sub>0</sub> bottles, triplicate 1-L polycarbonate bottles were rinsed twice and then gently filled by siphoning (to minimize physical damage to the grazers, viruses, and phytoplankton populations), ensuring that bottles are filled completely to avoid trapping air bubbles inside upon closure.
- After filling, triplicate experimental bottles are placed randomly into experimental conditions.



#### Note

For natural samples, if in situ incubation is not feasible, then the experimental environment should match the in situ temperature and light conditions (including lightdark period) as closely as possible. Therefore outdoor incubators should be temperature controlled and covered with neutral density screening to match natural light intensity at the sample depth. Laboratory assays should be set up as described above except that viral lysate capable of infecting the chosen host should be added (multiplicity of infection >1) to the cultured host before dilution (Baudoux et al. 2006) and experimental bottles incubated under the appropriate conditions for the chosen cultured virus-host system.

7 For determination of phytoplankton composition and abundance and virus abundance, triplicate subsamples (5 mL) are taken from the 10 L  $_{
m 0}$  dilution bottles and the 0.2– 0.45-μm and kDa diluents.

#### Note

However, an alternative method is to take samples out of each and every 1-L experimental bottle directly at  $t_0$ , and a making sure that there are no air bubbles in the bottles upon closing. Final time-point samples are also taken from every experimental bottle at the end of the experimental period, 24 h.

8 Initial (t<sub>0</sub>) and final (t<sub>24</sub>) phytoplankton composition and abundance estimates are typically determined by analysis of samples using flow cytometry (Evans et al. 2003; Baudoux et al. 2006; Kimmance et al. 2007).

#### Note

If flow cytometry is not available other cell-counting methods could be used (e.g., microscopy or COULTER COUNTER®). In the case of monoalgal blooms, chlorophyll measurements are a possible alternative (Evans et al. 2003). However, this is not recommended because of the highly specific nature of viral infection. It is more desirable to detect changes within specific phytoplankton groups, and with flow cytometry distinct groups can easily be discriminated by differences in fluorescence characteristics.

9 Apparent phytoplankton growth rates ( $\mu$ ,  $d^{-1}$ ) are calculated from each experimental bottle as the changes in abundance during the incubation using the equation:

# $\mu = \ln(P_t/P_o)/t$

#### Note

where  $P_t$  and  $P_0$  are the final and initial measured phytoplankton abundance, respectively, and t is the duration of the experiment.



The actual dilution rate is calculated by dividing the  $t_0$  phytoplankton abundance in each bottle by the averaged abundance of the replicate  $t_0$  100% counts (3–6 bottles).

#### Note

Model 1 linear regression analysis of apparent growth rates against fraction of whole water is applied to each of the dilution experimental series (0.2–0.45 µm and kDa) to estimate instantaneous growth and mortality due to grazing and/or viral lysis. As explained above, the regression coefficient of apparent growth rate versus fraction of whole water for the <0.45-µm dilution series represents the microzooplankton grazing rate  $[m_g=(\mu \text{ in the absence of grazing}-\text{observed net }\mu)=\text{slope}],$  whereas the regression coefficient from the kDa dilution series represents the combined mortality impact of both grazers and viruses  $(m_g+m_v)$ , where  $m_v=\text{viral mortality }(\underline{\textbf{Fig 2}})$ . The viral-induced phytoplankton mortality rate therefore is the difference between these two regression slopes, i.e.,  $m_v=[(m_v+m_g)-m_g]$ . However, a significant rate can be derived only when there is a significant difference between the two mortality slopes. To determine this, experiments that produce significant regression slopes for both the 0.2–0.45-µm and kDa dilution series are further analyzed, e.g., using an F-test. Only then can the difference between the two slopes be considered a significant viral lysis rate.