Sep 07, 2018

biomechanical testing of repaired tendon

DOI

dx.doi.org/10.17504/protocols.io.tadeia6

Susumu Yoneda¹

¹Washington University, Saint Louis

Susumu Yoneda





DOI: dx.doi.org/10.17504/protocols.io.tadeia6

Protocol Citation: Susumu Yoneda 2018. biomechanical testing of repaired tendon. protocols.io <u>https://dx.doi.org/10.17504/protocols.io.tadeia6</u>

License: This is an open access protocol distributed under the terms of the **<u>Creative Commons Attribution License</u>**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: September 05, 2018

Last Modified: September 07, 2018

Protocol Integer ID: 15397

Abstract

This protocol was used for biomechanical testing of repaired tendon.

 All tendon were kept wet with phosphate buffered saline during procedure. The proximal tendon segment was secured with sandpaper and a triangle-toothed grip and the distal phalanx was secured in a custom grip. Distance between grips was 60mm. A random texture was applied to the surface of the tendons using Verhoeff's stain for tracking.

Tendons were preconditioned to a preload of 1 N followed by five triangle wave cycles between 0-1 mm at a rate of 0.3mm/s. Then tendons were pulled in uniaxial tension at 0.3 mm/s until failure on a material testing machine (5866; Instron Corp., Norwood, MA), and strain was tracked optically using video (60Hz) of the test. Failure mode (i.e., suture pullout, suture breakage, or knot failure) was recorded.

Load and strain data were synchronized and used to generate load displacement and load strain curves. The maximum load (N), load to create a 2 mm gap (N), stiffness (N/mm) (the slope of the linear portion of the load–deformation curve), strain at 20 N (%), yield force (N) and a modified version of resilience (N strain) (the area under the load–strain curve up to the yield point) were determined using a custom-written code in MATLAB (Natick, MA). Load to create a gap of 2 mm (a threshold level that leads to decreased repair strength and increased adhesions) between tendon stump ends was calculated by tracking the displacement between points ± 2 mm from the repair site with code in MATLAB.