

Oct 05, 2022

Electrophysiology

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l274mjpg1y/v1>

gurvir.virdi¹

¹UCL Institute of Neurology



gurvir.virdi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l274mjpg1y/v1>

Protocol Citation: gurvirdi 2022. Electrophysiology. protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l274mjpg1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 05, 2022

Last Modified: May 31, 2024



Protocol Integer ID: 70867

Keywords: ASAPCRN, electrophysiology electrophysiology of ipsc, derived mda neuron, mda neuron, electrophysiology electrophysiology, ipsc

Abstract

Electrophysiology of iPSC-derived mDA neurons

Electrophysiology

- 1 Visualized patch-clamp recordings from cell cultures were performed using an infrared differential interference contrast imaging system and a Multipatch 700B amplifier controlled by pClamp 10.2 software package (Molecular Devices, USA).
 - 1.1 For the recordings, a neuronal culture on a glass coverslip was placed in a recording chamber mounted on the stage of an Olympus BX51WI upright microscope (Olympus, Japan).
 - 1.2 The perfusion solution contained the following (in mM): 119 NaCl, 2.5 KCl, 1.3 Na₂SO₄, 2.5 CaCl₂, 26.2 NaHCO₃, 1 NaH₂PO₄, 2 CaCl₂, 2 MgCl₂, 10 glucose (or 22 in some recordings) and was continuously bubbled with 95% O₂ and 5% CO₂, at a pH of 7.4.
 - 1.3 Whole-cell recordings were performed at  32–34 °C ; the patch-clamp pipette resistance was 3–7 MΩ depending on particular experimental conditions.
 - 1.4 Series resistance was monitored throughout experiments using a +5 mV step command, cells with very high series resistance (above 25 MΩ) or unstable holding current were rejected.
 - 1.5 The intracellular pipette solution for voltage-clamp experiments contained (in mM): 120.5 CsCl, 10 KOH-HEPES, 2 EGTA, 8 NaCl, 5 QX-314 Br[−] salt, 2 Na-ATP, 0.3 Na-GTP.
 - 1.6 For current-clamp experiments, the intracellular solution contained (in mM): 126 K-gluconate, 4 NaCl, 5 HEPES, 15 glucose, 1 K₂SO₄×7 H₂O, 2 BAPTA, 3 Na-ATP. The pH was adjusted to 7.2 and osmolarity adjusted to 295 mOsm.
 - 1.7 To isolate response of NMDA receptors we added to a perfusion solution: 50 mM picrotoxin, 20 mM NBQX, 1 mM strychnine, 1 mM CGP-55845, 100 mM MCPG, with zero Mg²⁺.
 - 1.8 To isolate response of GABA_A receptors, we added 50 mM APV, 20 mM NBQX, 1 mM strychnine, 1 mM CGP-55845, 100 mM MCPG. All chemicals were purchased from Tocris Bioscience.
- 2 In the whole-cell (immediately after membrane breakthrough), iPSC-derived neurons were recorded for the resting membrane potential (V_{rest}), membrane capacitance (C_m), the membrane time constant (τ_m), and input resistance (R_{in}), measured from the hyperpolarizing square current pulse steps in current mode
 - 2.1 To assess the firing capability of the cells, a series of sub- and supra-threshold rectangular current pulses were applied to elicit neuronal firing, with a stepwise-



increased stimulus intensity (an increment of 5–10 pA).

2.2 The V_{rest} was set at –60 mV to –70 mV, by injecting a hyperpolarizing bias current where required.

2.3 The analysis of the AP waveform was performed for the first AP only.

2.4 The parameters of individual APs recorded were: the spike amplitude (measured from the threshold to the peak), the threshold value, overshoot and the spike width (duration at half-maximal amplitude), the rates of depolarisation and repolarisation phases.