

Jun 14, 2019

# Golden Gate plasmids used for transfection of *Perkinsus marinus*

DOI

[dx.doi.org/10.17504/protocols.io.37egrje](https://dx.doi.org/10.17504/protocols.io.37egrje)



Elin Einarsson<sup>1</sup>

<sup>1</sup>Waller Lab, Department of Biochemistry, University of Cambridge

Protist Research to Opti...



Elin Einarsson

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.37egrje>

**Document Citation:** Elin Einarsson 2019. Golden Gate plasmids used for transfection of *Perkinsus marinus*. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.37egrje>

**License:** This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Created:** June 14, 2019

**Last Modified:** June 14, 2019

**Document Integer ID:** 24518

**Keywords:** plasmids for transfection, transfection of perkinsus marinus, golden gate plasmid, using golden gate cloning, golden gate cloning, parasite perkinsus marinus, transfected cell, created plasmid, transfection, moe membrane protein, perkinsus marinus, bleomycin, puromycin

## Abstract

We have created plasmids for transfection of the parasite Perkinsus marinus using Golden Gate cloning. We have explored the possibility to either FACS sort (GFP/mCherry fused to MOE membrane protein) or using drug selection (bleomycin, blasticidin S or puromycin) to enrich transfected cells.

## Attachments



GoldenGate plasmids ...

303KB

## Troubleshooting

